

Synoptic CB: Porewater DIC

September 2025 Samples

2025-10-21

Contents

| | | |
|------|--|---|
| 0.1 | Import Data Functions | 2 |
| 0.2 | Import Sample Data | 2 |
| 0.3 | Assessing Standard Curves | 3 |
| 0.4 | CRM Check - Don't run chunk if no CRMs run | 4 |
| 0.5 | Assess Check Standards | 4 |
| 0.6 | Assess Blanks | 5 |
| 0.7 | Assess Duplicates - No duplicates included on this run | 6 |
| 0.8 | Sample Flagging - Are samples Within the range of the curve? | 6 |
| 0.9 | Visualize Data by Plot | 7 |
| 0.10 | Convert data from mg/L to uMoles/L | 8 |
| 0.11 | Check to see if samples run match metadata & merge info | 8 |
| 0.12 | Export Processed Data | 8 |

```

##Setup - Change things here & write any notes

#identify section
cat("Setup Information")

## Setup Information

##### Run information - PLEASE CHANGE
Date_Run = "09/22/25" #Date that instrument was run
Run_by = "Stephanie J. Wilson" #Instrument user
Script_run_by = "Stephanie J. Wilson" #Code user
run_notes = " " #any notes from the run
samples <- c("GCW", "GWI", "MSM", "SWH") #whatever identifies your samples within the same names
samples_pattern <- paste(samples, collapse = "|")
#samples_pattern <- "GCW" #use this instead of the line above if you have only one site code
chks_name = "Chk_Std_" #what did you name your check standards?
crm_name = "CRM|crm" #what did you name your CRMS?

##### File Names - PLEASE CHANGE
#file path and name for raw summary data file
raw_file_name = "Raw Data/COMPASS_SynopticCB_PW_DIC_202509.txt"

#file path and name for raw all peaks file
raw_allpeaks_name = "Raw Data/COMPASS_SynopticCB_PW_DIC_202509_allpeaks.txt"

#file path and name of processed data file
processed_file_name = "Processed Data/COMPASS_SynopticCB_PW_Processed_DIC_202509.csv"

##### Log Files - PLEASE CHECK
#downloaded metadata csv - downloaded from Google drive as csv for this year
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2025.csv"

#qaqc log file path for this year
Log_path = "Raw Data/COMPASS_Synoptic_DIC_QAQClog_2025.csv"

##Set Up Code
##Read in metadata and create similar sample IDs for matching to samples

```

0.1 Import Data Functions

0.2 Import Sample Data

```

## Import Sample Data

## New names:
## * ` ` -> '...14'

## # A tibble: 6 x 3
##   sample_name          ic_raw run_datetime
##   <chr>                <dbl> <chr>

```

```

## 1 202509_GCW_TR_LysB_10cm    4.59 9/23/2025 2:34:07 AM
## 2 202509_GCW_WC_LysA_10cm 150    9/23/2025 2:49:08 AM
## 3 202509_GCW_WC_LysA_20cm 147.   9/23/2025 3:05:24 AM
## 4 202509_GCW_WC_LysA_45cm 182.   9/23/2025 3:26:43 AM
## 5 202509_GCW_WC_LysB_10cm 162.   9/23/2025 3:42:42 AM
## 6 202509_GCW_WC_LysB_20cm 158.   9/23/2025 4:01:13 AM

```


0.3 Assessing Standard Curves

```
## Assess the Standard Curves
```

```

## New names:
## `geom_smooth()` using formula = 'y ~ x'
## * `-->` '...18'

```




```

## Warning: Removed 4 rows containing missing values or values outside the scale range
## ('geom_point()').

## Warning: Removed 4 rows containing missing values or values outside the scale range
## ('geom_line()').

```

Slope Drift Assessment



```
## [1] "IC Curve r2 GOOD"
```

0.4 CRM Check - Don't run chunk if no CRMs run

```
## Assess the CRMs

## New names:
## * `` -> '...14'

## [1] "IC crm has a % Difference <25% of expected - PROCEED"


## Run mean = 22.68

## Expected = 22.19
```

0.5 Assess Check Standards

```
## Assess the Check Standards

## New names:
## * `` -> '...14'
```




```
## [1] ">60% of IC Check Standards are within range of expected concentration"
```

0.6 Assess Blanks

```
## Assess Blanks
```

```
## New names:  
## * ' ' -> '...14'
```

```
## [1] ">60% of Carbon Blank concentrations are lower 25% quartile of samples"
```



```
## carbon blanks:
```


```
## [1] 1.34069
```

0.7 Assess Duplicates - No duplicates included on this run

0.8 Sample Flagging - Are samples Within the range of the curve?

```
## Sample Flagging
```

C: Grey = Within Range of Curve




0.9 Visualize Data by Plot

```
## Visualize Data
```

```
## Warning in rbind(c("202509", "GCW", "TR", "LysB", "10cm"), c("202509", "GCW", :  
## number of columns of result is not a multiple of vector length (arg 11)
```

Samples: DIC



0.10 Convert data from mg/L to uMoles/L

0.11 Check to see if samples run match metadata & merge info

```
## Check Sample IDs with Metadata
## All sample IDs are present in metadata.
```

0.12 Export Processed Data

```
## Export Processed Data

## # A tibble: 6 x 18
##   Project      Region Site  Zone  Replicate Depth_cm Sample_ID Year Month Day
##   <chr>        <chr>  <chr> <fct> <chr>       <int> <chr>    <int> <int> <int>
## 1 COMPASS: Sy~ CB    MSM    UP     B          10 202509_M~ 2025    9    11
## 2 COMPASS: Sy~ CB    MSM    UP     C          10 202509_M~ 2025    9    11
## 3 COMPASS: Sy~ CB    MSM    UP     C          20 202509_M~ 2025    9    11
## 4 COMPASS: Sy~ CB    GWI    UP     A          20 202509_G~ 2025    9    17
## 5 COMPASS: Sy~ CB    GWI    UP     A          45 202509_G~ 2025    9    17
## 6 COMPASS: Sy~ CB    GWI    UP     B          10 202509_G~ 2025    9    17
## # i 8 more variables: Time <chr>, Time_Zone <chr>, ic_mgL <dbl>, ic_uM <dbl>,
## #   ic_flag <chr>, Analysis_runtime <chr>, Run_notes <chr>, Field_notes <chr>

#end
```