

Dionex_COMPASS_September2022

Stephanie J. Wilson

2023-01-23

Daily Set up

Read in and Format the raw data - change wd & file names

```
# SULFATE DATA:

## Read in raw data file from Dionex - copied and saved as a txt
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202305_S04.txt",sep='\t' , header=T, skip=3)
head(Sdat)

##      X          X.1          X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank Unknown 4.887  0.0084   22.74  0.0132   0.06
## 2 2 Lab Blank Unknown 4.907  0.0091   23.29  0.0143   0.06
## 3 3 Lab Blank Unknown 4.907  0.0094   32.59  0.0148   0.06
## 4 4 Lab Blank Unknown 4.890  0.0090   32.12  0.0142   0.06
## 5 5 Standard 1 Calibration Standard 4.893  0.4337    7.79  0.6859   3.45
## 6 6 Standard 2 Calibration Standard 4.900  1.0909   10.03  1.7251   6.45
##      IC.S04.5
## 1      BMB
## 2      BMB
## 3      BMB
## 4      BMB
## 5      Rd
## 6      M

## Only keep the columns that we need
Sdat <- Sdat[,c(2,5)] # dont need this here
head(Sdat)

##          X.1 IC.S04.1
## 1  Lab Blank  0.0084
## 2  Lab Blank  0.0091
## 3  Lab Blank  0.0094
## 4  Lab Blank  0.0090
## 5 Standard 1  0.4337
## 6 Standard 2  1.0909

## Name the columns correctly
colnames(Sdat) <- c( "Sample_ID", "S04_ppm")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID S04_ppm
## 1  Lab Blank   0.0084
## 2  Lab Blank   0.0091
## 3  Lab Blank   0.0094
## 4  Lab Blank   0.0090
## 5 Standard 1   0.4337
## 6 Standard 2   1.0909
```

#Chloride data

Read in raw data file from Dionex - copied and saved as a txt

```
Cldat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202305_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X          X.1          X.2  IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1  Lab Blank Unknown 3.890  0.0086   22.84  0.0133   0.05
## 2 2  Lab Blank Unknown 3.857  0.0093   23.54  0.0145   0.05
## 3 3  Lab Blank Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4  Lab Blank Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard 3.900  5.1547   90.74  7.9918  33.37
## 6 6 Standard 2 Calibration Standard 3.860  9.9035  89.27  15.3544  65.84
##      IC.Cl.5
## 1      MB
## 2      MB
## 3     n.a.
## 4     n.a.
## 5      M
## 6      M
```

Only keep the columns that we need

```
Cldat <- Cldat[,c(2,5)]
head(Cldat)
```

```
##          X.1  IC.Cl.1
## 1  Lab Blank 0.0086
## 2  Lab Blank 0.0093
## 3  Lab Blank    n.a.
## 4  Lab Blank    n.a.
## 5 Standard 1 5.1547
## 6 Standard 2 9.9035
```

Name the columns correctly

```
colnames(Cldat) <- c( "Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```

Cldat <- as.data.frame(Cldat)
head(Cldat)

##      Sample_ID Cl_ppm
## 1    Lab Blank  0.0086
## 2    Lab Blank  0.0093
## 3    Lab Blank     NA
## 4    Lab Blank     NA
## 5 Standard 1 5.1547
## 6 Standard 2 9.9035

## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)

##      Sample_ID S04_ppm   Cl_ppm
## 1           911.5208 8679.8805
## 2           911.5208 2977.5413
## 3           911.5208     0.0086
## 4           911.5208 2622.8924
## 5           911.5208       NA
## 6          265.1608 8679.8805

## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID=="") , ]
head(all_dat)

##                  Sample_ID S04_ppm   Cl_ppm
## 26      1_GCW_202305_UP_lyxA_20cm  6.9968 66.1749
## 27      10_GCW_202305_TR_LysB_45cm 10.4414 75.9797
## 28    10_GCW_202305_TR_LysB_45cm_dup 10.4558 75.9221
## 29    100_GWI_202305_WC_LysB_20cm 64.0087 276.9108
## 30  100_GWI_202305_WC_LysB_20cm_dup 33.9433 300.0573
## 31    101_GWI_202305_WC_LysB_45cm 25.5771 260.4622

all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)

##                  Sample_ID S04_ppm   Cl_ppm
## 26      1_GCW_202305_UP_lyxA_20cm  6.9968 66.1749
## 27      10_GCW_202305_TR_LysB_45cm 10.4414 75.9797
## 28    10_GCW_202305_TR_LysB_45cm_dup 10.4558 75.9221
## 29    100_GWI_202305_WC_LysB_20cm 64.0087 276.9108
## 30  100_GWI_202305_WC_LysB_20cm_dup 33.9433 300.0573
## 31    101_GWI_202305_WC_LysB_45cm 25.5771 260.4622

```

Pull out standards - could do some checks here if we want

```

stds <- all_dat[grep("Standard", all_dat$Sample_ID),]
head(stds)

##      Sample_ID S04_ppm Cl_ppm
## 459 Standard 1  0.4622 5.2745
## 460 Standard 1  0.4622 5.1604
## 461 Standard 1  0.4622 5.1547
## 462 Standard 1  0.4622 5.1937
## 463 Standard 1  0.4622 5.2025
## 464 Standard 1  0.4388 5.2745

stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(S04_ppm), sd = sd(S04_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv <2, 'YES', 'NO, rerun')
head(stds_chk_S)

## # A tibble: 5 x 5
##   Sample_ID     mean       sd     cv flag
##   <fct>     <dbl>    <dbl>  <dbl> <chr>
## 1 Standard 1  0.444  0.00987  2.22 NO, rerun
## 2 Standard 2  1.03   0.150    14.5  NO, rerun
## 3 Standard 3  2.06   0.0194   0.940 YES
## 4 Standard 4  9.97   0.165    1.66  YES
## 5 Standard 5 20.0    0.219    1.10  YES

stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv <2, 'YES', 'NO, rerun')
head(stds_chk_Cl)

## # A tibble: 5 x 5
##   Sample_ID     mean       sd     cv flag
##   <fct>     <dbl>    <dbl>  <dbl> <chr>
## 1 Standard 1  5.20   0.0437  0.841 YES
## 2 Standard 2 10.1    0.133   1.32  YES
## 3 Standard 3 20.1    0.188   0.934 YES
## 4 Standard 4 102.    0.860   0.845 YES
## 5 Standard 5 202.    1.58    0.785 YES

```

Calculate mmol/L concentrations

```

#remove standards from sample dataframe
sampledat <- all_dat[!grep("Standard", all_dat$Sample_ID),]
head(sampledat)

```

```

##                                     Sample_ID S04_ppm   Cl_ppm
## 26      1_GCW_202305_UP_lyxA_20cm  6.9968  66.1749
## 27      10_GCW_202305_TR_LysB_45cm 10.4414  75.9797
## 28     10_GCW_202305_TR_LysB_45cm_dup 10.4558  75.9221
## 29      100_GWI_202305_WC_LysB_20cm 64.0087 276.9108
## 30    100_GWI_202305_WC_LysB_20cm_dup 33.9433 300.0573
## 31      101_GWI_202305_WC_LysB_45cm 25.5771 260.4622

# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06       #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.807 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)

##                                     Sample_ID S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 26      1_GCW_202305_UP_lyxA_20cm  6.9968  66.1749  0.2182408 1.866711
## 27      10_GCW_202305_TR_LysB_45cm 10.4414  75.9797  0.3256831 2.143292
## 28     10_GCW_202305_TR_LysB_45cm_dup 10.4558  75.9221  0.3261323 2.141667
## 29      100_GWI_202305_WC_LysB_20cm 64.0087 276.9108  1.9965284 7.811306
## 30    100_GWI_202305_WC_LysB_20cm_dup 33.9433 300.0573  1.0587430 8.464240
## 31      101_GWI_202305_WC_LysB_45cm 25.5771 260.4622  0.7977885 7.347312
##      salinity
## 26 0.1196040
## 27 0.1373213
## 28 0.1372172
## 29 0.5004038
## 30 0.5422295
## 31 0.4706812

```

Pull out dups and check with percent difference

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##                                     Sample_ID S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 26      1_GCW_202305_UP_lyxA_20cm  6.9968  66.1749  0.2182408 1.866711
## 27      10_GCW_202305_TR_LysB_45cm 10.4414  75.9797  0.3256831 2.143292

```

```

## 28 10_GCW_202305_TR_LysB_45cm_dup 10.4558 75.9221 0.3261323 2.141667
## 29      100_GWI_202305_WC_LysB_20cm 64.0087 276.9108 1.9965284 7.811306
## 30 100_GWI_202305_WC_LysB_20cm_dup 33.9433 300.0573 1.0587430 8.464240
## 31      101_GWI_202305_WC_LysB_45cm 25.5771 260.4622 0.7977885 7.347312
##      salinity
## 26 0.1196040
## 27 0.1373213
## 28 0.1372172
## 29 0.5004038
## 30 0.5422295
## 31 0.4706812

```

```

#pull out any rows that have "dup" in the SampleID column
dups <- sampledat %>%
  filter(str_detect(Sample_ID, "dup"))      #have to change this to match data
head(dups)

```

```

##                               Sample_ID  SO4_ppm    Cl_ppm    SO4_mM    Cl_mM
## 1      10_GCW_202305_TR_LysB_45cm_dup 10.4558 75.9221 0.3261323 2.141667
## 2      100_GWI_202305_WC_LysB_20cm_dup 33.9433 300.0573 1.0587430 8.464240
## 3      110_SWH_202305_UP_LysB_10cm_dup 147.8545 2731.8051 4.6118060 77.060793
## 4     120_SWH_202305_UPCON_LysB_20cm_dup 133.1384 1467.6991 4.1527885 41.401949
## 5     130_SWH_202305_TR_LysB_45cm_dup 462.7866 3997.2233 14.4350156 112.756652
## 6     140_SWH_202305_WC_LysC_10cm_dup 409.7469 4046.0590 12.7806269 114.134245
##      salinity
## 1 0.1372172
## 2 0.5422295
## 3 4.9363978
## 4 2.6521583
## 5 7.2230085
## 6 7.3112546

```

```

#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledat %>%
  filter(!str_detect(Sample_ID, "dup")) %>%
  filter(!str_detect(Sample_ID, "spk"))
head(sampledat2)

```

```

##                               Sample_ID  SO4_ppm    Cl_ppm    SO4_mM    Cl_mM  salinity
## 1      1_GCW_202305_UP_lyxA_20cm 6.9968 66.1749 0.2182408 1.866711 0.1196040
## 2      10_GCW_202305_TR_LysB_45cm 10.4414 75.9797 0.3256831 2.143292 0.1373213
## 3     100_GWI_202305_WC_LysB_20cm 64.0087 276.9108 1.9965284 7.811306 0.5004038
## 4     101_GWI_202305_WC_LysB_45cm 25.5771 260.4622 0.7977885 7.347312 0.4706812
## 5     102_GWI_202305_WC_LysC_10cm 34.8612 272.9539 1.0873737 7.699687 0.4932537
## 6     103_GWI_202305_WC_LysC_20cm 39.7708 302.7670 1.2405115 8.540677 0.5471260

```

```

#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[, -c(2,3)]
colnames(dups) <- c('Sample_ID', 'SO4_mM_dup', 'Cl_mM_dup', "salinity_dup")
head(dups)

```

```

##                               Sample_ID S04_mM_dup  Cl_mM_dup salinity_dup
## 1      10_GCW_202305_TR_LysB_45cm  0.3261323  2.141667   0.1372172
## 2      100_GWI_202305_WC_LysB_20cm  1.0587430  8.464240   0.5422295
## 3      110_SWH_202305_UP_LysB_10cm  4.6118060  77.060793  4.9363978
## 4     120_SWH_202305_UPCON_LysB_20cm  4.1527885 41.401949  2.6521583
## 5      130_SWH_202305_TR_LysB_45cm 14.4350156 112.756652  7.2230085
## 6      140_SWH_202305_WC_LysC_10cm 12.7806269 114.134245  7.3112546

```

```

#put it back together with the old data set and look for duplicates
QAdups <- merge(sampledat2, dups)
head(QAdups)

```

```

##                               Sample_ID S04_ppm    Cl_ppm    S04_mM    Cl_mM
## 1      10_GCW_202305_TR_LysB_45cm 10.4414    75.9797  0.3256831  2.143292
## 2      100_GWI_202305_WC_LysB_20cm 64.0087   276.9108 1.9965284  7.811306
## 3      110_SWH_202305_UP_LysB_10cm 215.0081  2446.1470 6.7064286 69.002736
## 4     120_SWH_202305_UPCON_LysB_20cm 162.4131 2237.4438 5.0659108 63.115481
## 5      130_SWH_202305_TR_LysB_45cm 271.3318 3655.9974 8.4632502 103.131097
## 6      140_SWH_202305_WC_LysC_10cm  92.1882 4023.6572 2.8754897 113.502319
##   salinity S04_mM_dup  Cl_mM_dup salinity_dup
## 1 0.1373213  0.3261323  2.141667   0.1372172
## 2 0.5004038  1.0587430  8.464240   0.5422295
## 3 4.4202136  4.6118060  77.060793  4.9363978
## 4 4.0430869  4.1527885 41.401949  2.6521583
## 5 6.6064133 14.4350156 112.756652  7.2230085
## 6 7.2707746 12.7806269 114.134245  7.3112546

```

```

QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*100
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')

```

```

QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')

```

```

head(QAdups)

```

```

##                               Sample_ID S04_ppm    Cl_ppm    S04_mM    Cl_mM
## 1      10_GCW_202305_TR_LysB_45cm 10.4414    75.9797  0.3256831  2.143292
## 2      100_GWI_202305_WC_LysB_20cm 64.0087   276.9108 1.9965284  7.811306
## 3      110_SWH_202305_UP_LysB_10cm 215.0081  2446.1470 6.7064286 69.002736
## 4     120_SWH_202305_UPCON_LysB_20cm 162.4131 2237.4438 5.0659108 63.115481
## 5      130_SWH_202305_TR_LysB_45cm 271.3318 3655.9974 8.4632502 103.131097
## 6      140_SWH_202305_WC_LysC_10cm  92.1882 4023.6572 2.8754897 113.502319
##   salinity S04_mM_dup  Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 0.1373213  0.3261323  2.141667   0.1372172  0.1378175      YES
## 2 0.5004038  1.0587430  8.464240   0.5422295 61.3880268      NO, rerun
## 3 4.4202136  4.6118060  77.060793  4.9363978 37.0132386      NO, rerun
## 4 4.0430869  4.1527885 41.401949  2.6521583 19.8102192      NO, rerun
## 5 6.6064133 14.4350156 112.756652  7.2230085 52.1591068      NO, rerun
## 6 7.2707746 12.7806269 114.134245  7.3112546 126.5337690      NO, rerun
##   Cl_dups_chk Cl_dups_flag
## 1 0.07583847      YES
## 2 8.02349385      YES
## 3 11.03363239    NO, rerun

```

```

## 4 41.55006815      NO, rerun
## 5 8.91718437       YES
## 6 0.55520664       YES


#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between SO4 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call `lifecycle::last_lifecycle_warnings()` to see where this warning was
## generated.

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

	Flag	S_no_rows	Flag	Cl_no_rows	Total	S_Percent	Cl_Percent
## 1	NO, rerun	6	NO, rerun	2	14	42.85714	14.28571
## 2	YES	8	YES	12	14	57.14286	85.71429

Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2, 1, sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv < 11, 'YES', 'NO, rerun')

head(QAdups)

```

	Sample_ID	S04_ppm	Cl_ppm	S04_mM	Cl_mM
## 1	10_GCW_202305_TR_LysB_45cm	10.4414	75.9797	0.3256831	2.143292
## 2	100_GWI_202305_WC_LysB_20cm	64.0087	276.9108	1.9965284	7.811306
## 3	110_SWH_202305_UP_LysB_10cm	215.0081	2446.1470	6.7064286	69.002736
## 4	120_SWH_202305_UPCON_LysB_20cm	162.4131	2237.4438	5.0659108	63.115481
## 5	130_SWH_202305_TR_LysB_45cm	271.3318	3655.9974	8.4632502	103.131097
## 6	140_SWH_202305_WC_LysC_10cm	92.1882	4023.6572	2.8754897	113.502319
	salinity S04_mM_dup	Cl_mM_dup	salinity_dup	S04_dups_chk	S04_dups_flag
## 1	0.1373213	0.3261323	2.141667	0.1372172	0.1378175 YES
## 2	0.5004038	1.0587430	8.464240	0.5422295	61.3880268 NO, rerun
## 3	4.4202136	4.6118060	77.060793	4.9363978	37.0132386 NO, rerun
## 4	4.0430869	4.1527885	41.401949	2.6521583	19.8102192 NO, rerun
## 5	6.6064133	14.4350156	112.756652	7.2230085	52.1591068 NO, rerun
## 6	7.2707746	12.7806269	114.134245	7.3112546	126.5337690 NO, rerun
	Cl_dups_chk	Cl_dups_flag	S04_dups_cv	S04_dups_cv_flag	
## 1	0.07583847	YES	0.09745169	YES	

```


## 2 8.02349385      YES 43.40789003      NO, rerun
## 3 11.03363239     NO, rerun 26.17231202     NO, rerun
## 4 41.55006815     NO, rerun 14.00794033     NO, rerun
## 5 8.91718437      YES 36.88205809      NO, rerun
## 6 0.55520664      YES 89.47288611      NO, rerun

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

	Flag	S_no_rows	Flag	Cl_no_rows	Total	S_Percent	Cl_Percent
## 1	NO, rerun	6	NO, rerun	2	14	42.85714	14.28571
## 2	YES	8	YES	12	14	57.14286	85.71429

Pull out spikes and check

```

#Show me the data that we have from the calculations
head(sampledat)

```

	Sample_ID	S04_ppm	Cl_ppm	S04_mM	Cl_mM
## 26	1_GCW_202305_UP_lysA_20cm	6.9968	66.1749	0.2182408	1.866711
## 27	10_GCW_202305_TR_LysB_45cm	10.4414	75.9797	0.3256831	2.143292
## 28	10_GCW_202305_TR_LysB_45cm_dup	10.4558	75.9221	0.3261323	2.141667
## 29	100_GWI_202305_WC_LysB_20cm	64.0087	276.9108	1.9965284	7.811306
## 30	100_GWI_202305_WC_LysB_20cm_dup	33.9433	300.0573	1.0587430	8.464240
## 31	101_GWI_202305_WC_LysB_45cm	25.5771	260.4622	0.7977885	7.347312
##	salinity				
## 26	0.1196040				
## 27	0.1373213				
## 28	0.1372172				
## 29	0.5004038				
## 30	0.5422295				
## 31	0.4706812				

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk"))      #have to change this to match data
head(spks)

```

	Sample_ID	S04_ppm	Cl_ppm	S04_mM	Cl_mM
## 1	101_GWI_202305_WC_LysB_45cm_spk	71.1032	289.7112	2.2178166	8.172389
## 2	11_GCW_202305_TR_LysC_10cm_spk	8.4098	78.3821	0.2623144	2.211061
## 3	111_SWH_202305_UP_LysB_20cm_spk	201.8888	2235.9854	6.2972177	63.074341
## 4	121_SWH_202305_UPCON_LysB_45cm_spk	28.7040	473.5673	0.8953213	13.358739

```

## 5     131_SWH_202305_TR_LysC_10cm_spk 318.4301 3501.4783 9.9323175 98.772307
## 6     141_SWH_202305_WC_LysC_20cm_spk 109.9728 3997.5809 3.4302183 112.766739
##   salinity
## 1 0.5235341
## 2 0.1416625
## 3 4.0404516
## 4 0.8557621
## 5 6.3271973
## 6 7.2236547

```

```

#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[,-c(2,3, 5,6)]
colnames(spks) <- c('Sample_ID', 'SO4_mM_spk')
head(spks)

```

```

##                               Sample_ID SO4_mM_spk
## 1     101_GWI_202305_WC_LysB_45cm 2.2178166
## 2     11_GCW_202305_TR_LysC_10cm 0.2623144
## 3     111_SWH_202305_UP_LysB_20cm 6.2972177
## 4    121_SWH_202305_UPCON_LysB_45cm 0.8953213
## 5     131_SWH_202305_TR_LysC_10cm 9.9323175
## 6     141_SWH_202305_WC_LysC_20cm 3.4302183

```

```

#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat, spks)
head(QAspks)

```

```

##                               Sample_ID SO4_ppm   Cl_ppm   SO4_mM   Cl_mM
## 1     101_GWI_202305_WC_LysB_45cm 25.5771 260.4622 0.7977885 7.347312
## 2     11_GCW_202305_TR_LysC_10cm  8.6995  77.5035 0.2713506 2.186276
## 3     111_SWH_202305_UP_LysB_20cm 190.0090 2248.7667 5.9266687 63.434886
## 4    121_SWH_202305_UPCON_LysB_45cm 204.1464 2663.3140 6.3676357 75.128745
## 5     131_SWH_202305_TR_LysC_10cm 353.5336 3529.3181 11.0272489 99.557633
## 6     141_SWH_202305_WC_LysC_20cm 388.4101 3977.8297 12.1150998 112.209583
##   salinity SO4_mM_spk
## 1 0.4706812 2.2178166
## 2 0.1400748 0.2623144
## 3 4.0635474 6.2972177
## 4 4.8126344 0.8953213
## 5 6.3775038 9.9323175
## 6 7.1879643 3.4302183

```

```

#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvol <- 10                # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$SO4_spk_Conc <- (spkconc)*spkvol          # mmoles of SO4
head(QAspks)

```

```

##                               Sample_ID SO4_ppm   Cl_ppm   SO4_mM   Cl_mM

```

```

## 1    101_GWI_202305_WC_LysB_45cm 25.5771 260.4622 0.7977885 7.347312
## 2    11_GCW_202305_TR_LysC_10cm 8.6995 77.5035 0.2713506 2.186276
## 3    111_SWH_202305_UP_LysB_20cm 190.0090 2248.7667 5.9266687 63.434886
## 4 121_SWH_202305_UPCON_LysB_45cm 204.1464 2663.3140 6.3676357 75.128745
## 5    131_SWH_202305_TR_LysC_10cm 353.5336 3529.3181 11.0272489 99.557633
## 6    141_SWH_202305_WC_LysC_20cm 388.4101 3977.8297 12.1150998 112.209583
##      salinity S04_mM_spk S04_spk_Conc
## 1 0.4706812 2.2178166 7.797879e-05
## 2 0.1400748 0.2623144 7.797879e-05
## 3 4.0635474 6.2972177 7.797879e-05
## 4 4.8126344 0.8953213 7.797879e-05
## 5 6.3775038 9.9323175 7.797879e-05
## 6 7.1879643 3.4302183 7.797879e-05

#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 50, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 1501, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

```

```

##                               Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 1    101_GWI_202305_WC_LysB_45cm 25.5771 260.4622 0.7977885 7.347312
## 2    11_GCW_202305_TR_LysC_10cm 8.6995 77.5035 0.2713506 2.186276
## 3    111_SWH_202305_UP_LysB_20cm 190.0090 2248.7667 5.9266687 63.434886

```

```

## 4 121_SWH_202305_UPCON_LysB_45cm 204.1464 2663.3140 6.3676357 75.128745
## 5 131_SWH_202305_TR_LysC_10cm 353.5336 3529.3181 11.0272489 99.557633
## 6 141_SWH_202305_WC_LysC_20cm 388.4101 3977.8297 12.1150998 112.209583
## salinity S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 0.4706812 2.2178166 7.797879e-05 1 1e-06
## 2 0.1400748 0.2623144 7.797879e-05 1 1e-06
## 3 4.0635474 6.2972177 7.797879e-05 1 1e-06
## 4 4.8126344 0.8953213 7.797879e-05 1 1e-06
## 5 6.3775038 9.9323175 7.797879e-05 1 1e-06
## 6 7.1879643 3.4302183 7.797879e-05 1 1e-06

#gives us the total S04 in the sample in mmoles
QAspks$S04_Total_unspkd <- (QAspks$S04_mM/QAspks$Dilution)*(QAspks$SampleVol)

##total S04 in spiked sample in mmoles
QAspks$S04_Total_spkd <- (QAspks$S04_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)

QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES' , 'NO', rerun)


head(QAspks)

## Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 101_GWI_202305_WC_LysB_45cm 25.5771 260.4622 0.7977885 7.347312
## 2 11_GCW_202305_TR_LysC_10cm 8.6995 77.5035 0.2713506 2.186276
## 3 111_SWH_202305_UP_LysB_20cm 190.0090 2248.7667 5.9266687 63.434886
## 4 121_SWH_202305_UPCON_LysB_45cm 204.1464 2663.3140 6.3676357 75.128745
## 5 131_SWH_202305_TR_LysC_10cm 353.5336 3529.3181 11.0272489 99.557633
## 6 141_SWH_202305_WC_LysC_20cm 388.4101 3977.8297 12.1150998 112.209583
## salinity S04_mM_spk S04_spk_Conc Dilution SampleVol S04_Total_unspkd
## 1 0.4706812 2.2178166 7.797879e-05 1 1e-06 7.977885e-07
## 2 0.1400748 0.2623144 7.797879e-05 1 1e-06 2.713506e-07
## 3 4.0635474 6.2972177 7.797879e-05 1 1e-06 5.926669e-06
## 4 4.8126344 0.8953213 7.797879e-05 1 1e-06 6.367636e-06
## 5 6.3775038 9.9323175 7.797879e-05 1 1e-06 1.102725e-05
## 6 7.1879643 3.4302183 7.797879e-05 1 1e-06 1.211510e-05
## S04_Total_spkd S04_expctd_spkd spk_recovery S04_spks_flag
## 1 2.439598e-05 7.877658e-05 30.96857 NO, rerun
## 2 2.885459e-06 7.825014e-05 3.68748 NO, rerun
## 3 6.926939e-05 8.390546e-05 82.55648 YES
## 4 9.848534e-06 8.434643e-05 11.67629 NO, rerun
## 5 1.092555e-04 8.900604e-05 122.75065 NO, rerun
## 6 3.773240e-05 9.009389e-05 41.88120 NO, rerun

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)

```

spksbar



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(SO4_spks_flag) %>%
  summarise(no_rows = length(SO4_spks_flag))
Perc_spks$Total <- length(QAspks$SO4_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 2 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>          <int>  <int>   <dbl>
## 1 NO, rerun      13     14    92.9
## 2 YES            1     14     7.14
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names  
#pull the sample ID and separate it by the underscores  
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_', fixed=TRUE)))
```

```
## Warning in rbind(c("1", "GCW", "202305", "UP", "lysA", "20cm"), c("10", :  
## number of columns of result is not a multiple of vector length (arg 7)
```

```
colnames(IDs) <- c("Analysis_No" , "Date", "Site", "Zone", "Replicate", "Depth")  
head(IDs)
```

```
## Analysis_No Date Site Zone Replicate Depth  
## 1 1 GCW 202305 UP lysA 20cm  
## 2 10 GCW 202305 TR LysB 45cm  
## 3 100 GWI 202305 WC LysB 20cm  
## 4 101 GWI 202305 WC LysB 45cm  
## 5 102 GWI 202305 WC LysC 10cm  
## 6 103 GWI 202305 WC LysC 20cm
```

```
#rejoin them to the dataframe  
alldat <- cbind(IDs, sampledat2)  
head(alldat)
```

```
## Analysis_No Date Site Zone Replicate Depth Sample_ID  
## 1 1 GCW 202305 UP lysA 20cm 1_GCW_202305_UP_lysA_20cm  
## 2 10 GCW 202305 TR LysB 45cm 10_GCW_202305_TR_LysB_45cm  
## 3 100 GWI 202305 WC LysB 20cm 100_GWI_202305_WC_LysB_20cm  
## 4 101 GWI 202305 WC LysB 45cm 101_GWI_202305_WC_LysB_45cm  
## 5 102 GWI 202305 WC LysC 10cm 102_GWI_202305_WC_LysC_10cm  
## 6 103 GWI 202305 WC LysC 20cm 103_GWI_202305_WC_LysC_20cm  
## S04_ppm Cl_ppm S04_mM Cl_mM salinity  
## 1 6.9968 66.1749 0.2182408 1.866711 0.1196040  
## 2 10.4414 75.9797 0.3256831 2.143292 0.1373213  
## 3 64.0087 276.9108 1.9965284 7.811306 0.5004038  
## 4 25.5771 260.4622 0.7977885 7.347312 0.4706812  
## 5 34.8612 272.9539 1.0873737 7.699687 0.4932537  
## 6 39.7708 302.7670 1.2405115 8.540677 0.5471260
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_S04_202305.csv") #Change f
```

END