

# AQ300Test Parameters



Date & Time: 2022-07-19 13:50:25

Test Number: 11

Long Test Name: Vanadium NOx

Short Test Name: V-NOx

Units: mg N/L

Decimal Places: 4

Test Type: 1. End Point

## Main Parameters

|                 |     |                 |   |                   |            |
|-----------------|-----|-----------------|---|-------------------|------------|
| Sample Volume:  | 180 | Water Volume:   | 0 | Number of Mixes:  | 1          |
| Cuvette primes: | 2   | Cuvette Washes: | 2 | Baseline on Wash: | YES        |
| Cadmium Volume: | 0   | Reduction Time: | 0 | Diluent Location: | Reagent 18 |

## Math Parameters

|                     |        |                       |        |                         |        |
|---------------------|--------|-----------------------|--------|-------------------------|--------|
| Reaction Time:      | 1825   | Wavelength:           | 520nm  | Polynomial Order:       | 2      |
| Repeat Delta Check: | 0.250  | Auto Dilution Factor: | 5.000  | Manual Dilution Factor: | 0.000  |
| Linearity Low:      | 0.1766 | Linearity High:       | 0.9483 | Latest Slope:           | 1.0000 |
| Alert Low:          | 0.000  | Alert High:           | 0.000  | Latest Intercept:       | 0.0000 |

## Reagent Parameters

Number of Reagents: 2

| Reagent | Name                   | Volume | Delay | Replaced in Blank |
|---------|------------------------|--------|-------|-------------------|
| 1       | DI                     | 10     | 0     | NO                |
| 2       | Vanadium Color Reagent | 310    | 0     | NO                |
| Blank   | Water Bottle           | 10     |       |                   |

## Standards

Auto Std. Number: 7

Auto Std. Concentration : 1.0000

Exclude the Blank : NO

Auto Std. Position: Cup 1

Standard 1 (Blank) Position : Reagent 18

Correlation limit : 0.9950

| Standard | Value  | Standard | Value  | Auto Standard | Percentage | Value  |
|----------|--------|----------|--------|---------------|------------|--------|
| 1        | 0.0000 | 11       | 0.0000 | S90           | 2.2222     | 0.0222 |
| 2        | 0.0000 | 12       | 0.0000 | S91           | 5.0000     | 0.0500 |
| 3        | 0.0000 | 13       | 0.0000 | S92           | 10.0000    | 0.1000 |
| 4        | 0.0000 | 14       | 0.0000 | S93           | 25.0000    | 0.2500 |
| 5        | 0.0000 | 15       | 0.0000 | S94           | 50.0000    | 0.5000 |
| 6        | 0.0000 | 16       | 0.0000 | S95           | 75.0000    | 0.7500 |
| 7        | 0.0000 | 17       | 0.0000 | S96           | 100.0000   | 1.0000 |

## Controls

| Control | Low Value | High Value | Location | Frequency | Start Before |
|---------|-----------|------------|----------|-----------|--------------|
|---------|-----------|------------|----------|-----------|--------------|

## QC PRO

| Control | Low Value | High Value | Location     | Frequency  | Start Before |
|---------|-----------|------------|--------------|------------|--------------|
| CCV     | 0.4000    | 0.6000     | Reagent 10   | 10         | 1            |
| CCB     | -0.0400   | 0.0400     | Reagent 18   | 10         | 1            |
| ICV     | 0.0000    | 0.0000     | As scheduled | 0          | Off          |
| LCS     | 0.0000    | 0.0000     | As scheduled | 0          | Off          |
| ICB     | 0.0000    | 0.0000     | As scheduled | 0          | Off          |
| MB      | 0.0000    | 0.0000     | As scheduled | 0          | Off          |
| Control | % (+or-)  | Stock Conc | Spike Added  | Location   |              |
| DUP     | 10.0000   |            |              |            |              |
| SPKA    | 50.0000   | 5.0000     | 0.5000       | Reagent 11 |              |
| SPKM    | 0.0000    |            | 0.0000       |            |              |
| SPKD    | 0.0000    |            |              |            |              |