

# Dionex\_COMPASS\_June2023

Stephanie J. Wilson

2023-06-23

## Daily Set up

Read in and Format the raw data - change wd & file names

```
# SULFATE DATA:
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202309_S04.txt", sep='\t' , header=T, skip=3)
head(Sdat)
```

```
##      X          X.1                X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard  n.a.    n.a.    n.a.    n.a.    n.a.
## 6 6 Standard 2 Calibration Standard 4.753  0.9994    6.87    1.1465    7.60
```

```
## Only keep the columns that we need
```

```
Sdat <- Sdat[,c(2,5)] # dont need this here
head(Sdat)
```

```
##          X.1 IC.S04.1
## 1 Lab Blank    n.a.
## 2 Lab Blank    n.a.
## 3 Lab Blank    n.a.
## 4 Lab Blank    n.a.
## 5 Standard 1    n.a.
## 6 Standard 2 0.9994
```

```
## Name the columns correctly
```

```
colnames(Sdat) <- c( "Sample_ID", "S04_ppm")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID S04_ppm
## 1 Lab Blank      NA
## 2 Lab Blank      NA
## 3 Lab Blank      NA
## 4 Lab Blank      NA
## 5 Standard 1      NA
## 6 Standard 2  0.9994
```

```
#Chloride data
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Cldat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202309_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X      X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank      Unknown 3.700 0.0065 100.00 0.0102 0.07
## 2 2 Lab Blank      Unknown  n.a.   n.a.   n.a.   n.a.   n.a.
## 3 3 Lab Blank      Unknown  n.a.   n.a.   n.a.   n.a.   n.a.
## 4 4 Lab Blank      Unknown  n.a.   n.a.   n.a.   n.a.   n.a.
## 5 5 Standard 1 Calibration Standard 3.740 4.9411 93.38 7.7044 64.09
## 6 6 Standard 2 Calibration Standard 3.767 9.9736 93.13 15.5513 127.45
```

```
## Only keep the columns that we need
```

```
Cldat <- Cldat[ ,c(2,5)]
head(Cldat)
```

```
##      X.1 IC.Cl.1
## 1 Lab Blank 0.0065
## 2 Lab Blank  n.a.
## 3 Lab Blank  n.a.
## 4 Lab Blank  n.a.
## 5 Standard 1 4.9411
## 6 Standard 2 9.9736
```

```
## Name the columns correctly
```

```
colnames(Cldat) <- c( "Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##      Sample_ID Cl_ppm
## 1 Lab Blank 0.0065
## 2 Lab Blank  NA
## 3 Lab Blank  NA
## 4 Lab Blank  NA
## 5 Standard 1 4.9411
## 6 Standard 2 9.9736
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##      Sample_ID  S04_ppm  Cl_ppm
## 1          1205.7288 9437.8836
## 2          1205.7288 2765.0461
## 3          1205.7288   0.0065
## 4          1205.7288 3281.0669
## 5          1205.7288    NA
## 6          279.7333 9437.8836
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##              Sample_ID  S04_ppm  Cl_ppm
## 26      1_GCW_202309_UP_LysA_20cm  7.4516  3.3628
## 27     10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764
## 28 10_GCW_202309_WC_LysA_10cm_dup 596.0266 7030.0612
## 29     11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026
## 30 11_GCW_202309_WC_LysA_20cm_spk 490.6059 6428.3445
## 31     12_GCW_202309_WC_LysA_45cm 228.8569 5506.8946
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
```

```
## If there is an n.a. in the dataframe that is because there was no peak, which
# would mean there was no Cl or SO4 and we want to know that so make all n.a.'s into zeros
all_dat[is.na(all_dat)] <- 0
head(all_dat)
```

```
##              Sample_ID  S04_ppm  Cl_ppm
## 26      1_GCW_202309_UP_LysA_20cm  7.4516  3.3628
## 27     10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764
## 28 10_GCW_202309_WC_LysA_10cm_dup 596.0266 7030.0612
## 29     11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026
## 30 11_GCW_202309_WC_LysA_20cm_spk 490.6059 6428.3445
## 31     12_GCW_202309_WC_LysA_45cm 228.8569 5506.8946
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grepl("Standard", all_dat$Sample_ID),]
#stds <- stds[-c(17),] #this is if you need to remove one for any reason
head(stds)
```

```
##      Sample_ID S04_ppm Cl_ppm
## 340 Standard 1  0.4898 5.0593
## 341 Standard 1  0.4898 5.2609
## 342 Standard 1  0.4898 4.9411
```

```
## 343 Standard 1 0.4898 5.1478
## 344 Standard 1 0.5051 5.0593
## 345 Standard 1 0.5051 5.2609
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(SO4_ppm), sd = sd(SO4_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID   mean     sd    cv flag
##   <fct>       <dbl>  <dbl> <dbl> <chr>
## 1 Standard 1 0.373 0.222 59.6 NO, rerun
## 2 Standard 2 1.03 0.0276 2.67 YES
## 3 Standard 3 2.04 0.0371 1.82 YES
## 4 Standard 4 7.11 5.33 75.0 NO, rerun
## 5 Standard 5 14.0 10.5 75.0 NO, rerun
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID   mean     sd    cv flag
##   <fct>       <dbl>  <dbl> <dbl> <chr>
## 1 Standard 1 5.10 0.121 2.37 YES
## 2 Standard 2 10.3 0.288 2.79 YES
## 3 Standard 3 20.0 0.432 2.16 YES
## 4 Standard 4 102. 2.18 2.13 YES
## 5 Standard 5 206. 4.37 2.12 YES
```

## Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##           Sample_ID SO4_ppm  Cl_ppm
## 26 1_GCW_202309_UP_LysA_20cm 7.4516 3.3628
## 27 10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764
## 28 10_GCW_202309_WC_LysA_10cm_dup 596.0266 7030.0612
## 29 11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026
## 30 11_GCW_202309_WC_LysA_20cm_spk 490.6059 6428.3445
## 31 12_GCW_202309_WC_LysA_45cm 228.8569 5506.8946
```

```

# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06      #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$SO4_mM <- (sampledat$SO4_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)

```

```

##           Sample_ID  SO4_ppm  Cl_ppm  SO4_mM  Cl_mM
## 26      1_GCW_202309_UP_LysA_20cm  7.4516  3.3628  0.2324267  0.09486037
## 27     10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764 18.5110480 196.86252186
## 28 10_GCW_202309_WC_LysA_10cm_dup 596.0266 7030.0612 18.5909732 198.30920169
## 29     11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026  9.4869994 185.31742172
## 30 11_GCW_202309_WC_LysA_20cm_spk 490.6059 6428.3445 15.3027417 181.33552891
## 31     12_GCW_202309_WC_LysA_45cm 228.8569 5506.8946  7.1383936 155.34258392
##      salinity
## 26  0.00610258
## 27 12.61067495
## 28 12.70334659
## 29 11.87111720
## 30 11.61604451
## 31  9.95098454

```

Pull out dups and check with percent difference

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##           Sample_ID  SO4_ppm  Cl_ppm  SO4_mM  Cl_mM
## 26      1_GCW_202309_UP_LysA_20cm  7.4516  3.3628  0.2324267  0.09486037
## 27     10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764 18.5110480 196.86252186
## 28 10_GCW_202309_WC_LysA_10cm_dup 596.0266 7030.0612 18.5909732 198.30920169
## 29     11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026  9.4869994 185.31742172
## 30 11_GCW_202309_WC_LysA_20cm_spk 490.6059 6428.3445 15.3027417 181.33552891
## 31     12_GCW_202309_WC_LysA_45cm 228.8569 5506.8946  7.1383936 155.34258392
##      salinity
## 26  0.00610258
## 27 12.61067495
## 28 12.70334659

```

```
## 29 11.87111720
## 30 11.61604451
## 31 9.95098454
```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledats %>%
  filter(str_detect(Sample_ID, "dup")) #have to change this to match data
head(dups)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 10_GCW_202309_WC_LysA_10cm_dup 596.0266 7030.0612 18.5909732 1.983092e+02
## 2 20_GCW_202309_SW_B_dup 770.0896 5898.2462 24.0202620 1.663821e+02
## 3 30_MSM_202309_TR_LysA_45cm_dup 251.9186 1988.5976 7.8577230 5.609584e+01
## 4 40_MSM_202309_WC_LysB_10cm_dup 634.3614 8701.6746 19.7866937 2.454633e+02
## 5 50_SWH_202309_UPCON_LysA_20cm_dup 0.0000 0.0679 0.0000000 1.915374e-03
## 6 60_SWH_202309_UP_LysA_45cm_dup 5.4856 214.0801 0.1711042 6.038931e+00
##           salinity
## 1 1.270335e+01
## 2 1.065816e+01
## 3 3.593422e+00
## 4 1.572395e+01
## 5 1.486953e-04
## 6 3.868687e-01
```

```
#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledats %>%
  filter(!str_detect(Sample_ID, "dup")) %>%
  filter(!str_detect(Sample_ID, "spk"))
head(sampledat2)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 1_GCW_202309_UP_LysA_20cm 7.4516 3.3628 0.2324267 0.09486037
## 2 10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764 18.5110480 196.86252186
## 3 11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026 9.4869994 185.31742172
## 4 12_GCW_202309_WC_LysA_45cm 228.8569 5506.8946 7.1383936 155.34258392
## 5 13_GCW_202309_WC_LysB_10cm 293.2979 6840.1267 9.1484061 192.95138787
## 6 14_GCW_202309_WC_LysB_20cm 134.6521 6335.1128 4.2000031 178.70557969
##           salinity
## 1 0.00610258
## 2 12.61067495
## 3 11.87111720
## 4 9.95098454
## 5 12.36013495
## 6 11.44757483
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup", "", as.character(dups$Sample_ID))
dups <- dups[ , -c(2,3)]
colnames(dups) <- c('Sample_ID', 'S04_mM_dup', 'Cl_mM_dup', 'salinity_dup')
head(dups)
```

```
##           Sample_ID S04_mM_dup  Cl_mM_dup salinity_dup
```

```
## 1 10_GCW_202309_WC_LysA_10cm 18.5909732 1.983092e+02 1.270335e+01
## 2 20_GCW_202309_SW_B 24.0202620 1.663821e+02 1.065816e+01
## 3 30_MSM_202309_TR_LysA_45cm 7.8577230 5.609584e+01 3.593422e+00
## 4 40_MSM_202309_WC_LysB_10cm 19.7866937 2.454633e+02 1.572395e+01
## 5 50_SWH_202309_UPCON_LysA_20cm 0.0000000 1.915374e-03 1.486953e-04
## 6 60_SWH_202309_UP_LysA_45cm 0.1711042 6.038931e+00 3.868687e-01
```

*#put it back together with the old data set and look for duplicates*

```
QAdups <- merge(sampledat2, dups)
head(QAdups)
```

```
##           Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764 18.51104803 196.86252186
## 2 20_GCW_202309_SW_B 809.5272 6202.5891 25.25038054 174.96725247
## 3 30_MSM_202309_TR_LysA_45cm 503.5731 3988.6796 15.70720836 112.51564457
## 4 40_MSM_202309_WC_LysB_10cm 685.5478 9424.2428 21.38327511 265.84605924
## 5 50_SWH_202309_UPCON_LysA_20cm 0.1084 0.6823 0.00338116 0.01924683
## 6 60_SWH_202309_UP_LysA_45cm 7.1839 227.6073 0.22407673 6.42051622
## salinity S04_mM_dup Cl_mM_dup salinity_dup
## 1 12.610674955 18.5909732 1.983092e+02 1.270335e+01
## 2 11.208104504 24.0202620 1.663821e+02 1.065816e+01
## 3 7.207570037 7.8577230 5.609584e+01 3.593422e+00
## 4 17.029632740 19.7866937 2.454633e+02 1.572395e+01
## 5 0.001258916 0.0000000 1.915374e-03 1.486953e-04
## 6 0.411312391 0.1711042 6.038931e+00 3.868687e-01
```

```
QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*100
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')
```

```
QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')
```

```
head(QAdups)
```

```
##           Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764 18.51104803 196.86252186
## 2 20_GCW_202309_SW_B 809.5272 6202.5891 25.25038054 174.96725247
## 3 30_MSM_202309_TR_LysA_45cm 503.5731 3988.6796 15.70720836 112.51564457
## 4 40_MSM_202309_WC_LysB_10cm 685.5478 9424.2428 21.38327511 265.84605924
## 5 50_SWH_202309_UPCON_LysA_20cm 0.1084 0.6823 0.00338116 0.01924683
## 6 60_SWH_202309_UP_LysA_45cm 7.1839 227.6073 0.22407673 6.42051622
## salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 12.610674955 18.5909732 1.983092e+02 1.270335e+01 0.4308398 YES
## 2 11.208104504 24.0202620 1.663821e+02 1.065816e+01 4.9933123 YES
## 3 7.207570037 7.8577230 5.609584e+01 3.593422e+00 66.6200568 NO, rerun
## 4 17.029632740 19.7866937 2.454633e+02 1.572395e+01 7.7560487 YES
## 5 0.001258916 0.0000000 1.915374e-03 1.486953e-04 200.0000000 NO, rerun
## 6 0.411312391 0.1711042 6.038931e+00 3.868687e-01 26.8092663 NO, rerun
## Cl_dups_chk Cl_dups_flag
## 1 0.7321778 YES
## 2 5.0301139 YES
## 3 66.9228457 NO, rerun
## 4 7.9727628 YES
```

```
## 5 163.7963210    NO, rerun
## 6   6.1252370      YES
```


*#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this*

```
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)
```

```
ggarrange(Sdupsbar, Cldupsbar,ncol=2, nrow=1)
```





```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##          Flag S_no_rows      Flag Cl_no_rows Total S_Percent Cl_Percent
## 1 NO, rerun          3 NO, rerun          2      7  42.85714  28.57143
## 2      YES          4      YES          5      7  57.14286  71.42857

```

## Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)

```

```

##          Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764 18.51104803 196.86252186
## 2      20_GCW_202309_SW_B 809.5272 6202.5891 25.25038054 174.96725247
## 3 30_MSM_202309_TR_LysA_45cm 503.5731 3988.6796 15.70720836 112.51564457
## 4 40_MSM_202309_WC_LysB_10cm 685.5478 9424.2428 21.38327511 265.84605924
## 5 50_SWH_202309_UPCON_LysA_20cm 0.1084 0.6823 0.00338116 0.01924683
## 6 60_SWH_202309_UP_LysA_45cm 7.1839 227.6073 0.22407673 6.42051622
##          salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 12.610674955 18.5909732 1.983092e+02 1.270335e+01 0.4308398 YES
## 2 11.208104504 24.0202620 1.663821e+02 1.065816e+01 4.9933123 YES
## 3 7.207570037 7.8577230 5.609584e+01 3.593422e+00 66.6200568 NO, rerun
## 4 17.029632740 19.7866937 2.454633e+02 1.572395e+01 7.7560487 YES
## 5 0.001258916 0.0000000 1.915374e-03 1.486953e-04 200.0000000 NO, rerun
## 6 0.411312391 0.1711042 6.038931e+00 3.868687e-01 26.8092663 NO, rerun
##          Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1 0.7321778 YES 0.3046498 YES


```

## 2	5.0301139	YES	3.5308050	YES
## 3	66.9228457	NO, rerun	47.1074940	NO, rerun
## 4	7.9727628	YES	5.4843546	YES
## 5	163.7963210	NO, rerun	141.4213562	NO, rerun
## 6	6.1252370	YES	18.9570140	NO, rerun

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##           Flag S_no_rows           Flag Cl_no_rows Total S_Percent Cl_Percent
## 1 NO, rerun           3 NO, rerun           2      7  42.85714  28.57143
## 2      YES           4      YES           5      7  57.14286  71.42857

```

## Pull out spikes and check

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##           Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
## 26 1_GCW_202309_UP_LysA_20cm 7.4516 3.3628 0.2324267 0.09486037
## 27 10_GCW_202309_WC_LysA_10cm 593.4642 6978.7764 18.5110480 196.86252186
## 28 10_GCW_202309_WC_LysA_10cm_dup 596.0266 7030.0612 18.5909732 198.30920169
## 29 11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026 9.4869994 185.31742172
## 30 11_GCW_202309_WC_LysA_20cm_spk 490.6059 6428.3445 15.3027417 181.33552891
## 31 12_GCW_202309_WC_LysA_45cm 228.8569 5506.8946 7.1383936 155.34258392
##           salinity
## 26 0.00610258
## 27 12.61067495
## 28 12.70334659
## 29 11.87111720
## 30 11.61604451
## 31 9.95098454

```

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)

```

```

##           Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
## 1 11_GCW_202309_WC_LysA_20cm_spk 490.6059 6428.3445 15.3027417 181.335529
## 2 21_GCW_202309_SW_C_spk 958.1205 5884.6478 29.8852308 165.998528
## 3 31_MSM_202309_TR_LysB_10cm_spk 444.0675 3520.5436 13.8511385 99.310116
## 4 41_MSM_202309_WC_LysB_20cm_spk 922.5773 9076.6365 28.7765845 256.040522

```

```
## 5 51_SWH_202309_UPCON_LysA_45cm_spk 138.4937 108.5239 4.3198284 3.061323
## 6 61_SWH_202309_UP_LysB_10cm_spk 5.5928 165.5009 0.1744479 4.668573
## salinity
## 1 11.6160445
## 2 10.6335846
## 3 6.3616483
## 4 16.4015082
## 5 0.1961287
## 6 0.2990861
```

```
#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[ ,-c(2,3, 5,6)]
colnames(spks) <- c('Sample_ID', 'SO4_mM_spk')
head(spks)
```

```
## Sample_ID SO4_mM_spk
## 1 11_GCW_202309_WC_LysA_20cm 15.3027417
## 2 21_GCW_202309_SW_C 29.8852308
## 3 31_MSM_202309_TR_LysB_10cm 13.8511385
## 4 41_MSM_202309_WC_LysB_20cm 28.7765845
## 5 51_SWH_202309_UPCON_LysA_45cm 4.3198284
## 6 61_SWH_202309_UP_LysB_10cm 0.1744479
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat, spks)
head(QAspks)
```

```
## Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
## 1 11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026 9.4869994 185.317422
## 2 21_GCW_202309_SW_C 811.2022 6193.2976 25.3026263 174.705151
## 3 31_MSM_202309_TR_LysB_10cm 675.2835 6882.4336 21.0631160 194.144812
## 4 41_MSM_202309_WC_LysB_20cm 728.9486 9246.8867 22.7370119 260.843066
## 5 51_SWH_202309_UPCON_LysA_45cm 30.3570 117.3651 0.9468808 3.310722
## 6 61_SWH_202309_UP_LysB_10cm 48.4169 502.6916 1.5101965 14.180299
## salinity SO4_mM_spk
## 1 11.8711172 15.3027417
## 2 11.1913148 29.8852308
## 3 12.4365835 13.8511385
## 4 16.7091503 28.7765845
## 5 0.2121047 4.3198284
## 6 0.9083897 0.1744479
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw) # in mM
spkvol <- 10 # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$SO4_spk_Conc <- (spkconc)*spkvol # mmoles of SO4
head(QAspks)
```

```
## Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
```

```
## 1    11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026 9.4869994 185.317422
## 2          21_GCW_202309_SW_C 811.2022 6193.2976 25.3026263 174.705151
## 3    31_MSM_202309_TR_LysB_10cm 675.2835 6882.4336 21.0631160 194.144812
## 4    41_MSM_202309_WC_LysB_20cm 728.9486 9246.8867 22.7370119 260.843066
## 5 51_SWH_202309_UPCON_LysA_45cm 30.3570 117.3651 0.9468808 3.310722
## 6    61_SWH_202309_UP_LysB_10cm 48.4169 502.6916 1.5101965 14.180299
##      salinity S04_mM_spk S04_spk_Conc
## 1 11.8711172 15.3027417 7.797879e-05
## 2 11.1913148 29.8852308 7.797879e-05
## 3 12.4365835 13.8511385 7.797879e-05
## 4 16.7091503 28.7765845 7.797879e-05
## 5 0.2121047 4.3198284 7.797879e-05
## 6 0.9083897 0.1744479 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
```

```
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 50, QAspks$Dilution)
```

```
#Set Sample volumes in uL
```

```
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 1462, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 1501, QAspks$SampleVol)
```

```
#change sample volume to L
```

```
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)
```

```
##           Sample_ID S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 1    11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026 9.4869994 185.317422
## 2          21_GCW_202309_SW_C 811.2022 6193.2976 25.3026263 174.705151
## 3    31_MSM_202309_TR_LysB_10cm 675.2835 6882.4336 21.0631160 194.144812
```

```
## 4 41_MSM_202309_WC_LysB_20cm 728.9486 9246.8867 22.7370119 260.843066
## 5 51_SWH_202309_UPCON_LysA_45cm 30.3570 117.3651 0.9468808 3.310722
## 6 61_SWH_202309_UP_LysB_10cm 48.4169 502.6916 1.5101965 14.180299
## salinity S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 11.8711172 15.3027417 7.797879e-05 1 1e-06
## 2 11.1913148 29.8852308 7.797879e-05 1 1e-06
## 3 12.4365835 13.8511385 7.797879e-05 1 1e-06
## 4 16.7091503 28.7765845 7.797879e-05 1 1e-06
## 5 0.2121047 4.3198284 7.797879e-05 1 1e-06
## 6 0.9083897 0.1744479 7.797879e-05 1 1e-06
```

*#gives us the total SO4 in the sample in mmoles*

```
QAspks$S04_Total_unspkd <- (QAspks$S04_mM/QAspks$Dilution)*(QAspks$SampleVol)
```

*##total SO4 in spiked sample in mmoles*

```
QAspks$S04_Total_spkd <- (QAspks$S04_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)
```

```
QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
```

```
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun')
```


```
head(QAspks)
```

```
## Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 11_GCW_202309_WC_LysA_20cm 304.1532 6569.5026 9.4869994 185.317422
## 2 21_GCW_202309_SW_C 811.2022 6193.2976 25.3026263 174.705151
## 3 31_MSM_202309_TR_LysB_10cm 675.2835 6882.4336 21.0631160 194.144812
## 4 41_MSM_202309_WC_LysB_20cm 728.9486 9246.8867 22.7370119 260.843066
## 5 51_SWH_202309_UPCON_LysA_45cm 30.3570 117.3651 0.9468808 3.310722
## 6 61_SWH_202309_UP_LysB_10cm 48.4169 502.6916 1.5101965 14.180299
## salinity S04_mM_spk S04_spk_Conc Dilution SampleVol S04_Total_unspkd
## 1 11.8711172 15.3027417 7.797879e-05 1 1e-06 9.486999e-06
## 2 11.1913148 29.8852308 7.797879e-05 1 1e-06 2.530263e-05
## 3 12.4365835 13.8511385 7.797879e-05 1 1e-06 2.106312e-05
## 4 16.7091503 28.7765845 7.797879e-05 1 1e-06 2.273701e-05
## 5 0.2121047 4.3198284 7.797879e-05 1 1e-06 9.468808e-07
## 6 0.9083897 0.1744479 7.797879e-05 1 1e-06 1.510197e-06
## S04_Total_spkd S04_expctd_spkd spk_recovery S04_spks_flag
## 1 1.683302e-04 8.746579e-05 192.452570 NO, rerun
## 2 3.287375e-04 1.032814e-04 318.293021 NO, rerun
## 3 1.523625e-04 9.904191e-05 153.836421 NO, rerun
## 4 3.165424e-04 1.007158e-04 314.292717 NO, rerun
## 5 4.751811e-05 7.892567e-05 60.206157 NO, rerun
## 6 1.918927e-06 7.948899e-05 2.414079 NO, rerun
```

*#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on*

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

spksbar



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 2 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 NO, rerun           6     7   85.7
## 2 YES                 1     7   14.3
```

## Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_', fixed=TRUE)))
```

```
## Warning in rbind(c("1", "GCW", "202309", "UP", "LysA", "20cm"), c("10", :
## number of columns of result is not a multiple of vector length (arg 11)
```

```
colnames(IDs) <- c("Analysis_No", "Site", "Date", "Zone", "Replicate", "Depth") #, "RHZ", "RHZ_Rep")
head(IDs)
```

```
##   Analysis_No Site   Date Zone Replicate Depth
## 1           1  GCW 202309  UP      LysA  20cm
## 2          10  GCW 202309  WC      LysA  10cm
## 3          11  GCW 202309  WC      LysA  20cm
## 4          12  GCW 202309  WC      LysA  45cm
## 5          13  GCW 202309  WC      LysB  10cm
## 6          14  GCW 202309  WC      LysB  20cm
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampled2)
head(alldat)
```

```
##   Analysis_No Site   Date Zone Replicate Depth      Sample_ID
## 1           1  GCW 202309  UP      LysA  20cm 1_GCW_202309_UP_LysA_20cm
## 2          10  GCW 202309  WC      LysA  10cm 10_GCW_202309_WC_LysA_10cm
## 3          11  GCW 202309  WC      LysA  20cm 11_GCW_202309_WC_LysA_20cm
## 4          12  GCW 202309  WC      LysA  45cm 12_GCW_202309_WC_LysA_45cm
## 5          13  GCW 202309  WC      LysB  10cm 13_GCW_202309_WC_LysB_10cm
## 6          14  GCW 202309  WC      LysB  20cm 14_GCW_202309_WC_LysB_20cm
##   S04_ppm   Cl_ppm   S04_mM   Cl_mM   salinity
## 1   7.4516   3.3628  0.2324267  0.09486037  0.00610258
## 2 593.4642 6978.7764 18.5110480 196.86252186 12.61067495
## 3 304.1532 6569.5026  9.4869994 185.31742172 11.87111720
## 4 228.8569 5506.8946  7.1383936 155.34258392  9.95098454
## 5 293.2979 6840.1267  9.1484061 192.95138787 12.36013495
## 6 134.6521 6335.1128  4.2000031 178.70557969 11.44757483
```

## Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_S04_202309.csv")
```

```
#C
```

END