

Dionex_COMPASS_May2022

Stephanie J. Wilson

2023-01-23

Daily Set up

Read in and Format the raw data - change wd & file names

```
#setwd("S:/Biogeochemistry/People/Wilson (Steph)/Data/Dionex/Raw Data Files")

# SULFATE DATA:

## Read in raw data file from Dionex - copied and saved as a txt
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202205_S04.txt",sep='\t' , header=T, skip=3)
head(Sdat)

##      X      X.1          X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard 4.133   0.4711   6.47   0.4192   2.16
## 6 6 Standard 2 Calibration Standard 4.130   0.9084   6.22   0.8083   4.14

## Only keep the columns that we need
Sdat <- Sdat[,c(2,5)] # dont need this here
head(Sdat)

##          X.1 IC.S04.1
## 1  Lab Blank    n.a.
## 2  Lab Blank    n.a.
## 3  Lab Blank    n.a.
## 4  Lab Blank    n.a.
## 5 Standard 1  0.4711
## 6 Standard 2  0.9084

## Name the columns correctly
colnames(Sdat) <- c( "Sample_ID", "S04_ppm")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)

## Warning: NAs introduced by coercion
```

```

Sdat <- as.data.frame(Sdat)
head(Sdat)

##      Sample_ID SO4_ppm
## 1    Lab Blank     NA
## 2    Lab Blank     NA
## 3    Lab Blank     NA
## 4    Lab Blank     NA
## 5 Standard 1  0.4711
## 6 Standard 2  0.9084

#Chloride data
## Read in raw data file from Dionex - copied and saved as a txt
Cldat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202205_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)

##      X       X.1          X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard 3.243  5.0098  92.77  6.0066  38.67
## 6 6 Standard 2 Calibration Standard 3.243 10.0990  93.21 12.1084  77.34
##      IC.Cl.5
## 1      n.a.
## 2      n.a.
## 3      n.a.
## 4      n.a.
## 5      BMB
## 6      BMB

## Only keep the columns that we need
Cldat <- Cldat[,c(2,5)]
head(Cldat)

##      X.1 IC.Cl.1
## 1  Lab Blank   n.a.
## 2  Lab Blank   n.a.
## 3  Lab Blank   n.a.
## 4  Lab Blank   n.a.
## 5 Standard 1  5.0098
## 6 Standard 2 10.0990

## Name the columns correctly
colnames(Cldat) <- c( "Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)

## Warning: NAs introduced by coercion

```

```

Cldat <- as.data.frame(Cldat)
head(Cldat)

##      Sample_ID Cl_ppm
## 1    Lab Blank     NA
## 2    Lab Blank     NA
## 3    Lab Blank     NA
## 4    Lab Blank     NA
## 5 Standard 1  5.0098
## 6 Standard 2 10.0990

## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)

##      Sample_ID S04_ppm   Cl_ppm
## 1           1122.213 8909.0812
## 2           1122.213 3574.4197
## 3           1122.213     5.0098
## 4           1122.213 3462.5246
## 5           1122.213       NA
## 6           429.473 8909.0812

## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID=="") , ]
head(all_dat)

##                               Sample_ID S04_ppm   Cl_ppm
## 26      COMPASS_20220502_MSM_SipA_10cm 924.7531 5530.741
## 27      COMPASS_20220502_MSM_SipA_20cm 686.1893 5431.449
## 28      COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549
## 29      COMPASS_20220502_MSM_SipA_45cm_Dup 398.1761 6525.085
## 30 COMPASS_20220502_MSM_SipA_45cm_Spike 606.2450 6339.224
## 31      COMPASS_20220502_MSM_SipB_10cm 747.5236 5169.818

all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)

##                               Sample_ID S04_ppm   Cl_ppm
## 26      COMPASS_20220502_MSM_SipA_10cm 924.7531 5530.741
## 27      COMPASS_20220502_MSM_SipA_20cm 686.1893 5431.449
## 28      COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549
## 29      COMPASS_20220502_MSM_SipA_45cm_Dup 398.1761 6525.085
## 30 COMPASS_20220502_MSM_SipA_45cm_Spike 606.2450 6339.224
## 31      COMPASS_20220502_MSM_SipB_10cm 747.5236 5169.818

```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grep1("Standard", all_dat$Sample_ID),]
head(stds)
```

```
##           Sample_ID S04_ppm Cl_ppm
## 184 Standard 1   0.4979 5.2975
## 185 Standard 1   0.4979 5.0098
## 186 Standard 1   0.4979 5.1659
## 187 Standard 1   0.4979 5.2056
## 188 Standard 1   0.4711 5.2975
## 189 Standard 1   0.4711 5.0098
```

Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grep1("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##           Sample_ID S04_ppm Cl_ppm
## 26      COMPASS_20220502_MSM_SipA_10cm 924.7531 5530.741
## 27      COMPASS_20220502_MSM_SipA_20cm 686.1893 5431.449
## 28      COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549
## 29      COMPASS_20220502_MSM_SipA_45cm_Dup 398.1761 6525.085
## 30 COMPASS_20220502_MSM_SipA_45cm_Spike 606.2450 6339.224
## 31      COMPASS_20220502_MSM_SipB_10cm 747.5236 5169.818
```

```
# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06        #molecular weight of sulfur: 32.06
```

```
# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)
```

```
# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000
```

```
head(sampledat)
```

```
##           Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 26      COMPASS_20220502_MSM_SipA_10cm 924.7531 5530.741 28.84445 156.0153
## 27      COMPASS_20220502_MSM_SipA_20cm 686.1893 5431.449 21.40328 153.2144
## 28      COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753
## 29      COMPASS_20220502_MSM_SipA_45cm_Dup 398.1761 6525.085 12.41972 184.0645
```

```

## 30 COMPASS_20220502_MSM_SipA_45cm_Spike 606.2450 6339.224 18.90970 178.8215
## 31      COMPASS_20220502_MSM_SipB_10cm 747.5236 5169.818 23.31639 145.8341
##   salinity
## 26  9.994075
## 27  9.814655
## 28 11.945288
## 29 11.790855
## 30 11.455003
## 31  9.341887

```

Pull out dups and check with percent difference

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##                               Sample_ID  SO4_ppm   Cl_ppm   SO4_mM   Cl_mM
## 26      COMPASS_20220502_MSM_SipA_10cm 924.7531 5530.741 28.84445 156.0153
## 27      COMPASS_20220502_MSM_SipA_20cm 686.1893 5431.449 21.40328 153.2144
## 28      COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753
## 29      COMPASS_20220502_MSM_SipA_45cm_Dup 398.1761 6525.085 12.41972 184.0645
## 30 COMPASS_20220502_MSM_SipA_45cm_Spike 606.2450 6339.224 18.90970 178.8215
## 31      COMPASS_20220502_MSM_SipB_10cm 747.5236 5169.818 23.31639 145.8341
##   salinity
## 26  9.994075
## 27  9.814655
## 28 11.945288
## 29 11.790855
## 30 11.455003
## 31  9.341887

```

```

#pull out any rows that have "dup" in the SampleID column
dups <- sampledat %>%
  filter(str_detect(Sample_ID, "Dup"))      #have to change this to match data
head(dups)

```

```

##                               Sample_ID  SO4_ppm   Cl_ppm   SO4_mM   Cl_mM
## 1 COMPASS_20220502_MSM_SipA_45cm_Dup 398.1761 6525.085 12.41972 184.0645
## 2 COMPASS_20220505_GWI_SipA_20cm_Dup 1069.4710 8909.081 33.35842 251.3140
##   salinity
## 1 11.79085
## 2 16.09874

```

```

#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledat %>%
  filter(!str_detect(Sample_ID, "Dup")) %>%
  filter(!str_detect(Sample_ID, "Spike"))
head(sampledat2)

```

```

##                               Sample_ID  SO4_ppm   Cl_ppm   SO4_mM   Cl_mM  salinity
## 1 COMPASS_20220502_MSM_SipA_10cm 924.7531 5530.741 28.84445 156.0153  9.994075

```

```

## 2 COMPASS_20220502_MSM_SipA_20cm 686.1893 5431.449 21.40328 153.2144 9.814655
## 3 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.945288
## 4 COMPASS_20220502_MSM_SipB_10cm 747.5236 5169.818 23.31639 145.8341 9.341887
## 5 COMPASS_20220502_MSM_SipB_20cm 706.2678 5130.657 22.02956 144.7294 9.271123
## 6 COMPASS_20220502_MSM_SipC_10cm 869.9081 5414.163 27.13375 152.7267 9.783419

#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_Dup","",as.character(dups$Sample_ID))
dups <- dups[ , -c(2,3)]
colnames(dups) <- c('Sample_ID', 'SO4_mM_dup', "Cl_mM_dup", "salinity_dup")
head(dups)

##                               Sample_ID SO4_mM_dup Cl_mM_dup salinity_dup
## 1 COMPASS_20220502_MSM_SipA_45cm    12.41972   184.0645     11.79085
## 2 COMPASS_20220505_GWI_SipA_20cm    33.35842   251.3140     16.09874

#put it back together with the old data set and look for duplicates
QAdups <- merge(sampledat2, dups)
head(QAdups)

##                               Sample_ID SO4_ppm   Cl_ppm   SO4_mM   Cl_mM salinity
## 1 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.94529
## 2 COMPASS_20220505_GWI_SipA_20cm 1068.0465 8895.024 33.31399 250.9175 16.07333
##   SO4_mM_dup Cl_mM_dup salinity_dup
## 1    12.41972   184.0645     11.79085
## 2    33.35842   251.3140     16.09874

QAdups$SO4_dups_chk <- ((abs(QAdups$SO4_mM-QAdups$SO4_mM_dup))/((QAdups$SO4_mM+QAdups$SO4_mM_dup)/2))*100
QAdups$SO4_dups_flag <- ifelse(QAdups$SO4_dups_chk <10, 'YES', 'NO, rerun')

QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')

head(QAdups)

##                               Sample_ID SO4_ppm   Cl_ppm   SO4_mM   Cl_mM salinity
## 1 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.94529
## 2 COMPASS_20220505_GWI_SipA_20cm 1068.0465 8895.024 33.31399 250.9175 16.07333
##   SO4_mM_dup Cl_mM_dup salinity_dup SO4_dups_chk SO4_dups_flag Cl_dups_chk
## 1    12.41972   184.0645     11.79085      0.3972243        YES    1.3012528
## 2    33.35842   251.3140     16.09874      0.1332855        YES    0.1579119
##   Cl_dups_flag
## 1           YES
## 2           YES

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = SO4_dups_chk, fill=SO4_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between SO4 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

```

```

## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.


```

```

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

```

```

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES        2  YES        2      2     100      100

```

Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)


##                               Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM salinity
## 1 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.94529
## 2 COMPASS_20220505_GWI_SipA_20cm 1068.0465 8895.024 33.31399 250.9175 16.07333
##   S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag Cl_dups_chk
## 1    12.41972 184.0645    11.79085  0.3972243      YES  1.3012528
## 2    33.35842 251.3140    16.09874  0.1332855      YES  0.1579119
##   Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1           YES  0.28088003      YES
## 2           YES  0.09424705      YES

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x = "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x = "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES        2    YES       2     2      100      100
```

Pull out spikes and check

```
#Show me the data that we have from the calculations
head(sampledat)
```

```

##                                     Sample_ID S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 26      COMPASS_20220502_MSM_SipA_10cm 924.7531 5530.741 28.84445 156.0153
## 27      COMPASS_20220502_MSM_SipA_20cm 686.1893 5431.449 21.40328 153.2144
## 28      COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753
## 29      COMPASS_20220502_MSM_SipA_45cm_Dup 398.1761 6525.085 12.41972 184.0645
## 30 COMPASS_20220502_MSM_SipA_45cm_Spike 606.2450 6339.224 18.90970 178.8215
## 31      COMPASS_20220502_MSM_SipB_10cm 747.5236 5169.818 23.31639 145.8341
##   salinity
## 26  9.994075
## 27  9.814655
## 28 11.945288
## 29 11.790855
## 30 11.455003
## 31  9.341887

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "Spike"))      #have to change this to match data
head(spks)

##                                     Sample_ID S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 1 COMPASS_20220502_MSM_SipA_45cm_Spike 606.245 6339.224 18.9097 178.8215
##   salinity
## 1    11.455

#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_Spike","",as.character(spks$Sample_ID))
spks <- spks[,-c(2,3, 5,6)]
colnames(spks) <- c('Sample_ID', 'S04_mM_spk')
head(spks)

##                                     Sample_ID S04_mM_spk
## 1 COMPASS_20220502_MSM_SipA_45cm     18.9097

#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat, spks)
head(QAspks)

##                                     Sample_ID S04_ppm   Cl_ppm   S04_mM   Cl_mM salinity
## 1 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.94529
##   S04_mM_spk
## 1    18.9097

#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvol <- 10                 # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvol      # mmoles of S04
head(QAspks)

```

```

##                               Sample_ID  SO4_ppm   Cl_ppm   SO4_mM     Cl_mM salinity
## 1 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.94529
##   SO4_mM_spk SO4_spk_Conc
## 1      18.9097 7.797879e-05

#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

##                               Sample_ID  SO4_ppm   Cl_ppm   SO4_mM     Cl_mM salinity
## 1 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.94529
##   SO4_mM_spk SO4_spk_Conc Dilution SampleVol
## 1      18.9097 7.797879e-05      1    1e-06

#gives us the total SO4 in the sample in mmoles
QAspks$SO4_Total_unspkd <- (QAspks$SO4_mM/QAspks$Dilution)*(QAspks$SampleVol)

##total SO4 in spiked sample in mmoles
QAspks$SO4_Total_spkd <- (QAspks$SO4_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)

QAspks$SO4_expctd_spkd <- (QAspks$SO4_Total_unspkd + QAspks$SO4_spk_Conc)

```

```

QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES' , 'NO', rerun


head(QAspks)

##                               Sample_ID  S04_ppm   Cl_ppm   S04_mM   Cl_mM salinity
## 1 COMPASS_20220502_MSM_SipA_45cm 399.7609 6610.549 12.46915 186.4753 11.94529
##   S04_mM_spk S04_spk_Conc Dilution SampleVol S04_Total_unspkd S04_Total_spkd
## 1     18.9097 7.797879e-05         1      1e-06     1.246915e-05  0.0002080067
##   S04_expctd_spkd spk_recovery S04_spks_flag
## 1     9.044794e-05        229.974      NO, rerun

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)

spksbar

```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)

## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>          <int>   <int>    <dbl>
## 1 NO, rerun      1       1     100

```

Make final dataframe with IDs

```

#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_', fixed=TRUE)))
colnames(IDs) <- c("Analysis_No", "Site", "Zone", "Replicate", "Depth")
head(IDs)

## Analysis_No Site Zone Replicate Depth
## 1 COMPASS 20220502 MSM SipA 10cm
## 2 COMPASS 20220502 MSM SipA 20cm
## 3 COMPASS 20220502 MSM SipA 45cm
## 4 COMPASS 20220502 MSM SipB 10cm
## 5 COMPASS 20220502 MSM SipB 20cm
## 6 COMPASS 20220502 MSM SipC 10cm

#rejoin them to the dataframe
alldat <- cbind(IDs, sampledat2)
head(alldat)

## Analysis_No Site Zone Replicate Depth Sample_ID
## 1 COMPASS 20220502 MSM SipA 10cm COMPASS_20220502_MSM_SipA_10cm
## 2 COMPASS 20220502 MSM SipA 20cm COMPASS_20220502_MSM_SipA_20cm
## 3 COMPASS 20220502 MSM SipA 45cm COMPASS_20220502_MSM_SipA_45cm
## 4 COMPASS 20220502 MSM SipB 10cm COMPASS_20220502_MSM_SipB_10cm
## 5 COMPASS 20220502 MSM SipB 20cm COMPASS_20220502_MSM_SipB_20cm
## 6 COMPASS 20220502 MSM SipC 10cm COMPASS_20220502_MSM_SipC_10cm
##   S04_ppm Cl_ppm S04_mM Cl_mM salinity
## 1 924.7531 5530.741 28.84445 156.0153 9.994075
## 2 686.1893 5431.449 21.40328 153.2144 9.814655
## 3 399.7609 6610.549 12.46915 186.4753 11.945288
## 4 747.5236 5169.818 23.31639 145.8341 9.341887
## 5 706.2678 5130.657 22.02956 144.7294 9.271123
## 6 869.9081 5414.163 27.13375 152.7267 9.783419

```

Make final dataframe with IDs

```
#setwd("S:/Biogeochemistry/People/Wilson (Steph)/Data/Dionex/Final Data Files")      #Change wd  
write.csv(alldat, file="Processed Data/COMPASS_SynopticCB_PW_ProCESSED_Cl_S04_202205.csv")    #Change
```

END