

# Dionex\_COMPASS\_June2023

Stephanie J. Wilson

2023-06-23

## Daily Set up

Read in and Format the raw data - change wd & file names

```
# SULFATE DATA:
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202306_S04.txt", sep='\t' , header=T, skip=3)
head(Sdat)
```

```
##      X      X.1      X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank      Unknown 4.473  0.0159    3.37   0.0154    0.09
## 2 2 Lab Blank      Unknown 4.490  0.0155    3.53   0.0151    0.08
## 3 3 Lab Blank      Unknown 4.493  0.0156    3.53   0.0151    0.08
## 4 4 Lab Blank      Unknown 4.497  0.0163    3.62   0.0158    0.09
## 5 5 Standard 1 Calibration Standard 4.503  0.7135    9.66   0.6906    5.08
## 6 6 Standard 2 Calibration Standard 4.497  1.2518    8.87   1.2118    8.59
```

```
## Only keep the columns that we need
```

```
Sdat <- Sdat[,c(2,5)] # dont need this here
head(Sdat)
```

```
##      X.1 IC.S04.1
## 1 Lab Blank 0.0159
## 2 Lab Blank 0.0155
## 3 Lab Blank 0.0156
## 4 Lab Blank 0.0163
## 5 Standard 1 0.7135
## 6 Standard 2 1.2518
```

```
## Name the columns correctly
```

```
colnames(Sdat) <- c("Sample_ID", "S04_ppm")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID S04_ppm
## 1 Lab Blank  0.0159
## 2 Lab Blank  0.0155
## 3 Lab Blank  0.0156
## 4 Lab Blank  0.0163
## 5 Standard 1  0.7135
## 6 Standard 2  1.2518
```

```
#Chloride data
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Cldat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202306_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X      X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank      Unknown 3.543 0.0226 6.28 0.0287 0.14
## 2 2 Lab Blank      Unknown 3.540 0.0096 2.86 0.0122 0.10
## 3 3 Lab Blank      Unknown 3.543 0.0095 2.81 0.0120 0.10
## 4 4 Lab Blank      Unknown 3.543 0.0099 2.88 0.0126 0.11
## 5 5 Standard 1 Calibration Standard 3.550 5.0880 90.34 6.4594 56.89
## 6 6 Standard 2 Calibration Standard 3.550 9.8032 91.13 12.4456 108.66
```

```
## Only keep the columns that we need
```

```
Cldat <- Cldat[,c(2,5)]
head(Cldat)
```

```
##      X.1 IC.Cl.1
## 1 Lab Blank 0.0226
## 2 Lab Blank 0.0096
## 3 Lab Blank 0.0095
## 4 Lab Blank 0.0099
## 5 Standard 1 5.0880
## 6 Standard 2 9.8032
```

```
## Name the columns correctly
```

```
colnames(Cldat) <- c("Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##      Sample_ID Cl_ppm
## 1 Lab Blank 0.0226
## 2 Lab Blank 0.0096
## 3 Lab Blank 0.0095
## 4 Lab Blank 0.0099
## 5 Standard 1 5.0880
## 6 Standard 2 9.8032
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##   Sample_ID   S04_ppm   Cl_ppm
## 1         1790.3015 13351.2461
## 2         1790.3015  4162.4856
## 3         1790.3015    0.0056
## 4         1790.3015 4096.1384
## 5         1790.3015      NA
## 6         462.8695 13351.2461
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##           Sample_ID   S04_ppm   Cl_ppm
## 26  1_GCW_202306_UP_LysA_20cm    5.2274    7.1807
## 27  10_GCW_202306_TR_LysB_20cm    4.9649    6.2568
## 28  10_GCW_202306_TR_LysB_20cm_dup    5.0878    6.5878
## 29  100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880
## 30  100_GWI_202306_WC_LysA_10cm_dup 1192.1473 12081.0024
## 31  101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)
```

```
##           Sample_ID   S04_ppm   Cl_ppm
## 26  1_GCW_202306_UP_LysA_20cm    5.2274    7.1807
## 27  10_GCW_202306_TR_LysB_20cm    4.9649    6.2568
## 28  10_GCW_202306_TR_LysB_20cm_dup    5.0878    6.5878
## 29  100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880
## 30  100_GWI_202306_WC_LysA_10cm_dup 1192.1473 12081.0024
## 31  101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142
```

```
## If there is an n.a. in the dataframe that is because there was no peak, which
# would mean there was no Cl or S04 and we want to know that so make all n.a.'s into zeros
all_dat[is.na(all_dat)] <- 0
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grepl("Standard", all_dat$Sample_ID),]
stds <- stds[-c(17),]
head(stds)
```

```
##   Sample_ID S04_ppm Cl_ppm
## 621 Standard 1  0.7831 5.4746
## 622 Standard 1  0.7831 5.0880
## 623 Standard 1  0.7831 5.5265
```

```
## 624 Standard 1 0.7831 5.1245
## 625 Standard 1 0.7135 5.4746
## 626 Standard 1 0.7135 5.0880
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(SO4_ppm), sd = sd(SO4_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean    sd    cv flag
##   <fct>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1 0.757 0.0352 4.66 YES
## 2 Standard 2 1.33 0.0523 3.94 YES
## 3 Standard 3 2.33 0.0867 3.73 YES
## 4 Standard 4 10.7 0.293 2.74 YES
## 5 Standard 5 20.4 0.430 2.11 YES
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean    sd    cv flag
##   <fct>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1 5.30 0.205 3.86 YES
## 2 Standard 2 10.3 0.382 3.71 YES
## 3 Standard 3 20.9 0.792 3.79 YES
## 4 Standard 4 104. 3.04 2.93 YES
## 5 Standard 5 206. 4.84 2.35 YES
```

## Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##           Sample_ID SO4_ppm Cl_ppm
## 26 1_GCW_202306_UP_LysA_20cm 5.2274 7.1807
## 27 10_GCW_202306_TR_LysB_20cm 4.9649 6.2568
## 28 10_GCW_202306_TR_LysB_20cm_dup 5.0878 6.5878
## 29 100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880
## 30 100_GWI_202306_WC_LysA_10cm_dup 1192.1473 12081.0024
## 31 101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142
```

```

# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06      #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)

```

```

##              Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 26      1_GCW_202306_UP_LysA_20cm    5.2274    7.1807  0.1630505  0.2025585
## 27     10_GCW_202306_TR_LysB_20cm    4.9649    6.2568  0.1548628  0.1764965
## 28    10_GCW_202306_TR_LysB_20cm_dup    5.0878    6.5878  0.1586962  0.1858336
## 29     100_GWI_202306_WC_LysA_10cm  1136.1780  11561.9880  35.4391142  326.1491678
## 30    100_GWI_202306_WC_LysA_10cm_dup  1192.1473  12081.0024  37.1848815  340.7899126
## 31     101_GWI_202306_WC_LysA_20cm  1325.5768  10940.2142  41.3467498  308.6097094
##      salinity
## 26  0.01300152
## 27  0.01133204
## 28  0.01193015
## 29 20.89253832
## 30 21.83039734
## 31 19.76899306

```

Pull out dups and check with percent difference

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##              Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 26      1_GCW_202306_UP_LysA_20cm    5.2274    7.1807  0.1630505  0.2025585
## 27     10_GCW_202306_TR_LysB_20cm    4.9649    6.2568  0.1548628  0.1764965
## 28    10_GCW_202306_TR_LysB_20cm_dup    5.0878    6.5878  0.1586962  0.1858336
## 29     100_GWI_202306_WC_LysA_10cm  1136.1780  11561.9880  35.4391142  326.1491678
## 30    100_GWI_202306_WC_LysA_10cm_dup  1192.1473  12081.0024  37.1848815  340.7899126
## 31     101_GWI_202306_WC_LysA_20cm  1325.5768  10940.2142  41.3467498  308.6097094
##      salinity
## 26  0.01300152
## 27  0.01133204
## 28  0.01193015

```

```
## 29 20.89253832
## 30 21.83039734
## 31 19.76899306
```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledats %>%
  filter(str_detect(Sample_ID, "dup")) #have to change this to match data
head(dups)
```

```
##           Sample_ID    S04_ppm    Cl_ppm    S04_mM    Cl_mM
## 1 10_GCW_202306_TR_LysB_20cm_dup    5.0878    6.5878  0.1586962  0.1858336
## 2 100_GWI_202306_WC_LysA_10cm_dup 1192.1473 12081.0024 37.1848815 340.7899126
## 3 110_GWI_202306_SW_C_dup    1716.1138 12692.8628 53.5281909 358.0497264
## 4 120_MSM_202306_WC_RHZ_Col1_dup 1093.4677  9805.4513 34.1069152 276.5994725
## 5 130_MSM_202306_UP_PPR_3_dup    77.8144   975.9884  2.4271491  27.5314076
## 6 140_MSM_202306_TR_PPR_1_dup   158.3756  1642.2791  4.9399750  46.3266319
##           salinity
## 1 0.01193015
## 2 21.83039734
## 3 22.93602908
## 4 17.71847650
## 5 1.76363704
## 6 2.96762433
```

```
#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledats %>%
  filter(!str_detect(Sample_ID, "dup")) %>%
  filter(!str_detect(Sample_ID, "spk"))
head(sampledat2)
```

```
##           Sample_ID    S04_ppm    Cl_ppm    S04_mM    Cl_mM
## 1 1_GCW_202306_UP_LysA_20cm    5.2274    7.1807  0.1630505  0.2025585
## 2 10_GCW_202306_TR_LysB_20cm    4.9649    6.2568  0.1548628  0.1764965
## 3 100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880 35.4391142 326.1491678
## 4 101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142 41.3467498 308.6097094
## 5 102_GWI_202306_WC_LysA_45cm 1374.3504 11492.1490 42.8680724 324.1790973
## 6 103_GWI_202306_WC_LysB_10cm 1309.7705 12131.0664 40.8537274 342.2021551
##           salinity
## 1 0.01300152
## 2 0.01133204
## 3 20.89253832
## 4 19.76899306
## 5 20.76633924
## 6 21.92086298
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup", "", as.character(dups$Sample_ID))
dups <- dups[, -c(2,3)]
colnames(dups) <- c('Sample_ID', 'S04_mM_dup', 'Cl_mM_dup', 'salinity_dup')
head(dups)
```

```
##           Sample_ID S04_mM_dup  Cl_mM_dup salinity_dup
```

```
## 1 10_GCW_202306_TR_LysB_20cm 0.1586962 0.1858336 0.01193015
## 2 100_GWI_202306_WC_LysA_10cm 37.1848815 340.7899126 21.83039734
## 3 110_GWI_202306_SW_C 53.5281909 358.0497264 22.93602908
## 4 120_MSM_202306_WC_RHZ_Col1 34.1069152 276.5994725 17.71847650
## 5 130_MSM_202306_UP_PPR_3 2.4271491 27.5314076 1.76363704
## 6 140_MSM_202306_TR_PPR_1 4.9399750 46.3266319 2.96762433
```

*#put it back together with the old data set and look for duplicates*

```
QAdups <- merge(sampled2, dups)
head(QAdups)
```

```
##          Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 10_GCW_202306_TR_LysB_20cm 4.9649 6.2568 0.1548628 0.1764965
## 2 100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880 35.4391142 326.1491678
## 3 120_MSM_202306_WC_RHZ_Col1 1094.7124 9713.3353 34.1457392 274.0009958
## 4 140_MSM_202306_TR_PPR_1 164.8393 1705.2825 5.1415876 48.1038787
## 5 150_MSM_202306_TR_PPR_11 820.7486 6672.6120 25.6003930 188.2260085
## 6 160_MSM_202306_WC_PPR_9 1026.1569 8888.6427 32.0073893 250.7374528
## salinity S04_mM_dup Cl_mM_dup salinity_dup
## 1 0.01133204 0.1586962 0.1858336 0.01193015
## 2 20.89253832 37.1848815 340.7899126 21.83039734
## 3 17.55202289 34.1069152 276.5994725 17.71847650
## 4 3.08147148 4.9399750 46.3266319 2.96762433
## 5 12.05743588 25.5830724 188.8710635 12.09875691
## 6 16.06180336 33.0801029 256.5228039 16.43240225
```

```
QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*100
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')
```

```
QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')
```

```
head(QAdups)
```

```
##          Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 10_GCW_202306_TR_LysB_20cm 4.9649 6.2568 0.1548628 0.1764965
## 2 100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880 35.4391142 326.1491678
## 3 120_MSM_202306_WC_RHZ_Col1 1094.7124 9713.3353 34.1457392 274.0009958
## 4 140_MSM_202306_TR_PPR_1 164.8393 1705.2825 5.1415876 48.1038787
## 5 150_MSM_202306_TR_PPR_11 820.7486 6672.6120 25.6003930 188.2260085
## 6 160_MSM_202306_WC_PPR_9 1026.1569 8888.6427 32.0073893 250.7374528
## salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 0.01133204 0.1586962 0.1858336 0.01193015 2.44511425 YES
## 2 20.89253832 37.1848815 340.7899126 21.83039734 4.80768731 YES
## 3 17.55202289 34.1069152 276.5994725 17.71847650 0.11376577 YES
## 4 3.08147148 4.9399750 46.3266319 2.96762433 3.99962997 YES
## 5 12.05743588 25.5830724 188.8710635 12.09875691 0.06768064 YES
## 6 16.06180336 33.0801029 256.5228039 16.43240225 3.29622060 YES
## Cl_dups_chk Cl_dups_flag
## 1 5.1539168 YES
## 2 4.3904294 YES
## 3 0.9438701 YES
## 4 3.7641369 YES
```

```
## 5 0.3421162 YES
## 6 2.2810189 YES
```


*#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this*

```
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)
```

```
ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```





```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         14  YES         14     14         100        100

```

## Pull out dups and check with cv

```

#the cv
# calculate the sd between the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)

```

```

##           Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 1 10_GCW_202306_TR_LysB_20cm    4.9649    6.2568  0.1548628  0.1764965
## 2 100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880 35.4391142 326.1491678
## 3 120_MSM_202306_WC_RHZ_Col1 1094.7124  9713.3353 34.1457392 274.0009958
## 4  140_MSM_202306_TR_PPR_1   164.8393  1705.2825  5.1415876  48.1038787
## 5  150_MSM_202306_TR_PPR_11  820.7486  6672.6120 25.6003930 188.2260085
## 6  160_MSM_202306_WC_PPR_9 1026.1569  8888.6427 32.0073893 250.7374528
##      salinity S04_mM_dup   Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1  0.01133204  0.1586962   0.1858336   0.01193015   2.44511425         YES
## 2 20.89253832 37.1848815 340.7899126   21.83039734   4.80768731         YES
## 3 17.55202289 34.1069152 276.5994725   17.71847650   0.11376577         YES
## 4  3.08147148  4.9399750  46.3266319    2.96762433   3.99962997         YES
## 5 12.05743588 25.5830724 188.8710635   12.09875691   0.06768064         YES
## 6 16.06180336 33.0801029 256.5228039   16.43240225   3.29622060         YES
##      Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1    5.1539168         YES  1.72895687         YES
## 2    4.3904294         YES  3.39954830         YES

```


```
## 3 0.9438701 YES 0.08044455 YES
## 4 3.7641369 YES 2.82816547 YES
## 5 0.3421162 YES 0.04785744 YES
## 6 2.2810189 YES 2.33077994 YES
```

*#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this*

```
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar,ncol=2, nrow=1)
```



*#check for percent of no, reruns to see if it would warrant reruns*

```
Perc_dups <- QAdups %>%
```

```

group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES          14  YES          14    14         100         100

```

## Pull out spikes and check

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 26  1_GCW_202306_UP_LysA_20cm  5.2274  7.1807  0.1630505  0.2025585
## 27  10_GCW_202306_TR_LysB_20cm  4.9649  6.2568  0.1548628  0.1764965
## 28  10_GCW_202306_TR_LysB_20cm_dup  5.0878  6.5878  0.1586962  0.1858336
## 29  100_GWI_202306_WC_LysA_10cm 1136.1780 11561.9880 35.4391142 326.1491678
## 30 100_GWI_202306_WC_LysA_10cm_dup 1192.1473 12081.0024 37.1848815 340.7899126
## 31  101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142 41.3467498 308.6097094
##      salinity
## 26  0.01300152
## 27  0.01133204
## 28  0.01193015
## 29 20.89253832
## 30 21.83039734
## 31 19.76899306

```

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)

```

```

##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 101_GWI_202306_WC_LysA_20cm_spk 1614.9188 10400.4899 50.371765 293.3847645
## 2  11_GCW_202306_TR_LysB_45cm_spk   89.1234   4.2323  2.779894  0.1193879
## 3  111_MSM_202306_UP_RHZ_SF6_spk  222.0675  1072.2663  6.926622  30.2472863
## 4 121_MSM_202306_WC_RHZ_Co12_spk 1065.4168  9059.6735 33.231965 255.5620169
## 5   131_MSM_202306_UP_PPR_4_spk  186.3889   977.1367  5.813752  27.5637997
## 6   141_MSM_202306_TR_PPR_2_spk  324.6354  2338.5877 10.125870  65.9686234
##      salinity

```

```
## 1 18.793711249
## 2 0.007673766
## 3 1.937611204
## 4 16.370856014
## 5 1.765712017
## 6 4.225853974
```

```
#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[ , -c(2,3, 5,6)]
colnames(spks) <- c('Sample_ID', 'SO4_mM_spk')
head(spks)
```

```
##           Sample_ID SO4_mM_spk
## 1 101_GWI_202306_WC_LysA_20cm 50.371765
## 2 11_GCW_202306_TR_LysB_45cm 2.779894
## 3 111_MSM_202306_UP_RHZ_SF6 6.926622
## 4 121_MSM_202306_WC_RHZ_Co12 33.231965
## 5 131_MSM_202306_UP_PPR_4 5.813752
## 6 141_MSM_202306_TR_PPR_2 10.125870
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat, spks)
head(QAspks)
```

```
##           Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
## 1 101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142 41.3467498 308.609709
## 2 11_GCW_202306_TR_LysB_45cm 5.5728 4.2406 0.1738241 0.119622
## 3 111_MSM_202306_UP_RHZ_SF6 132.2880 1071.7023 4.1262633 30.231377
## 4 121_MSM_202306_WC_RHZ_Co12 875.2224 8937.4618 27.2995134 252.114578
## 5 131_MSM_202306_UP_PPR_4 99.0902 1030.9977 3.0907735 29.083151
## 6 141_MSM_202306_TR_PPR_2 237.8114 2340.9398 7.4176981 66.034973
##           salinity SO4_mM_spk
## 1 19.768993059 50.371765
## 2 0.007688764 2.779894
## 3 1.936592056 6.926622
## 4 16.150019473 33.231965
## 5 1.863038844 5.813752
## 6 4.230104219 10.125870
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw) # in mM
spkvol <- 10 # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$SO4_spk_Conc <- (spkconc)*spkvol # mmoles of SO4
head(QAspks)
```

```
##           Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
## 1 101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142 41.3467498 308.609709
## 2 11_GCW_202306_TR_LysB_45cm 5.5728 4.2406 0.1738241 0.119622
## 3 111_MSM_202306_UP_RHZ_SF6 132.2880 1071.7023 4.1262633 30.231377
```

```
## 4 121_MSM_202306_WC_RHZ_Co12 875.2224 8937.4618 27.2995134 252.114578
## 5 131_MSM_202306_UP_PPR_4 99.0902 1030.9977 3.0907735 29.083151
## 6 141_MSM_202306_TR_PPR_2 237.8114 2340.9398 7.4176981 66.034973
## salinity S04_mM_spk S04_spk_Conc
## 1 19.768993059 50.371765 7.797879e-05
## 2 0.007688764 2.779894 7.797879e-05
## 3 1.936592056 6.926622 7.797879e-05
## 4 16.150019473 33.231965 7.797879e-05
## 5 1.863038844 5.813752 7.797879e-05
## 6 4.230104219 10.125870 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
```

```
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 50, QAspks$Dilution)
```

```
#Set Sample volumes in uL
```

```
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 1462, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 1501, QAspks$SampleVol)
```

```
#change sample volume to L
```

```
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)
```

```
## Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142 41.3467498 308.609709
## 2 11_GCW_202306_TR_LysB_45cm 5.5728 4.2406 0.1738241 0.119622
## 3 111_MSM_202306_UP_RHZ_SF6 132.2880 1071.7023 4.1262633 30.231377
## 4 121_MSM_202306_WC_RHZ_Co12 875.2224 8937.4618 27.2995134 252.114578
## 5 131_MSM_202306_UP_PPR_4 99.0902 1030.9977 3.0907735 29.083151
## 6 141_MSM_202306_TR_PPR_2 237.8114 2340.9398 7.4176981 66.034973
```

```
##      salinity S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 19.768993059 50.371765 7.797879e-05      200 0.001462
## 2 0.007688764 2.779894 7.797879e-05        1 0.000001
## 3 1.936592056 6.926622 7.797879e-05       50 0.001501
## 4 16.150019473 33.231965 7.797879e-05      100 0.001475
## 5 1.863038844 5.813752 7.797879e-05       50 0.001501
## 6 4.230104219 10.125870 7.797879e-05       50 0.001501
```

*#gives us the total SO4 in the sample in mmoles*

```
QAspks$S04_Total_unspkd <- (QAspks$S04_mM/QAspks$Dilution)*(QAspks$SampleVol)
```

*##total SO4 in spiked sample in mmoles*

```
QAspks$S04_Total_spkd <- (QAspks$S04_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)
```

```
QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
```

```
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun')
```


```
head(QAspks)
```

```
##      Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 101_GWI_202306_WC_LysA_20cm 1325.5768 10940.2142 41.3467498 308.609709
## 2 11_GCW_202306_TR_LysB_45cm 5.5728 4.2406 0.1738241 0.119622
## 3 111_MSM_202306_UP_RHZ_SF6 132.2880 1071.7023 4.1262633 30.231377
## 4 121_MSM_202306_WC_RHZ_Co12 875.2224 8937.4618 27.2995134 252.114578
## 5 131_MSM_202306_UP_PPR_4 99.0902 1030.9977 3.0907735 29.083151
## 6 141_MSM_202306_TR_PPR_2 237.8114 2340.9398 7.4176981 66.034973
##      salinity S04_mM_spk S04_spk_Conc Dilution SampleVol S04_Total_unspkd
## 1 19.768993059 50.371765 7.797879e-05      200 0.001462 3.022447e-04
## 2 0.007688764 2.779894 7.797879e-05        1 0.000001 1.738241e-07
## 3 1.936592056 6.926622 7.797879e-05       50 0.001501 1.238704e-04
## 4 16.150019473 33.231965 7.797879e-05      100 0.001475 4.026678e-04
## 5 1.863038844 5.813752 7.797879e-05       50 0.001501 9.278502e-05
## 6 4.230104219 10.125870 7.797879e-05       50 0.001501 2.226793e-04
##      S04_Total_spkd S04_expctd_spkd spk_recovery S04_spks_flag
## 1 3.707362e-04 3.802235e-04 97.50480 YES
## 2 3.057883e-05 7.815261e-05 39.12708 NO, rerun
## 3 2.093225e-04 2.018492e-04 103.70242 YES
## 4 4.934947e-04 4.806466e-04 102.67308 YES
## 5 1.756916e-04 1.707638e-04 102.88573 YES
## 6 3.060038e-04 3.006581e-04 101.77800 YES
```

*#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on*

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

```
spksbar
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 2 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 NO, rerun      3     16    18.8
## 2 YES           13     16    81.2
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_ ', fixed=TRUE)))
```

```
## Warning in rbind(c("1", "GCW", "202306", "UP", "LysA", "20cm"), c("10", :
## number of columns of result is not a multiple of vector length (arg 11)
```

```
colnames(IDs) <- c("Analysis_No", "Site", "Date", "Zone", "Replicate", "Depth")
head(IDs)
```

```
##   Analysis_No Site   Date Zone Replicate Depth
## 1           1  GCW 202306  UP      LysA   20cm
## 2          10  GCW 202306  TR      LysB   20cm
## 3         100  GWI 202306  WC      LysA   10cm
## 4         101  GWI 202306  WC      LysA   20cm
## 5         102  GWI 202036  WC      LysA   45cm
## 6         103  GWI 202306  WC      LysB   10cm
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampled2)
head(alldat)
```

```
##   Analysis_No Site   Date Zone Replicate Depth      Sample_ID
## 1           1  GCW 202306  UP      LysA   20cm  1_GCW_202306_UP_LysA_20cm
## 2          10  GCW 202306  TR      LysB   20cm 10_GCW_202306_TR_LysB_20cm
## 3         100  GWI 202306  WC      LysA   10cm 100_GWI_202306_WC_LysA_10cm
## 4         101  GWI 202306  WC      LysA   20cm 101_GWI_202306_WC_LysA_20cm
## 5         102  GWI 202036  WC      LysA   45cm 102_GWI_202036_WC_LysA_45cm
## 6         103  GWI 202306  WC      LysB   10cm 103_GWI_202306_WC_LysB_10cm
##   S04_ppm   Cl_ppm   S04_mM   Cl_mM   salinity
## 1   5.2274   7.1807  0.1630505  0.2025585  0.01300152
## 2   4.9649   6.2568  0.1548628  0.1764965  0.01133204
## 3 1136.1780 11561.9880 35.4391142 326.1491678 20.89253832
## 4 1325.5768 10940.2142 41.3467498 308.6097094 19.76899306
## 5 1374.3504 11492.1490 42.8680724 324.1790973 20.76633924
## 6 1309.7705 12131.0664 40.8537274 342.2021551 21.92086298
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_S04_202306.csv")
```

*#C*

END