

# Synoptic\_CB\_Nutrients\_202305

May 2023

2025-07-01

## Contents

0.1	Import Data & Clean . . . . .	2
0.2	Assessing standard Curves . . . . .	2
0.3	Dilution Corrections - ensure the latest dilution is kept . . . . .	7
0.4	Temporary Chunk . . . . .	7
0.5	Performance Check . . . . .	7
0.6	Analyze the Check Standards . . . . .	9
0.7	Analyze Blanks . . . . .	10
0.8	Analyze Duplicates . . . . .	11
0.9	Spikes . . . . .	12
0.10	Matrix Effects . . . . .	13
0.11	Unit Converted Data Column Added (mg/L to uM ) . . . . .	13
0.12	Sample Flagging - Within range of standard curve . . . . .	13
0.13	Pull out sample id information . . . . .	13
0.14	Pulling Out Rhizon Samples . . . . .	13
0.15	Check to see if samples run match metadata & merge info . . . . .	13
0.16	Visualize Data . . . . .	14
0.17	Export Processed Data . . . . .	16

##Run Information

```
cat("Run Information: Input by User") #lets you know what section you're in
```

```
## Run Information: Input by User
```

```
#set the run date & user name
run_date <- "20240114"
sample_year <- 2023
sample_month <- 05
user <- "Stephanie Wilson"
```

```
#identify the files you want to read in
```

```

#read in as a list to accommodate multiple runs in a month
NOx_files <- c("Raw Data/SEAL_COMPASS_Synoptic_NOx_May2023_1.csv",
               "Raw Data/SEAL_COMPASS_Synoptic_NOx_May2023_2_redo.csv",
               "Raw Data/SEAL_COMPASS_Synoptic_NOx_May2023_3_redoA.csv")
NH3_PO4_files <- c("Raw Data/SEAL_COMPASS_Synoptic_NH3_PO4_May2023_1redopo4.csv",
                   "Raw Data/SEAL_COMPASS_Synoptic_NH3_PO4_May2023_3.csv")

# Define the file path for QAQC log file - NO Need to change just check year
file_path <- "Raw Data/SEAL_COMPASS_Synoptic_QAQC_Log_2023.csv"
final_path <- "Processed Data/COMPASS_Synoptic_Nutrients_202305.csv"

#record any notes about the run or anything other info here:
run_notes <- "Autospikes for PO4 not included in this run. "

#duplicate sample names to be changed
#list the sample iDs that are messed up and create a list
#with run number as well so that we can change them below
wrong_names <- c("GCW_202304_TR_LysC_45cm", "GCW_202304_TR_LysA_20cm_8",
                 "GWI_202304_UP_LysA_20cm", "GWI_202304_UP_LysA_20cm")
wrong_nums <- c(20, 16, 46, 44)
correct_names <- c("GCW_202304_TR_LysB_45cm", "GCW_202304_TR_LysA_20cm",
                  "GWI_202304_UP_LysA_10cm", "GWI_202304_UP_LysA_10cm")

#can't determine from metadata - for now unsure
remove_names <- c("GCW_202304_TR_LysA_20cm", "GCW_202304_TR_LysA_20cm",
                  "GCW_202304_TR_LysB_20cm_13", "GCW_202304_TR_LysB_20cm_13")
#couldn't tell which one this is from the metadata, no A_10cm which is what we thought
#marked on the sheet, need to check sample vials in freezer
#to see if we have a A_10cm from GCW_TR to be sure
remove_nums <- c(15, 13, 21, 19 )

#Set up file path for metadata
#downloaded metadata csv - downloaded from Google drive as csv for this year
#https://docs.google.com/spreadsheets/d/1HCANO_q6y17x0RUXVzID09hVal-RfwWc/edit?usp=sharing&ouid=10899
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2023.csv"

```

```
##Setup
```

```
##Read in metadata and create similar sample IDs for matching to samples
```

## 0.1 Import Data & Clean


## 0.2 Assessing standard Curves

```
#Pull out standards data
```

```
## Assess Standard Curves
```


```
#Plot standards data
```

```
## Assess Standard Curves
```




```
## 'geom_smooth()' using formula = 'y ~ x'
```

## NH3 Standard Curve



```
## 'geom_smooth()' using formula = 'y ~ x'
```

## PO4 Standard Curve



```
## [1] "NOx Curve r2 GOOD - PROCEED"
```

```
## [1] "NH3 Curve r2 is below cutoff! - REASSESS"
```

```
## [1] "PO4 Curve r2 is below cutoff! - REASSESS"
```

```
## [1] "QAQC log file exists and has been read into the code."
```




Table 1: Average Slope by Analyte

Test	avg_slope
NH3	1.899
NOx	1.261
PO4	2.353

### 0.3 Dilution Corrections - ensure the latest dilution is kept

```
## Dilution Corrections
```

```
## Duplicated samples: GCW_202305_UP_LysA_20cm, GCW_202305_UP_LysB_10cm, GCW_202305_UP_LysB_20cm, GCW_202305_UP_LysB_30cm
```

```
## Duplicated samples with NO dilution present (possible input error): GCW_202305_SW_A, GCW_202305_SW_B
```

### 0.4 Temporary Chunk

### 0.5 Performance Check

```
## [1] "NOx pe Check has a % Difference <10% - PROCEED"
```

```
## Run mean = 0.712214
```

```
## Expected = 0.706
```

```
## [1] "NH3 pe Check has a % Difference >10% - REASSESS"
```

```
## Run mean = 1.437577
```

```
## Expected = 0.948
```


```
## [1] "PO4 pe Check has a % Difference <10% - PROCEED"
```

```
## Run mean = 0.893199
```

```
## Expected = 0.818
```

```
#Check NOx Reduction Efficiency
```

```
## Assess Reduction Efficiency
```



```
## [1] "Mean NOx Reduction Efficiency <95% - REASSESS"
```

```
## [1] 93.2968
```




## 0.6 Analyze the Check Standards

```
## Analyze Check Standards
```

```
## [1] "NOx Check Standard RSD within Range - PROCEED"
```

```
## [1] "NH3 CHECK STANDARD RSD TOO HIGH - REASSESS"
```

```
## [1] "PO4 CHECK STANDARD RSD TOO HIGH - REASSESS"
```



```
## [1] ">60% of NOx Check Standards are within range of expected concentration - PROCEED"
```

```
## [1] ">60% of NH3 Check Standards are within range of expected concentration - PROCEED"
```

```
## [1] ">60% of PO4 Check Standards are within range of expected concentration - PROCEED"
```

## 0.7 Analyze Blanks

```
## Assess Blanks
```

```
## [1] ">60% of NOx Blank concentrations are lower than the lower 25% quartile of samples - PROCEED"
```

```
## [1] ">60% of NH3 Blank concentrations are lower than the lower 25% quartile of samples - PROCEED"
```

```
## [1] ">60% of PO4 Blank concentrations are lower than the lower 25% quartile of samples- PROCEED"
```




Table 2: Mean Concentration of Blanks

Test	Blank_Mean_Conc
NOx	0.0017
NH3	-0.0171
PO4	-0.0221

## 0.8 Analyze Duplicates

```
## Analyze Duplicates
```


```
## Warning: Returning more (or less) than 1 row per 'summarise()' group was deprecated in
## dplyr 1.1.0.
## i Please use 'reframe()' instead.
## i When switching from 'summarise()' to 'reframe()', remember that 'reframe()'
## always returns an ungrouped data frame and adjust accordingly.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
## [1] "<60% of NOx Duplicates have a CV <10% - REASSESS"
```

```
## [1] "<60% of NH3 Duplicates have a CV <10% - REASSESS"
```

```
## [1] "<60% of PO4 Duplicates have a CV <10% - REASSESS"
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```




## 0.9 Spikes

## [1] ">60% of Spikes have a CV <50% - PROCEED"

## [1] ">60% of Spikes have a CV <50% - PROCEED"

## [1] ">60% of Spikes have a CV <50% - PROCEED"



## 0.10 Matrix Effects

```
## [1] "NO NOx Matrix Effect, PROCEED"
```

```
## [1] "NO NH3 Matrix Effect, PROCEED"
```

```
## [1] "NO PO4 Matrix Effect, PROCEED"
```

## 0.11 Unit Converted Data Column Added (mg/L to uM )

## 0.12 Sample Flagging - Within range of standard curve

```
## Sample Flagging
```

## 0.13 Pull out sample id information

```
## Sample Processing
```

```
## Warning: Expected 5 pieces. Missing pieces filled with 'NA' in 32 rows [19, 20, 21, 46,  
## 47, 48, 83, 84, 108, 109, 110, 129, 130, 131, 150, 151, 152, 177, 178, 179,  
## ...].
```

## 0.14 Pulling Out Rhizon Samples

## 0.15 Check to see if samples run match metadata & merge info


```
## Check Sample IDs with Metadata
```


```
## Some sample IDs are missing from metadata.
```


```
## [1] "GCW_202305_UP_LysA_20cm" "GCW_202305_UP_LysB_10cm"  
## [3] "GCW_202305_UP_LysB_20cm" "MSM_202305_UP_LysC_45cm"  
## [5] "MSM_202305_TR_LysC_45cm" "MSM_202305_WC_LysA_45cm"  
## [7] "SWH_202305_UPCON_LysA_10cm" "SWH_202305_UPCON_LysA_20cm"  
## [9] "SWH_202305_UPCON_LysA_45cm" "SWH_202305_UPCON_LysB_10cm"  
## [11] "SWH_202305_UPCON_LysB_20cm" "SWH_202305_UPCON_LysB_45cm"  
## [13] "SWH_202305_UPCON_LysC_20cm" "SWH_202305_UPCON_LysC_45cm"  
## [15] "SWH_202305_UP_LysB_10cm" "SWH_202305_WC_LysA_10cm"  
## [17] "SWH_202305_WC_LysA_20cm" "SWH_202305_WC_LysA_45cm"  
## [19] "SWH_202305_WC_LysB_10cm" "SWH_202305_WC_LysB_20cm"  
## [21] "GWI_202305_TR_LysB_45cm"
```

```
## Warning: Expected 5 pieces. Missing pieces filled with 'NA' in 32 rows [19, 20, 21, 46,  
## 47, 48, 83, 84, 108, 109, 110, 129, 130, 131, 150, 151, 152, 177, 178, 179,  
## ...].
```

## ## Visualize Data







## 0.17 Export Processed Data

#end