

Dionex_COMPASS_June2023

Stephanie J. Wilson

2023-06-23

Daily Set up

Read in and Format the raw data - change wd & file names

```
#setwd("S:/Biogeochemistry/People/Wilson (Steph)/Data/Dionex/Raw Data Files")

# SULFATE DATA:

## Read in raw data file from Dionex - copied and saved as a txt
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202304_S04.txt",sep='\t' , header=T, skip=3)
head(Sdat)
```

```
##      X      X.1      X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank      Unknown   n.a.      n.a.      n.a.      n.a.      n.a.
## 2 2 Lab Blank      Unknown   n.a.      n.a.      n.a.      n.a.      n.a.
## 3 3 Lab Blank      Unknown   n.a.      n.a.      n.a.      n.a.      n.a.
## 4 4 Lab Blank      Unknown   n.a.      n.a.      n.a.      n.a.      n.a.
## 5 5 Standard 1 Calibration Standard 4.680  0.4466  6.56  0.6608  2.79
## 6 6 Standard 2 Calibration Standard 4.677  0.9155  7.16  1.3545  5.57
##      IC.S04.5
## 1      n.a.
## 2      n.a.
## 3      n.a.
## 4      n.a.
## 5      M
## 6      M
```

```
## Only keep the columns that we need
Sdat <- Sdat[ ,c(2,5)] # dont need this here
head(Sdat)
```

```
##      X.1 IC.S04.1
## 1 Lab Blank   n.a.
## 2 Lab Blank   n.a.
## 3 Lab Blank   n.a.
## 4 Lab Blank   n.a.
## 5 Standard 1 0.4466
## 6 Standard 2 0.9155
```

```
## Name the columns correctly
colnames(Sdat) <- c( "Sample_ID", "SO4_ppm")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$SO4_ppm <- as.numeric(Sdat$SO4_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID SO4_ppm
## 1 Lab Blank      NA
## 2 Lab Blank      NA
## 3 Lab Blank      NA
## 4 Lab Blank      NA
## 5 Standard 1  0.4466
## 6 Standard 2  0.9155
```

```
#Chloride data
## Read in raw data file from Dionex - copied and saved as a txt
Cldat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202304_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X      X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard 3.710  6.0596  92.44  9.3112  42.63
## 6 6 Standard 2 Calibration Standard 3.710  11.3033  91.83  17.3688  78.65
##      IC.Cl.5
## 1      n.a.
## 2      n.a.
## 3      n.a.
## 4      n.a.
## 5      M
## 6      M
```

```
## Only keep the columns that we need
Cldat <- Cldat[ ,c(2,5)]
head(Cldat)
```

```
##      X.1 IC.Cl.1
## 1 Lab Blank  n.a.
## 2 Lab Blank  n.a.
## 3 Lab Blank  n.a.
## 4 Lab Blank  n.a.
## 5 Standard 1  6.0596
## 6 Standard 2 11.3033
```

```
## Name the columns correctly
colnames(Cldat) <- c("Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##   Sample_ID Cl_ppm
## 1 Lab Blank    NA
## 2 Lab Blank    NA
## 3 Lab Blank    NA
## 4 Lab Blank    NA
## 5 Standard 1  6.0596
## 6 Standard 2 11.3033
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##   Sample_ID S04_ppm Cl_ppm
## 1          2556.7130 21806.2530
## 2          2556.7130  6022.9406
## 3          2556.7130    0.0017
## 4          2556.7130  5600.4666
## 5          2556.7130         NA
## 6          597.3037 21806.2530
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##           Sample_ID S04_ppm Cl_ppm
## 26 1_GCW_202304_UP_LysA_20cm  8.2441  0.0017
## 27 10_GCW_202304_TR_LysB_20cm 24.9499  73.6636
## 28 10_GCW_202304_TR_LysB_20cm_dup 25.0390  73.6625
## 29 100_MSM_202304_PPR_TR_4 741.9526 5993.3481
## 30 100_MSM_202304_PPR_TR_4 741.9526         NA
## 31 100_MSM_202304_PPR_TR_4 741.9526         NA
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
```

```
## If there is an n.a. in the dataframe that is because there was no peak, which
# would mean there was no Cl or SO4 and we want to know that so make all n.a.'s into zeros
all_dat[is.na(all_dat)] <- 0
head(all_dat)
```

```
##           Sample_ID S04_ppm Cl_ppm
```

```
## 26      1_GCW_202304_UP_LysA_20cm  8.2441    0.0017
## 27     10_GCW_202304_TR_LysB_20cm 24.9499    73.6636
## 28 10_GCW_202304_TR_LysB_20cm_dup 25.0390    73.6625
## 29      100_MSM_202304_PPR_TR_4 741.9526 5993.3481
## 30      100_MSM_202304_PPR_TR_4 741.9526    0.0000
## 31      100_MSM_202304_PPR_TR_4 741.9526    0.0000
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grepl("Standard", all_dat$Sample_ID),]
stds <- stds[-c(3, 8, 13, 18, 23),] #this is if you need to remove one for any reason
head(stds)
```

```
##      Sample_ID S04_ppm Cl_ppm
## 638 Standard 1  0.6281 6.4219
## 639 Standard 1  0.6281 6.0596
## 641 Standard 1  0.6281 5.9994
## 642 Standard 1  0.6281 6.2815
## 643 Standard 1  0.4466 6.4219
## 644 Standard 1  0.4466 6.0596
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(S04_ppm), sd = sd(S04_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean    sd    cv flag
##   <fct>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1  5.40 10.1  186. NO, rerun
## 2 Standard 2  1.01  0.137 13.6 NO, rerun
## 3 Standard 3  2.20  0.284 12.9 NO, rerun
## 4 Standard 4 10.7  1.25  11.7 NO, rerun
## 5 Standard 5 21.7  1.91   8.80 NO, rerun
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean    sd    cv flag
##   <fct>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1  6.19 0.174  2.82 YES
## 2 Standard 2 11.7  0.231  1.99 YES
## 3 Standard 3 22.9  0.533  2.33 YES
## 4 Standard 4 103.  2.59  2.51 YES
## 5 Standard 5 206.  4.47  2.17 YES
```

Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##                Sample_ID  S04_ppm  Cl_ppm
## 26      1_GCW_202304_UP_LysA_20cm    8.2441    0.0017
## 27     10_GCW_202304_TR_LysB_20cm   24.9499   73.6636
## 28 10_GCW_202304_TR_LysB_20cm_dup   25.0390   73.6625
## 29      100_MSM_202304_PPR_TR_4  741.9526 5993.3481
## 30      100_MSM_202304_PPR_TR_4  741.9526    0.0000
## 31      100_MSM_202304_PPR_TR_4  741.9526    0.0000
```

```
# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06       #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)
```

```
##                Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 26      1_GCW_202304_UP_LysA_20cm    8.2441    0.0017  0.2571460 4.795487e-05
## 27     10_GCW_202304_TR_LysB_20cm   24.9499   73.6636  0.7782252 2.077958e+00
## 28 10_GCW_202304_TR_LysB_20cm_dup   25.0390   73.6625  0.7810044 2.077927e+00
## 29      100_MSM_202304_PPR_TR_4  741.9526 5993.3481 23.1426263 1.690648e+02
## 30      100_MSM_202304_PPR_TR_4  741.9526    0.0000 23.1426263 0.000000e+00
## 31      100_MSM_202304_PPR_TR_4  741.9526    0.0000 23.1426263 0.000000e+00
##      salinity
## 26 2.907190e-05
## 27 1.331361e-01
## 28 1.331341e-01
## 29 1.083001e+01
## 30 2.600000e-05
## 31 2.600000e-05
```

Pull out dups and check with percent difference

```
#Show me the data that we have from the calculations  
head(sampledat)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM  
## 26      1_GCW_202304_UP_LysA_20cm  8.2441  0.0017  0.2571460 4.795487e-05  
## 27      10_GCW_202304_TR_LysB_20cm 24.9499  73.6636  0.7782252 2.077958e+00  
## 28 10_GCW_202304_TR_LysB_20cm_dup 25.0390  73.6625  0.7810044 2.077927e+00  
## 29      100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 1.690648e+02  
## 30      100_MSM_202304_PPR_TR_4 741.9526  0.0000 23.1426263 0.000000e+00  
## 31      100_MSM_202304_PPR_TR_4 741.9526  0.0000 23.1426263 0.000000e+00  
##      salinity  
## 26 2.907190e-05  
## 27 1.331361e-01  
## 28 1.331341e-01  
## 29 1.083001e+01  
## 30 2.600000e-05  
## 31 2.600000e-05
```

```
#pull out any rows that have "dup" in the SampleID column  
dups <- sampledat %>%  
  filter(str_detect(Sample_ID, "dup")) #have to change this to match data  
head(dups)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM  
## 1 10_GCW_202304_TR_LysB_20cm_dup 25.0390  73.6625  0.7810044  2.077927  
## 2   100_MSM_202304_PPR_TR_4_dup 739.1719 6000.1002 23.0558921 169.255295  
## 3   100_MSM_202304_PPR_TR_4_dup 739.1719  0.0000 23.0558921  0.000000  
## 4   100_MSM_202304_PPR_TR_4_dup 739.1719  0.0000 23.0558921  0.000000  
## 5   100_MSM_202304_PPR_TR_4_dup 11.4443 6000.1002  0.3569651 169.255295  
## 6   100_MSM_202304_PPR_TR_4_dup 11.4443  0.0000  0.3569651  0.000000  
##      salinity  
## 1  0.1331341  
## 2 10.8422071  
## 3  0.0000260  
## 4  0.0000260  
## 5 10.8422071  
## 6  0.0000260
```

```
#remove these from sample dataframe in a new dataframe  
sampledat2 <- sampledat %>%  
  filter(!str_detect(Sample_ID, "dup")) %>%  
  filter(!str_detect(Sample_ID, "spk"))  
head(sampledat2)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM  
## 1  1_GCW_202304_UP_LysA_20cm  8.2441  0.0017  0.2571460 4.795487e-05  
## 2 10_GCW_202304_TR_LysB_20cm 24.9499  73.6636  0.7782252 2.077958e+00  
## 3   100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 1.690648e+02  
## 4   100_MSM_202304_PPR_TR_4 741.9526  0.0000 23.1426263 0.000000e+00
```

```
## 5      100_MSM_202304_PPR_TR_4 741.9526      0.0000 23.1426263 0.000000e+00
## 6      100_MSM_202304_PPR_TR_4   9.2977 5993.3481  0.2900094 1.690648e+02
##          salinity
## 1 2.907190e-05
## 2 1.331361e-01
## 3 1.083001e+01
## 4 2.600000e-05
## 5 2.600000e-05
## 6 1.083001e+01
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[ , -c(2,3)]
colnames(dups) <- c('Sample_ID', 'S04_mM_dup', 'Cl_mM_dup', "salinity_dup")
head(dups)
```

```
##              Sample_ID S04_mM_dup Cl_mM_dup salinity_dup
## 1 10_GCW_202304_TR_LysB_20cm 0.7810044  2.077927  0.1331341
## 2      100_MSM_202304_PPR_TR_4 23.0558921 169.255295  10.8422071
## 3      100_MSM_202304_PPR_TR_4 23.0558921  0.000000  0.0000260
## 4      100_MSM_202304_PPR_TR_4 23.0558921  0.000000  0.0000260
## 5      100_MSM_202304_PPR_TR_4  0.3569651 169.255295  10.8422071
## 6      100_MSM_202304_PPR_TR_4  0.3569651  0.000000  0.0000260
```

```
#put it back together with the old data set and look for duplicates
QAdups <- merge(sampldat2, dups)
head(QAdups)
```

```
##              Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
## 1 10_GCW_202304_TR_LysB_20cm 24.9499  73.6636 0.7782252  2.077958
## 2      100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 3      100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 4      100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 5      100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 6      100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
##          salinity S04_mM_dup Cl_mM_dup salinity_dup
## 1 0.1331361 0.7810044  2.077927  0.1331341
## 2 10.8300060 23.0558921 169.255295  10.8422071
## 3 10.8300060 23.0558921  0.000000  0.0000260
## 4 10.8300060 23.0558921  0.000000  0.0000260
## 5 10.8300060  0.3569651 169.255295  10.8422071
## 6 10.8300060  0.3569651  0.000000  0.0000260
```

```
QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*100
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')

QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')

head(QAdups)
```

```
##              Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM
```

```
## 1 10_GCW_202304_TR_LysB_20cm 24.9499 73.6636 0.7782252 2.077958
## 2 100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 3 100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 4 100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 5 100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## 6 100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 169.064827
## salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 0.1331361 0.7810044 2.077927 0.1331341 0.3564791 YES
## 2 10.8300060 23.0558921 169.255295 10.8422071 0.3754850 YES
## 3 10.8300060 23.0558921 0.000000 0.0000260 0.3754850 YES
## 4 10.8300060 23.0558921 0.000000 0.0000260 0.3754850 YES
## 5 10.8300060 0.3569651 169.255295 10.8422071 193.9238932 NO, rerun
## 6 10.8300060 0.3569651 0.000000 0.0000260 193.9238932 NO, rerun
## Cl_dups_chk Cl_dups_flag
## 1 1.493286e-03 YES
## 2 1.125965e-01 YES
## 3 2.000000e+02 NO, rerun
## 4 2.000000e+02 NO, rerun
## 5 1.125965e-01 YES
## 6 2.000000e+02 NO, rerun
```

```
QAdups <- QAdups %>%
  mutate(row_number = row_number()) %>%
  filter(!Sample_ID == "100_MSM_202304_PPR_TR_4")
```


#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this

```
Sdupsbar <- ggplot(data = QAdups, aes(x = row_number, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
Cldupsbar <- ggplot(data = QAdups, aes(x = row_number, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```

```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups1 <- Perc_dups1[-c(3),]
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         8  YES         8      8      100      100
```

Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)

```

```

##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 10_GCW_202304_TR_LysB_20cm 24.9499  73.6636 0.7782252 2.077958
## 2 110_MSM_202304_PPR_WC_2 728.5031 13337.4514 22.7231160 376.232762
## 3 20_GCW_202304_WC_LysB_45cm 440.9736 7417.9642 13.7546351 209.251458
## 4 30_MSM_202304_UP_LysB_45cm 222.6503 2656.9852 6.9448004 74.950217
## 5 40_MSM_202304_TR_LysC_10cm 579.3436 5440.5188 18.0706051 153.470206
## 6 60_GWI_202304_WC_LysB_20cm 2050.6464 17258.5113 63.9627698 486.840939
##      salinity S04_mM_dup  Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 0.1331361 0.7810044 2.077927 0.1331341 0.3564791 YES
## 2 24.1008007 22.3924579 373.019901 23.8949908 1.4658270 YES
## 3 13.4042873 13.5096163 205.038812 13.1344330 1.7973625 YES
## 4 4.8011983 6.7942452 73.354116 4.6989549 2.1916400 YES
## 5 9.8310435 18.3222333 151.588496 9.7105046 1.3828446 YES
## 6 31.1861559 59.0870056 454.059289 29.0862241 7.9248648 YES
##      Cl_dups_chk Cl_dups_flag row_number S04_dups_cv S04_dups_cv_flag
## 1 0.001493286 YES 1 0.2520688 YES
## 2 0.857617337 YES 83 1.0364962 YES
## 3 2.033668796 YES 84 1.2709272 YES
## 4 2.152467862 YES 85 1.5497235 YES
## 5 1.233670395 YES 86 0.9778188 YES
## 6 6.968145871 YES 87 5.6037257 YES

```

```

QAdups <- QAdups %>%
  mutate(row_number = row_number()) %>%
  filter(!Sample_ID == "100_MSM_202304_PPR_TR_4")

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = row_number, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = row_number, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +


```

```

theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
      color = "black", size=1)

ggarrange(Sdupsbar, Cl_dupsbar, ncol=2, nrow=1)

```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups1 <- Perc_dups1[,-c(3),]
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         8  YES         8      8      100      100

```

Pull out spikes and check

#Show me the data that we have from the calculations

```
head(sampledat)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 26 1_GCW_202304_UP_LysA_20cm  8.2441  0.0017  0.2571460 4.795487e-05
## 27 10_GCW_202304_TR_LysB_20cm 24.9499 73.6636 0.7782252 2.077958e+00
## 28 10_GCW_202304_TR_LysB_20cm_dup 25.0390 73.6625 0.7810044 2.077927e+00
## 29 100_MSM_202304_PPR_TR_4 741.9526 5993.3481 23.1426263 1.690648e+02
## 30 100_MSM_202304_PPR_TR_4 741.9526 0.0000 23.1426263 0.000000e+00
## 31 100_MSM_202304_PPR_TR_4 741.9526 0.0000 23.1426263 0.000000e+00
##           salinity
## 26 2.907190e-05
## 27 1.331361e-01
## 28 1.331341e-01
## 29 1.083001e+01
## 30 2.600000e-05
## 31 2.600000e-05
```

#pull out any rows that have "spk" in the SampleID column

```
spks <- sampledat %>%
```

```
  filter(str_detect(Sample_ID, "spk"))
```

#have to change this to match data

```
head(spks)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 101_MSM_202304_PPR_TR_5_spk  5.7280  0.0000  0.178665  0.000000
## 2 101_MSM_202304_PPR_TR_5_spk  5.7280 6351.8704  0.178665 179.178291
## 3 101_MSM_202304_PPR_TR_5_spk 428.2802  0.0000 13.358709  0.000000
## 4 101_MSM_202304_PPR_TR_5_spk 428.2802 6351.8704 13.358709 179.178291
## 5 11_GCW_202304_TR_LysB_45cm_spk 80.9716  98.3241  2.525627  2.773599
## 6 111_MSM_202304_PPR_WC_3_spk 1822.7867 12621.9929 56.855480 356.050575
##           salinity
## 1 0.0000260
## 2 11.4778558
## 3 0.0000260
## 4 11.4778558
## 5 0.1776976
## 6 22.8079672
```

#remove the dup from these IDs so we will have duplicate sample names

```
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
```

```
spks <- spks[ ,-c(2,3, 5,6)]
```

```
colnames(spks) <- c('Sample_ID', 'S04_mM_spk')
```

```
head(spks)
```

```
##           Sample_ID S04_mM_spk
## 1 101_MSM_202304_PPR_TR_5 0.178665
## 2 101_MSM_202304_PPR_TR_5 0.178665
## 3 101_MSM_202304_PPR_TR_5 13.358709
## 4 101_MSM_202304_PPR_TR_5 13.358709
## 5 11_GCW_202304_TR_LysB_45cm 2.525627
## 6 111_MSM_202304_PPR_WC_3 56.855480
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledats, spks)
head(QAspks)
```

```
##           Sample_ID  S04_ppm Cl_ppm  S04_mM Cl_mM salinity S04_mM_spk
## 1 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 2 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 3 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 4 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 5 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665
## 6 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvol <- 10              # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvol      # mmoles of S04
head(QAspks)
```

```
##           Sample_ID  S04_ppm Cl_ppm  S04_mM Cl_mM salinity S04_mM_spk
## 1 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 2 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 3 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 4 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 5 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665
## 6 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665
##      S04_spk_Conc
## 1 7.797879e-05
## 2 7.797879e-05
## 3 7.797879e-05
## 4 7.797879e-05
## 5 7.797879e-05
## 6 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 50, QAspks$Dilution)
```

```

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 1462, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 1501, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

```

```

##           Sample_ID  S04_ppm Cl_ppm  S04_mM Cl_mM salinity S04_mM_spk
## 1 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 2 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 3 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 4 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 5 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665
## 6 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665
##      S04_spk_Conc Dilution SampleVol
## 1 7.797879e-05      1      1e-06
## 2 7.797879e-05      1      1e-06
## 3 7.797879e-05      1      1e-06
## 4 7.797879e-05      1      1e-06
## 5 7.797879e-05      1      1e-06
## 6 7.797879e-05      1      1e-06

```

```

#gives us the total S04 in the sample in mmoles
QAspks$S04_Total_unspkd <- (QAspks$S04_mM/QAspks$Dilution)*(QAspks$SampleVol)

##total S04 in spiked sample in mmoles
QAspks$S04_Total_spkd <- (QAspks$S04_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)

QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO', rerun=TRUE)

head(QAspks)

```

```

##           Sample_ID  S04_ppm Cl_ppm  S04_mM Cl_mM salinity S04_mM_spk
## 1 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 2 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05  0.178665
## 3 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 4 101_MSM_202304_PPR_TR_5  0.0000      0 0.00000      0 2.6e-05 13.358709
## 5 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665
## 6 101_MSM_202304_PPR_TR_5 783.8743      0 24.45023      0 2.6e-05  0.178665


```

```
##   S04_spk_Conc Dilution SampleVol S04_Total_unspkd S04_Total_spkd
## 1 7.797879e-05      1      1e-06      0.000000e+00      1.965315e-06
## 2 7.797879e-05      1      1e-06      0.000000e+00      1.965315e-06
## 3 7.797879e-05      1      1e-06      0.000000e+00      1.469458e-04
## 4 7.797879e-05      1      1e-06      0.000000e+00      1.469458e-04
## 5 7.797879e-05      1      1e-06      2.445023e-05      1.965315e-06
## 6 7.797879e-05      1      1e-06      2.445023e-05      1.965315e-06
##   S04_expctd_spkd spk_recovery S04_spks_flag
## 1 7.797879e-05      2.520320      NO, rerun
## 2 7.797879e-05      2.520320      NO, rerun
## 3 7.797879e-05     188.443288      NO, rerun
## 4 7.797879e-05     188.443288      NO, rerun
## 5 1.024290e-04      1.918709      NO, rerun
## 6 1.024290e-04      1.918709      NO, rerun
```

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

spksbar



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 2 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 NO, rerun      43    44    97.7
## 2 YES            1    44     2.27
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_ ', fixed=TRUE)))
```

```
## Warning in rbind(c("1", "GCW", "202304", "UP", "LysA", "20cm"), c("10", :
## number of columns of result is not a multiple of vector length (arg 1)
```

```
colnames(IDs) <- c("Analysis_No", "Site", "Date", "Zone", "Replicate", "Depth") #, "RHZ", "RHZ_Rep")
head(IDs)
```

```
##   Analysis_No Site   Date Zone Replicate Depth  NA
## 1           1  GCW 202304  UP      LysA  20cm   1
## 2          10  GCW 202304  TR      LysB  20cm  10
## 3         100  MSM 202304  PPR       TR    4 100
## 4         100  MSM 202304  PPR       TR    4 100
## 5         100  MSM 202304  PPR       TR    4 100
## 6         100  MSM 202304  PPR       TR    4 100
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampledat2)
head(alldat)
```

```
##   Analysis_No Site   Date Zone Replicate Depth  NA      Sample_ID
## 1           1  GCW 202304  UP      LysA  20cm   1  1_GCW_202304_UP_LysA_20cm
## 2          10  GCW 202304  TR      LysB  20cm  10 10_GCW_202304_TR_LysB_20cm
## 3         100  MSM 202304  PPR       TR    4 100  100_MSM_202304_PPR_TR_4
## 4         100  MSM 202304  PPR       TR    4 100  100_MSM_202304_PPR_TR_4
## 5         100  MSM 202304  PPR       TR    4 100  100_MSM_202304_PPR_TR_4
## 6         100  MSM 202304  PPR       TR    4 100  100_MSM_202304_PPR_TR_4
##   S04_ppm  Cl_ppm    S04_mM    Cl_mM    salinity
## 1   8.2441   0.0017 0.2571460 4.795487e-05 2.907190e-05
## 2  24.9499  73.6636 0.7782252 2.077958e+00 1.331361e-01
## 3 741.9526 5993.3481 23.1426263 1.690648e+02 1.083001e+01
## 4 741.9526   0.0000 23.1426263 0.000000e+00 2.600000e-05
## 5 741.9526   0.0000 23.1426263 0.000000e+00 2.600000e-05
## 6   9.2977 5993.3481 0.2900094 1.690648e+02 1.083001e+01
```


Make final dataframe with IDs

```
#setwd("S:/Biogeochemistry/People/Wilson (Steph)/Data/Dionex/Final Data Files")      #Change wd  
write.csv(alldat, file="Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_S04_202304.csv")  #C
```

END