

# Dionex\_COMPASS\_June2023

Stephanie J. Wilson

2023-06-23

## Daily Set up

Read in and Format the raw data - change wd & file names

```
# SULFATE DATA:
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202310_S04.txt", sep='\t', header=T, skip=3)
head(Sdat)
```

```
##      X      X.1      X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank      Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
## 2 2 Lab Blank      Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
## 3 3 Lab Blank      Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
## 4 4 Lab Blank      Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
## 5 5 Standard 1 Calibration Standard 5.097  0.4573  6.48  0.5739  3.63
## 6 6 Standard 2 Calibration Standard 5.150  0.9857  6.85  1.2371  7.50
```

```
## Only keep the columns that we need
```

```
Sdat <- Sdat[,c(2,5)] # dont need this here
head(Sdat)
```

```
##      X.1 IC.S04.1
## 1 Lab Blank      n.a.
## 2 Lab Blank      n.a.
## 3 Lab Blank      n.a.
## 4 Lab Blank      n.a.
## 5 Standard 1 0.4573
## 6 Standard 2 0.9857
```

```
## Name the columns correctly
```

```
colnames(Sdat) <- c("Sample_ID", "S04_ppm")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID S04_ppm
## 1 Lab Blank      NA
## 2 Lab Blank      NA
## 3 Lab Blank      NA
## 4 Lab Blank      NA
## 5 Standard 1  0.4573
## 6 Standard 2  0.9857
```

```
#Chloride data
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Cldat <- read.table("raw Data/COMPASS_Synoptic_CB_MonMon_202310_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X      X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard 4.040  4.9368  93.52  8.2828  64.04
## 6 6 Standard 2 Calibration Standard 4.030  10.0322  93.15  16.8317  127.41
```

```
## Only keep the columns that we need
```

```
Cldat <- Cldat[,c(2,5)]
head(Cldat)
```

```
##      X.1 IC.Cl.1
## 1 Lab Blank  n.a.
## 2 Lab Blank  n.a.
## 3 Lab Blank  n.a.
## 4 Lab Blank  n.a.
## 5 Standard 1  4.9368
## 6 Standard 2 10.0322
```

```
## Name the columns correctly
```

```
colnames(Cldat) <- c("Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##      Sample_ID Cl_ppm
## 1 Lab Blank      NA
## 2 Lab Blank      NA
## 3 Lab Blank      NA
## 4 Lab Blank      NA
## 5 Standard 1  4.9368
## 6 Standard 2 10.0322
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##      Sample_ID   S04_ppm   Cl_ppm
## 1           2322.4204 33352.2387
## 2           2322.4204  5714.1618
## 3           2322.4204    0.0061
## 4           2322.4204  5966.0950
## 5           2322.4204      NA
## 6           592.7974 33352.2387
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##              Sample_ID   S04_ppm   Cl_ppm
## 26      1_MSM_202310_UP_LysA_20cm 156.4865 2514.586
## 27     10_MSM_202310_TR_LysA_45cm 543.5656 4326.371
## 28 10_MSM_202310_TR_LysA_45cm_dup 533.3889 4242.405
## 29    100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979
## 30 100_GWI_202310_TR_LysC_45cm_dup 1151.9205 9252.682
## 31    101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
```

```
## If there is an n.a. in the dataframe that is because there was no peak, which
# would mean there was no Cl or S04 and we want to know that so make all n.a.'s into zeros
all_dat[is.na(all_dat)] <- 0
head(all_dat)
```

```
##              Sample_ID   S04_ppm   Cl_ppm
## 26      1_MSM_202310_UP_LysA_20cm 156.4865 2514.586
## 27     10_MSM_202310_TR_LysA_45cm 543.5656 4326.371
## 28 10_MSM_202310_TR_LysA_45cm_dup 533.3889 4242.405
## 29    100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979
## 30 100_GWI_202310_TR_LysC_45cm_dup 1151.9205 9252.682
## 31    101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grepl("Standard", all_dat$Sample_ID),]
#stds <- stds[-c(80),] #this is if you need to remove one for any reason
head(stds)
```

```
##      Sample_ID S04_ppm Cl_ppm
## 691 Standard 1  0.5460 5.6530
## 692 Standard 1  0.5460 4.9368
## 693 Standard 1  0.5460 5.1270
```

```
## 694 Standard 1 0.5460 5.6753
## 695 Standard 1 0.5460 5.5434
## 696 Standard 1 0.4573 5.6530
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(SO4_ppm), sd = sd(SO4_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID    mean      sd    cv flag
##   <fct>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1 0.516 0.0355 6.87 NO, rerun
## 2 Standard 2 1.09 0.0604 5.55 NO, rerun
## 3 Standard 3 2.14 0.107 5.03 NO, rerun
## 4 Standard 4 11.0 0.504 4.57 YES
## 5 Standard 5 21.4 0.977 4.55 YES
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 5, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID    mean      sd    cv flag
##   <fct>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1 5.39 0.306 5.67 NO, rerun
## 2 Standard 2 11.0 0.538 4.91 YES
## 3 Standard 3 21.4 1.05 4.92 YES
## 4 Standard 4 89.4 45.6 51.1 NO, rerun
## 5 Standard 5 216. 10.1 4.70 YES
```

## Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##           Sample_ID    SO4_ppm    Cl_ppm
## 26 1_MSM_202310_UP_LysA_20cm 156.4865 2514.586
## 27 10_MSM_202310_TR_LysA_45cm 543.5656 4326.371
## 28 10_MSM_202310_TR_LysA_45cm_dup 533.3889 4242.405
## 29 100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979
## 30 100_GWI_202310_TR_LysC_45cm_dup 1151.9205 9252.682
## 31 101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559
```

```

# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06      #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$SO4_mM <- (sampledat$SO4_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)

```

```

##              Sample_ID   SO4_ppm   Cl_ppm   SO4_mM   Cl_mM
## 26      1_MSM_202310_UP_LysA_20cm 156.4865 2514.586  4.881051  70.93331
## 27      10_MSM_202310_TR_LysA_45cm 543.5656 4326.371 16.954635 122.04149
## 28 10_MSM_202310_TR_LysA_45cm_dup 533.3889 4242.405 16.637208 119.67292
## 29      100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979 36.468936 266.14891
## 30 100_GWI_202310_TR_LysC_45cm_dup 1151.9205 9252.682 35.930147 261.00655
## 31      101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559 56.052024 411.01717
##      salinity
## 26 4.543882
## 27 7.817778
## 28 7.666052
## 29 17.049033
## 30 16.719623
## 31 26.329026

```

Pull out dups and check with percent difference

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##              Sample_ID   SO4_ppm   Cl_ppm   SO4_mM   Cl_mM
## 26      1_MSM_202310_UP_LysA_20cm 156.4865 2514.586  4.881051  70.93331
## 27      10_MSM_202310_TR_LysA_45cm 543.5656 4326.371 16.954635 122.04149
## 28 10_MSM_202310_TR_LysA_45cm_dup 533.3889 4242.405 16.637208 119.67292
## 29      100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979 36.468936 266.14891
## 30 100_GWI_202310_TR_LysC_45cm_dup 1151.9205 9252.682 35.930147 261.00655
## 31      101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559 56.052024 411.01717
##      salinity
## 26 4.543882
## 27 7.817778
## 28 7.666052

```

```
## 29 17.049033
## 30 16.719623
## 31 26.329026
```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledats %>%
  filter(str_detect(Sample_ID, "dup")) #have to change this to match data
head(dups)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 10_MSM_202310_TR_LysA_45cm_dup 533.3889 4242.405 16.63721 119.6729
## 2 100_GWI_202310_TR_LysC_45cm_dup 1151.9205 9252.682 35.93015 261.0065
## 3 110_GWI_202310_SW_B_dup 1924.4210 13794.211 60.02561 389.1174
## 4 120_MSM_202310_TR_RHZ_SF_Tree_3_dup 938.2737 8470.767 29.26618 238.9497
## 5 130_MSM_202310_WC_RHZ_SF_Collar_5_dup 1047.0756 9863.028 32.65988 278.2236
## 6 20_MSM_202310_WC_LysB_10cm_dup 872.9288 9105.025 27.22797 256.8413
## salinity
## 1 7.666052
## 2 16.719623
## 3 24.926166
## 4 15.306701
## 5 17.822518
## 6 16.452806
```

```
#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledats %>%
  filter(!str_detect(Sample_ID, "dup")) %>%
  filter(!str_detect(Sample_ID, "spk"))
head(sampledat2)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM  salinity
## 1 1_MSM_202310_UP_LysA_20cm 156.4865 2514.586 4.881051 70.93331 4.543882
## 2 10_MSM_202310_TR_LysA_45cm 543.5656 4326.371 16.954635 122.04149 7.817778
## 3 100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979 36.468936 266.14891 17.049033
## 4 101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559 56.052024 411.01717 26.329026
## 5 102_GWI_202310_WC_LysA_20cm 1448.6028 13511.788 45.184117 381.15059 24.415828
## 6 103_GWI_202310_WC_LysA_45cm 1413.1583 12517.120 44.078550 353.09223 22.618461
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[ ,-c(2,3)]
colnames(dups) <- c('Sample_ID', 'S04_mM_dup', 'Cl_mM_dup', "salinity_dup")
head(dups)
```

```
##           Sample_ID S04_mM_dup Cl_mM_dup salinity_dup
## 1 10_MSM_202310_TR_LysA_45cm 16.63721 119.6729 7.666052
## 2 100_GWI_202310_TR_LysC_45cm 35.93015 261.0065 16.719623
## 3 110_GWI_202310_SW_B 60.02561 389.1174 24.926166
## 4 120_MSM_202310_TR_RHZ_SF_Tree_3 29.26618 238.9497 15.306701
## 5 130_MSM_202310_WC_RHZ_SF_Collar_5 32.65988 278.2236 17.822518
## 6 20_MSM_202310_WC_LysB_10cm 27.22797 256.8413 16.452806
```

```
#put it back together with the old data set and look for duplicates
```

```
QAdups <- merge(sampledat2, dups)
```

```
head(QAdups)
```

```
##           Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 1      10_MSM_202310_TR_LysA_45cm 543.5656 4326.3707 16.95464 122.041487
## 2      100_GWI_202310_TR_LysC_45cm 1169.1941 9434.9790 36.46894 266.148914
## 3           110_GWI_202310_SW_B 1925.0477  108.1539 60.04516   3.050886
## 4      120_MSM_202310_TR_RHZ_SF_Tree_3 904.1450 8163.1000 28.20165 230.270804
## 5 130_MSM_202310_WC_RHZ_SF_Collar_5 1056.5185 9928.2328 32.95441 280.062984
## 6       20_MSM_202310_WC_LysB_10cm 888.9918 9296.0155 27.72900 262.228928
##      salinity S04_mM_dup Cl_mM_dup salinity_dup
## 1  7.8177779   16.63721  119.6729    7.666052
## 2 17.0490331   35.93015  261.0065   16.719623
## 3  0.1954601   60.02561  389.1174   24.926166
## 4 14.7507477   29.26618  238.9497   15.306701
## 5 17.9403427   32.65988  278.2236   17.822518
## 6 16.7979260   27.22797  256.8413   16.452806
```

```
QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*100
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')
```

```
QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')
```

```
head(QAdups)
```

```
##           Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 1      10_MSM_202310_TR_LysA_45cm 543.5656 4326.3707 16.95464 122.041487
## 2      100_GWI_202310_TR_LysC_45cm 1169.1941 9434.9790 36.46894 266.148914
## 3           110_GWI_202310_SW_B 1925.0477  108.1539 60.04516   3.050886
## 4      120_MSM_202310_TR_RHZ_SF_Tree_3 904.1450 8163.1000 28.20165 230.270804
## 5 130_MSM_202310_WC_RHZ_SF_Collar_5 1056.5185 9928.2328 32.95441 280.062984
## 6       20_MSM_202310_WC_LysB_10cm 888.9918 9296.0155 27.72900 262.228928
##      salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1  7.8177779   16.63721  119.6729    7.666052   1.88990343         YES
## 2 17.0490331   35.93015  261.0065   16.719623   1.48838838         YES
## 3  0.1954601   60.02561  389.1174   24.926166   0.03256034         YES
## 4 14.7507477   29.26618  238.9497   15.306701   3.70477134         YES
## 5 17.9403427   32.65988  278.2236   17.822518   0.89778727         YES
## 6 16.7979260   27.22797  256.8413   16.452806   1.82335118         YES
##      Cl_dups_chk Cl_dups_flag
## 1   1.9598039         YES
## 2   1.9509868         YES
## 3 196.8881871    NO, rerun
## 4   3.6992794         YES
## 5   0.6589252         YES
## 6   2.0758663         YES
```

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
```


```
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
```

```
scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
geom_bar(stat = 'identity') +
scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
```



```

Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows      Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         15 NO, rerun         4     15         100    26.66667
## 2  YES         15     YES         11     15         100    73.33333

```

## Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)

```

```

##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1  10_MSM_202310_TR_LysA_45cm  543.5656  4326.3707  16.95464  122.041487
## 2  100_GWI_202310_TR_LysC_45cm 1169.1941  9434.9790  36.46894  266.148914
## 3      110_GWI_202310_SW_B 1925.0477  108.1539  60.04516   3.050886
## 4  120_MSM_202310_TR_RHZ_SF_Tree_3  904.1450  8163.1000  28.20165  230.270804
## 5 130_MSM_202310_WC_RHZ_SF_Collar_5 1056.5185  9928.2328  32.95441  280.062984
## 6   20_MSM_202310_WC_LysB_10cm  888.9918  9296.0155  27.72900  262.228928
##   salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1  7.8177779  16.63721  119.6729   7.666052  1.88990343         YES
## 2 17.0490331  35.93015  261.0065  16.719623  1.48838838         YES
## 3  0.1954601  60.02561  389.1174   24.926166  0.03256034         YES
## 4 14.7507477  29.26618  238.9497   15.306701  3.70477134         YES
## 5 17.9403427  32.65988  278.2236   17.822518  0.89778727         YES
## 6 16.7979260  27.22797  256.8413   16.452806  1.82335118         YES
##   Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1  1.9598039         YES  1.33636353         YES
## 2  1.9509868         YES  1.05244952         YES
## 3 196.8881871    NO, rerun  0.02302364         YES
## 4  3.6992794         YES  2.61966894         YES
## 5  0.6589252         YES  0.63483147         YES


```

```
## 6      2.0758663      YES      1.28930398      YES
```

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar,ncol=2, nrow=1)
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
```

```

group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows      Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         15 NO, rerun         4     15         100    26.66667
## 2  YES         15    YES         11     15         100    73.33333

```

## Pull out spikes and check

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 26  1_MSM_202310_UP_LysA_20cm 156.4865 2514.586 4.881051 70.93331
## 27  10_MSM_202310_TR_LysA_45cm 543.5656 4326.371 16.954635 122.04149
## 28  10_MSM_202310_TR_LysA_45cm_dup 533.3889 4242.405 16.637208 119.67292
## 29  100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979 36.468936 266.14891
## 30  100_GWI_202310_TR_LysC_45cm_dup 1151.9205 9252.682 35.930147 261.00655
## 31  101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559 56.052024 411.01717
##      salinity
## 26  4.543882
## 27  7.817778
## 28  7.666052
## 29 17.049033
## 30 16.719623
## 31 26.329026

```

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)

```

```

##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1  101_GWI_202310_WC_LysA_10cm_spk 2123.4532 13985.840 66.23372 394.5230
## 2   11_MSM_202310_TR_LysB_10cm_spk  824.2115  6953.347 25.70841 196.1452
## 3   111_GWI_202310_SW_C_spk 2322.4204 13800.892 72.43981 389.3058
## 4 121_MSM_202310_TR_RHZ_SF_Tree_4_spk 1015.5222  8245.100 31.67568 232.5839
## 5   131_MSM_202307_WC_RHZ_LysA_spk 1420.6403  9776.295 44.31192 275.7770
## 6   21_MSM_202310_WC_LysB_20cm_spk  931.4789  9421.024 29.05424 265.7553
##      salinity
## 1 25.27244
## 2 12.56472

```

```
## 3 24.93824
## 4 14.89892
## 5 17.66579
## 6 17.02382
```

```
#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[ , -c(2,3, 5,6)]
colnames(spks) <- c('Sample_ID', 'SO4_mM_spk')
head(spks)
```

```
##              Sample_ID SO4_mM_spk
## 1 101_GWI_202310_WC_LysA_10cm 66.23372
## 2 11_MSM_202310_TR_LysB_10cm 25.70841
## 3 111_GWI_202310_SW_C 72.43981
## 4 121_MSM_202310_TR_RHZ_SF_Tree_4 31.67568
## 5 131_MSM_202307_WC_RHZ_LysA 44.31192
## 6 21_MSM_202310_WC_LysB_20cm 29.05424
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampldat, spks)
head(QAspks)
```

```
##              Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
## 1 101_GWI_202310_WC_LysA_10cm 1797.0279 14570.5587 56.052024 411.01717
## 2 11_MSM_202310_TR_LysB_10cm 744.7543 7094.4554 23.230016 200.12568
## 3 111_GWI_202310_SW_C 1946.6882 13991.2256 60.720156 394.67491
## 4 121_MSM_202310_TR_RHZ_SF_Tree_4 894.6568 8047.9702 27.905702 227.02314
## 5 21_MSM_202310_WC_LysB_20cm 747.1359 9510.8725 23.304301 268.28977
## 6 31_GCW_202310_TR_LysB_20cm 50.6452 776.0658 1.579701 21.89184
## salinity SO4_mM_spk
## 1 26.329026 66.233724
## 2 12.819707 25.708406
## 3 25.282171 72.439813
## 4 14.542708 31.675677
## 5 17.186173 29.054239
## 6 1.402377 4.434775
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw) # in mM
spkvol <- 10 # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$SO4_spk_Conc <- (spkconc)*spkvol # mmoles of SO4
head(QAspks)
```

```
##              Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM
## 1 101_GWI_202310_WC_LysA_10cm 1797.0279 14570.5587 56.052024 411.01717
## 2 11_MSM_202310_TR_LysB_10cm 744.7543 7094.4554 23.230016 200.12568
## 3 111_GWI_202310_SW_C 1946.6882 13991.2256 60.720156 394.67491
## 4 121_MSM_202310_TR_RHZ_SF_Tree_4 894.6568 8047.9702 27.905702 227.02314
## 5 21_MSM_202310_WC_LysB_20cm 747.1359 9510.8725 23.304301 268.28977
```

```
## 6      31_GCW_202310_TR_LysB_20cm  50.6452  776.0658  1.579701  21.89184
##      salinity S04_mM_spk S04_spk_Conc
## 1 26.329026  66.233724 7.797879e-05
## 2 12.819707  25.708406 7.797879e-05
## 3 25.282171  72.439813 7.797879e-05
## 4 14.542708  31.675677 7.797879e-05
## 5 17.186173  29.054239 7.797879e-05
## 6  1.402377   4.434775 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
```

```
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 50, QAspks$Dilution)
```

```
#Set Sample volumes in uL
```

```
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 1462, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 1501, QAspks$SampleVol)
```

```
#change sample volume to L
```

```
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)
```

```
##      Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 101_GWI_202310_WC_LysA_10cm 1797.0279 14570.5587 56.052024 411.01717
## 2  11_MSM_202310_TR_LysB_10cm  744.7543  7094.4554 23.230016 200.12568
## 3      111_GWI_202310_SW_C 1946.6882 13991.2256 60.720156 394.67491
## 4 121_MSM_202310_TR_RHZ_SF_Tree_4  894.6568  8047.9702 27.905702 227.02314
## 5   21_MSM_202310_WC_LysB_20cm  747.1359  9510.8725 23.304301 268.28977
## 6   31_GCW_202310_TR_LysB_20cm  50.6452  776.0658  1.579701  21.89184
##      salinity S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 26.329026  66.233724 7.797879e-05      1      1e-06
```

```
## 2 12.819707 25.708406 7.797879e-05 1 1e-06
## 3 25.282171 72.439813 7.797879e-05 1 1e-06
## 4 14.542708 31.675677 7.797879e-05 1 1e-06
## 5 17.186173 29.054239 7.797879e-05 1 1e-06
## 6 1.402377 4.434775 7.797879e-05 1 1e-06
```

*#gives us the total SO4 in the sample in mmoles*

```
QAspks$SO4_Total_unspkd <- (QAspks$SO4_mM/QAspks$Dilution)*(QAspks$SampleVol)
```

*##total SO4 in spiked sample in mmoles*

```
QAspks$SO4_Total_spkd <- (QAspks$SO4_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)
```

```
QAspks$SO4_expctd_spkd <- (QAspks$SO4_Total_unspkd + QAspks$SO4_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$SO4_Total_spkd/QAspks$SO4_expctd_spkd)*100
```

```
QAspks$SO4_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun')
```


```
head(QAspks)
```

```
##           Sample_ID  SO4_ppm  Cl_ppm  SO4_mM  Cl_mM
## 1 101_GWI_202310_WC_LysA_10cm 1797.0279 14570.5587 56.052024 411.01717
## 2 11_MSM_202310_TR_LysB_10cm 744.7543 7094.4554 23.230016 200.12568
## 3 111_GWI_202310_SW_C 1946.6882 13991.2256 60.720156 394.67491
## 4 121_MSM_202310_TR_RHZ_SF_Tree_4 894.6568 8047.9702 27.905702 227.02314
## 5 21_MSM_202310_WC_LysB_20cm 747.1359 9510.8725 23.304301 268.28977
## 6 31_GCW_202310_TR_LysB_20cm 50.6452 776.0658 1.579701 21.89184
## salinity SO4_mM_spk SO4_spk_Conc Dilution SampleVol SO4_Total_unspkd
## 1 26.329026 66.233724 7.797879e-05 1 1e-06 5.605202e-05
## 2 12.819707 25.708406 7.797879e-05 1 1e-06 2.323002e-05
## 3 25.282171 72.439813 7.797879e-05 1 1e-06 6.072016e-05
## 4 14.542708 31.675677 7.797879e-05 1 1e-06 2.790570e-05
## 5 17.186173 29.054239 7.797879e-05 1 1e-06 2.330430e-05
## 6 1.402377 4.434775 7.797879e-05 1 1e-06 1.579701e-06
## SO4_Total_spkd SO4_expctd_spkd spk_recovery SO4_spks_flag
## 1 7.285710e-04 1.340308e-04 543.58468 NO, rerun
## 2 2.827925e-04 1.012088e-04 279.41489 NO, rerun
## 3 7.968379e-04 1.386989e-04 574.50901 NO, rerun
## 4 3.484324e-04 1.058845e-04 329.06844 NO, rerun
## 5 3.195966e-04 1.012831e-04 315.54786 NO, rerun
## 6 4.878253e-05 7.955849e-05 61.31656 NO, rerun
```

*#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on*

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

```
spksbar
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(SO4_spks_flag) %>%
  summarise(no_rows = length(SO4_spks_flag))
Perc_spks$Total <- length(QAspks$SO4_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 NO, rerun          13    13     100
```

## Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampled2$Sample_ID), '_ ', fixed=TRUE)))

## Warning in rbind(c("1", "MSM", "202310", "UP", "LysA", "20cm"), c("10", :
## number of columns of result is not a multiple of vector length (arg 1)
```

```
colnames(IDs) <- c("Analysis_No", "Site", "Date", "Zone", "Replicate", "Depth", "RHZ", "RHZ_Rep")
head(IDs)
```

```
##   Analysis_No Site   Date Zone Replicate Depth RHZ RHZ_Rep
## 1           1  MSM 202310  UP      LysA  20cm   1    MSM
## 2          10  MSM 202310  TR      LysA  45cm  10    MSM
## 3         100  GWI 202310  TR      LysC  45cm 100    GWI
## 4         101  GWI 202310  WC      LysA  10cm 101    GWI
## 5         102  GWI 202310  WC      LysA  20cm 102    GWI
## 6         103  GWI 202310  WC      LysA  45cm 103    GWI
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampled2)
head(alldat)
```

```
##   Analysis_No Site   Date Zone Replicate Depth RHZ RHZ_Rep
## 1           1  MSM 202310  UP      LysA  20cm   1    MSM
## 2          10  MSM 202310  TR      LysA  45cm  10    MSM
## 3         100  GWI 202310  TR      LysC  45cm 100    GWI
## 4         101  GWI 202310  WC      LysA  10cm 101    GWI
## 5         102  GWI 202310  WC      LysA  20cm 102    GWI
## 6         103  GWI 202310  WC      LysA  45cm 103    GWI
##           Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM   salinity
## 1  1_MSM_202310_UP_LysA_20cm 156.4865 2514.586 4.881051 70.93331 4.543882
## 2 10_MSM_202310_TR_LysA_45cm 543.5656 4326.371 16.954635 122.04149 7.817778
## 3 100_GWI_202310_TR_LysC_45cm 1169.1941 9434.979 36.468936 266.14891 17.049033
## 4 101_GWI_202310_WC_LysA_10cm 1797.0279 14570.559 56.052024 411.01717 26.329026
## 5 102_GWI_202310_WC_LysA_20cm 1448.6028 13511.788 45.184117 381.15059 24.415828
## 6 103_GWI_202310_WC_LysA_45cm 1413.1583 12517.120 44.078550 353.09223 22.618461
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_Synoptic_CB_Dionex_ClSo4_Final_202310.csv")
```

#C

END