

# Dionex\_COMPASS\_June2023

Stephanie J. Wilson

2023-06-23

## Daily Set up

Read in and Format the raw data - change wd & file names

```
# SULFATE DATA:

## Read in raw data file from Dionex - copied and saved as a txt
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202311_S04.txt",sep='\t' , header=T, skip=3)
head(Sdat)

##      X          X.1          X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Lab Blank Unknown   n.a.    n.a.    n.a.    n.a.    n.a.
## 6 6 Standard 1 Calibration Standard 5.797  0.4644    6.54    0.5648    3.18

## Only keep the columns that we need
Sdat <- Sdat[,c(2,5)] # dont need this here
head(Sdat)

##          X.1 IC.S04.1
## 1 1 Lab Blank   n.a.
## 2 2 Lab Blank   n.a.
## 3 3 Lab Blank   n.a.
## 4 4 Lab Blank   n.a.
## 5 5 Lab Blank   n.a.
## 6 6 Standard 1  0.4644

## Name the columns correctly
colnames(Sdat) <- c( "Sample_ID", "S04_ppm")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)

## Warning: NAs introduced by coercion
```

```

Sdat <- as.data.frame(Sdat)
head(Sdat)

##      Sample_ID SO4_ppm
## 1    Lab Blank     NA
## 2    Lab Blank     NA
## 3    Lab Blank     NA
## 4    Lab Blank     NA
## 5    Lab Blank     NA
## 6 Standard 1 0.4644

#Chloride data
## Read in raw data file from Dionex - copied and saved as a txt
Cldat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202311_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)

##      X          X.1          X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank Unknown   n.a.   n.a.   n.a.   n.a.   n.a.
## 2 2 Lab Blank Unknown   n.a.   n.a.   n.a.   n.a.   n.a.
## 3 3 Lab Blank Unknown   n.a.   n.a.   n.a.   n.a.   n.a.
## 4 4 Lab Blank Unknown   n.a.   n.a.   n.a.   n.a.   n.a.
## 5 5 Lab Blank Unknown   n.a.   n.a.   n.a.   n.a.   n.a.
## 6 6 Standard 1 Calibration Standard 4.567  5.0079  93.30  8.0534  55.78

## Only keep the columns that we need
Cldat <- Cldat[,c(2,5)]
head(Cldat)

##          X.1 IC.Cl.1
## 1    Lab Blank   n.a.
## 2    Lab Blank   n.a.
## 3    Lab Blank   n.a.
## 4    Lab Blank   n.a.
## 5    Lab Blank   n.a.
## 6 Standard 1 5.0079

## Name the columns correctly
colnames(Cldat) <- c("Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)

## Warning: NAs introduced by coercion

Cldat <- as.data.frame(Cldat)
head(Cldat)

##      Sample_ID Cl_ppm
## 1    Lab Blank     NA
## 2    Lab Blank     NA
## 3    Lab Blank     NA
## 4    Lab Blank     NA
## 5    Lab Blank     NA
## 6 Standard 1 5.0079

```

```

## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)

##   Sample_ID   S04_ppm     Cl_ppm
## 1          1821.9714 14254.5661
## 2          1821.9714  4583.2258
## 3          1821.9714      0.0054
## 4          1821.9714  4220.2294
## 5          1821.9714       NA
## 6          503.9588 14254.5661

## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID=="") , ]
head(all_dat)

##                               Sample_ID   S04_ppm     Cl_ppm
## 26      1_MSM_202311_UP_LysA_20cm 200.7171 2967.900
## 27      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421
## 28    10_MSM_202311_TR_LysA_45cm_dup 534.2280 4378.928
## 29      100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329
## 30    100_GWI_202311_WC_LysA_10cm_dup 1544.0406 12622.756
## 31      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407

all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])

## If there is an n.a. in the dataframe that is because there was no peak, which
# would mean there was no Cl or SO4 and we want to know that so make all n.a.'s into zeros
all_dat[is.na(all_dat)] <- 0
head(all_dat)

##                               Sample_ID   S04_ppm     Cl_ppm
## 26      1_MSM_202311_UP_LysA_20cm 200.7171 2967.900
## 27      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421
## 28    10_MSM_202311_TR_LysA_45cm_dup 534.2280 4378.928
## 29      100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329
## 30    100_GWI_202311_WC_LysA_10cm_dup 1544.0406 12622.756
## 31      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407

```

Pull out standards - could do some checks here if we want

```

stds <- all_dat[grep("Standard", all_dat$Sample_ID),]
#stds <- stds[-c(80),]  #this is if you need to remove one for any reason
head(stds)

```

```

##   Sample_ID S04_ppm Cl_ppm
## 625 Standard 1  0.5148 5.1969
## 626 Standard 1  0.5148 5.2412
## 627 Standard 1  0.5148 5.2076

```

```

## 628 Standard 1 0.5148 5.1488
## 629 Standard 1 0.5148 5.0079
## 630 Standard 1 0.5053 5.1969

stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(SO4_ppm), sd = sd(SO4_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv <5, 'YES', 'NO, rerun')
head(stds_chk_S)

## # A tibble: 5 x 5
##   Sample_ID     mean      sd      cv flag
##   <fct>       <dbl>    <dbl>    <dbl> <chr>
## 1 Standard 1  0.496  0.0184  3.72  YES
## 2 Standard 2  1.02   0.0212  2.08  YES
## 3 Standard 3  2.04   0.0263  1.29  YES
## 4 Standard 4 10.2   0.0990  0.969 YES
## 5 Standard 5 20.3   0.225   1.11  YES

stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv <5, 'YES', 'NO, rerun')
head(stds_chk_Cl)

## # A tibble: 5 x 5
##   Sample_ID     mean      sd      cv flag
##   <fct>       <dbl>    <dbl>    <dbl> <chr>
## 1 Standard 1  5.16   0.0835  1.62  YES
## 2 Standard 2 10.4   0.134   1.28  YES
## 3 Standard 3 20.6   0.248   1.20  YES
## 4 Standard 4 102.    1.17    1.14  YES
## 5 Standard 5 203.    2.13    1.05  YES

```

## Calculate mmol/L concentrations

```

#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)

```

```

##                               Sample_ID    SO4_ppm    Cl_ppm
## 26      1_MSM_202311_UP_LysA_20cm 200.7171 2967.900
## 27      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421
## 28  10_MSM_202311_TR_LysA_45cm_dup 534.2280 4378.928
## 29      100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329
## 30  100_GWI_202311_WC_LysA_10cm_dup 1544.0406 12622.756
## 31      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407

```

```

# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06        #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$SO4_mM <- (sampledat$SO4_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)

```

```

##                               Sample_ID    SO4_ppm     Cl_ppm    SO4_mM    Cl_mM
## 26      1_MSM_202311_UP_LysA_20cm 200.7171 2967.900  6.260671 83.72074
## 27      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.619885 122.74811
## 28     10_MSM_202311_TR_LysA_45cm_dup 534.2280 4378.928 16.663381 123.52407
## 29      100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.866254 361.64539
## 30     100_GWI_202311_WC_LysA_10cm_dup 1544.0406 12622.756 48.160967 356.07211
## 31      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623
##      salinity
## 26      5.363022
## 27      7.863043
## 28      7.912749
## 29      23.166360
## 30     22.809347
## 31     21.977495

```

## Pull out dups and check with percent difference

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##                               Sample_ID    SO4_ppm     Cl_ppm    SO4_mM    Cl_mM
## 26      1_MSM_202311_UP_LysA_20cm 200.7171 2967.900  6.260671 83.72074
## 27      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.619885 122.74811
## 28     10_MSM_202311_TR_LysA_45cm_dup 534.2280 4378.928 16.663381 123.52407
## 29      100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.866254 361.64539
## 30     100_GWI_202311_WC_LysA_10cm_dup 1544.0406 12622.756 48.160967 356.07211
## 31      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623
##      salinity
## 26      5.363022
## 27      7.863043
## 28      7.912749

```

```

## 29 23.166360
## 30 22.809347
## 31 21.977495

#pull out any rows that have "dup" in the SampleID column
dups <- sampledat %>%
  filter(str_detect(Sample_ID, "dup"))      #have to change this to match data
head(dups)

##                               Sample_ID    SO4_ppm     Cl_ppm    SO4_mM     Cl_mM
## 1 10_MSM_202311_TR_LysA_45cm_dup 534.2280 4378.928 16.66338 123.5241
## 2 100_GWI_202311_WC_LysA_10cm_dup 1544.0406 12622.756 48.16097 356.0721
## 3 110_GWI_202311_SW_C_dup        1736.4688 12657.213 54.16309 357.0441
## 4 120_MSM_202311_TR_RHZ_SF_Tree_3_dup 877.2742 8017.520 27.36351 226.1642
## 5 130_MSM_202311_WC_RHZ_SF_Collar_5_dup 1036.8625 9369.887 32.34131 264.3127
## 6 20_MSM_202311_WC_LysB_20cm_dup   821.2223 9048.720 25.61517 255.2530
##   salinity
## 1 7.912749
## 2 22.809347
## 3 22.871611
## 4 14.487685
## 5 16.931411
## 6 16.351064

#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledat %>%
  filter(!str_detect(Sample_ID, "dup")) %>%
  filter(!str_detect(Sample_ID, "spk"))
head(sampledat2)

##                               Sample_ID    SO4_ppm     Cl_ppm    SO4_mM     Cl_mM  salinity
## 1 1_MSM_202311_UP_LysA_20cm 200.7171 2967.900  6.260671  83.72074 5.363022
## 2 10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.619885 122.74811 7.863043
## 3 100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.866254 361.64539 23.166360
## 4 101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623 21.977495
## 5 102_GWI_202311_WC_LysA_45cm 1199.1944 10824.812 37.404691 305.35435 19.560461
## 6 103_GWI_202311_WC_LysB_10cm 1592.8158 13596.121 49.682339 383.52951 24.568217

#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID <- gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[ , -c(2,3)]
colnames(dups) <- c('Sample_ID', 'SO4_mM_dup', "Cl_mM_dup", "salinity_dup")
head(dups)

##                               Sample_ID SO4_mM_dup Cl_mM_dup salinity_dup
## 1 10_MSM_202311_TR_LysA_45cm 16.66338 123.5241    7.912749
## 2 100_GWI_202311_WC_LysA_10cm 48.16097 356.0721   22.809347
## 3 110_GWI_202311_SW_C         54.16309 357.0441   22.871611
## 4 120_MSM_202311_TR_RHZ_SF_Tree_3 27.36351 226.1642   14.487685
## 5 130_MSM_202311_WC_RHZ_SF_Collar_5 32.34131 264.3127   16.931411
## 6 20_MSM_202311_WC_LysB_20cm   25.61517 255.2530   16.351064

```

```
#put it back together with the old data set and look for duplicates
```

```
QAdups <- merge(sampledat2, dups)
```

```
head(QAdups)
```

```
##                               Sample_ID    SO4_ppm    Cl_ppm    SO4_mM    Cl_mM
## 1      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.61988 122.7481
## 2      100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.86625 361.6454
## 3      110_GWI_202311_SW_C        1766.9177 12862.314 55.11284 362.8297
## 4      120_MSM_202311_TR_RHZ_SF_Tree_3 858.8025 7845.456 26.78735 221.3105
## 5      130_MSM_202311_WC_RHZ_SF_Collar_5 1046.7724 9517.950 32.65042 268.4894
## 6      20_MSM_202311_WC_LysB_20cm    842.8573 9269.786 26.29000 261.4890
##   salinity SO4_mM_dup Cl_mM_dup salinity_dup
## 1 7.863043 16.66338 123.5241    7.912749
## 2 23.166360 48.16097 356.0721   22.809347
## 3 23.242227 54.16309 357.0441   22.871611
## 4 14.176765 27.36351 226.1642   14.487685
## 5 17.198962 32.34131 264.3127   16.931411
## 6 16.750529 25.61517 255.2530   16.351064
```

```
QAdups$SO4_dups_chk <- ((abs(QAdups$SO4_mM-QAdups$SO4_mM_dup))/((QAdups$SO4_mM+QAdups$SO4_mM_dup)/2))*100
QAdups$SO4_dups_flag <- ifelse(QAdups$SO4_dups_chk <10, 'YES', 'NO, rerun')
```

```
QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')
```

```
head(QAdups)
```

```
##                               Sample_ID    SO4_ppm    Cl_ppm    SO4_mM    Cl_mM
## 1      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.61988 122.7481
## 2      100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.86625 361.6454
## 3      110_GWI_202311_SW_C        1766.9177 12862.314 55.11284 362.8297
## 4      120_MSM_202311_TR_RHZ_SF_Tree_3 858.8025 7845.456 26.78735 221.3105
## 5      130_MSM_202311_WC_RHZ_SF_Collar_5 1046.7724 9517.950 32.65042 268.4894
## 6      20_MSM_202311_WC_LysB_20cm    842.8573 9269.786 26.29000 261.4890
##   salinity SO4_mM_dup Cl_mM_dup salinity_dup SO4_dups_chk SO4_dups_flag
## 1 7.863043 16.66338 123.5241    7.912749 0.2613720      YES
## 2 23.166360 48.16097 356.0721   22.809347 1.4537919      YES
## 3 23.242227 54.16309 357.0441   22.871611 1.7382553      YES
## 4 14.176765 27.36351 226.1642   14.487685 2.1279820      YES
## 5 17.198962 32.34131 264.3127   16.931411 0.9512127      YES
## 6 16.750529 25.61517 255.2530   16.351064 2.6002362      YES
##   Cl_dups_chk Cl_dups_flag
## 1 0.6301604      YES
## 2 1.5530546      YES
## 3 1.6073981      YES
## 4 2.1693861      YES
## 5 1.5678185      YES
## 6 2.4135756      YES
```

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = SO4_dups_chk, fill=SO4_dups_chk)) +
```

```
geom_bar(stat = 'identity') +
```

```


scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
theme_classic() + labs(x= "Sample ID", y="Difference Between SO4 Duplicates (%)") +
theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
color = "black", size=1)

## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call `lifecycle::last_lifecycle_warnings()` to see where this warning was
## generated.

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
geom_bar(stat = 'identity') +
scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))

```

```

Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES      13  YES      13     13      100       100

```

Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2, 1, sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)

##                               Sample_ID    S04_ppm    Cl_ppm    S04_mM    Cl_mM
## 1 10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.61988 122.7481
## 2 100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.86625 361.6454
## 3 110_GWI_202311_SW_C          1766.9177 12862.314 55.11284 362.8297
## 4 120_MSM_202311_TR_RHZ_SF_Tree_3 858.8025 7845.456 26.78735 221.3105
## 5 130_MSM_202311_WC_RHZ_SF_Collar_5 1046.7724 9517.950 32.65042 268.4894
## 6 20_MSM_202311_WC_LysB_20cm    842.8573 9269.786 26.29000 261.4890
##   salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 7.863043 16.66338 123.5241    7.912749  0.2613720      YES
## 2 23.166360 48.16097 356.0721   22.809347  1.4537919      YES
## 3 23.242227 54.16309 357.0441   22.871611  1.7382553      YES
## 4 14.176765 27.36351 226.1642   14.487685  2.1279820      YES
## 5 17.198962 32.34131 264.3127   16.931411  0.9512127      YES
## 6 16.750529 25.61517 255.2530   16.351064  2.6002362      YES
##   Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1 0.6301604      YES  0.1848179      YES
## 2 1.5530546      YES  1.0279861      YES
## 3 1.6073981      YES  1.2291321      YES
## 4 2.1693861      YES  1.5047105      YES
## 5 1.5678185      YES  0.6726090      YES
## 6 2.4135756      YES  1.8386446      YES

```


```

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between SO4 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))

```

```

colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1 YES      13   YES      13     13      100      100

```

## Pull out spikes and check

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##                               Sample_ID    SO4_ppm    Cl_ppm    SO4_mM    Cl_mM
## 26      1_MSM_202311_UP_LysA_20cm 200.7171 2967.900  6.260671 83.72074
## 27      10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.619885 122.74811
## 28     10_MSM_202311_TR_LysA_45cm_dup 534.2280 4378.928 16.663381 123.52407
## 29     100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.866254 361.64539
## 30    100_GWI_202311_WC_LysA_10cm_dup 1544.0406 12622.756 48.160967 356.07211
## 31     101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623
##   salinity
## 26  5.363022
## 27  7.863043
## 28  7.912749
## 29 23.166360
## 30 22.809347
## 31 21.977495

```

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk"))      #have to change this to match data
head(spks)

```

```

##                               Sample_ID    SO4_ppm    Cl_ppm    SO4_mM    Cl_mM
## 1     101_GWI_202311_WC_LysA_20cm_spk 1651.5201 11892.932 51.51342 335.48469
## 2     11_MSM_202311_TR_LysB_10cm_spk  803.5885  7175.359 25.06514 202.40786
## 3 111_MSM_202311_UP_RHZ_SF_Tree_1_spk 398.1390  2820.719 12.41856 79.56895
## 4 121_MSM_202311_TR_RHZ_SF_Tree_4_spk 999.5127  7901.435 31.17632 222.88955
## 5     131_MSM_202311_WC_RHZ_LysA_spk 1481.2258  9900.823 46.20168 279.28978
## 6     21_MSM_202311_WC_LysB_45cm_spk  328.0532  8615.336 10.23248 243.02780
##   salinity
## 1 21.490554
## 2 12.965899
## 3 5.097066
## 4 14.277918
## 5 17.890813
## 6 15.567937

```

```

#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[ ,-c(2,3, 5,6)]
colnames(spks) <- c('Sample_ID', 'S04_mM_spk')
head(spks)

##                                     Sample_ID S04_mM_spk
## 1      101_GWI_202311_WC_LysA_20cm   51.51342
## 2      11_MSM_202311_TR_LysB_10cm   25.06514
## 3 111_MSM_202311_UP_RHZ_SF_Tree_1  12.41856
## 4 121_MSM_202311_TR_RHZ_SF_Tree_4  31.17632
## 5      131_MSM_202311_WC_RHZ_LysA  46.20168
## 6      21_MSM_202311_WC_LysB_45cm  10.23248

#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat, spks)
head(QAspks)

##                                     Sample_ID S04_ppm    Cl_ppm    S04_mM    Cl_mM
## 1      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623
## 2      11_MSM_202311_TR_LysB_10cm   709.2730  7178.410 22.123300 202.49394
## 3 111_MSM_202311_UP_RHZ_SF_Tree_1  310.6557  2855.731  9.689822 80.55658
## 4 121_MSM_202311_TR_RHZ_SF_Tree_4  961.2419  8335.104 29.982592 235.12283
## 5      131_MSM_202311_WC_RHZ_LysA 1320.9987 10110.144 41.203952 285.19447
## 6      21_MSM_202311_WC_LysB_45cm   145.5572  8790.947  4.540150 247.98159
##     salinity S04_mM_spk
## 1 21.977495   51.51342
## 2 12.971413   25.06514
## 3 5.160331    12.41856
## 4 15.061560   31.17632
## 5 18.269056   46.20168
## 6 15.885268   10.23248

#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvolt <- 10              # in uL
spkvolt <- spkvolt/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvolt          # mmoles of S04
head(QAspks)

##                                     Sample_ID S04_ppm    Cl_ppm    S04_mM    Cl_mM
## 1      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623
## 2      11_MSM_202311_TR_LysB_10cm   709.2730  7178.410 22.123300 202.49394
## 3 111_MSM_202311_UP_RHZ_SF_Tree_1  310.6557  2855.731  9.689822 80.55658
## 4 121_MSM_202311_TR_RHZ_SF_Tree_4  961.2419  8335.104 29.982592 235.12283
## 5      131_MSM_202311_WC_RHZ_LysA 1320.9987 10110.144 41.203952 285.19447
## 6      21_MSM_202311_WC_LysB_45cm   145.5572  8790.947  4.540150 247.98159
##     salinity S04_mM_spk S04_spk_Conc
## 1 21.977495   51.51342 7.797879e-05
## 2 12.971413   25.06514 7.797879e-05
## 3 5.160331    12.41856 7.797879e-05

```

```

## 4 15.061560 31.17632 7.797879e-05
## 5 18.269056 46.20168 7.797879e-05
## 6 15.885268 10.23248 7.797879e-05

#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...

QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 50, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_202306_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_202306_WC"), 1462, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_202306_WC"), 1501, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

##                               Sample_ID    SO4_ppm    Cl_ppm    SO4_mM    Cl_mM
## 1 101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623
## 2 11_MSM_202311_TR_LysB_10cm  709.2730  7178.410 22.123300 202.49394
## 3 111_MSM_202311_UP_RHZ_SF_Tree_1 310.6557 2855.731 9.689822 80.55658
## 4 121_MSM_202311_TR_RHZ_SF_Tree_4 961.2419 8335.104 29.982592 235.12283
## 5 131_MSM_202311_WC_RHZ_LysA 1320.9987 10110.144 41.203952 285.19447
## 6 21_MSM_202311_WC_LysB_45cm  145.5572  8790.947  4.540150 247.98159
##   salinity SO4_mM_spk SO4_spk_Conc Dilution SampleVol
## 1 21.977495 51.51342 7.797879e-05      1 1e-06
## 2 12.971413 25.06514 7.797879e-05      1 1e-06
## 3 5.160331 12.41856 7.797879e-05      1 1e-06
## 4 15.061560 31.17632 7.797879e-05      1 1e-06
## 5 18.269056 46.20168 7.797879e-05      1 1e-06
## 6 15.885268 10.23248 7.797879e-05      1 1e-06

```

```

#gives us the total SO4 in the sample in mmoles
QAspks$SO4_Total_unspkd <- (QAspks$SO4_mM/QAspks$Dilution)*(QAspks$SampleVol)

##total SO4 in spiked sample in mmoles
QAspks$SO4_Total_spkd <- (QAspks$SO4_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)

QAspks$SO4_expctd_spkd <- (QAspks$SO4_Total_unspkd + QAspks$SO4_spk_Conc)
QAspks$spk_recovery <- (QAspks$SO4_Total_spkd/QAspks$SO4_expctd_spkd)*100
QAspks$SO4_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES' , 'NO, rerun'


head(QAspks)

##                                         Sample_ID    SO4_ppm      Cl_ppm      SO4_mM      Cl_mM
## 1      101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623
## 2      11_MSM_202311_TR_LysB_10cm   709.2730  7178.410 22.123300 202.49394
## 3  111_MSM_202311_UP_RHZ_SF_Tree_1 310.6557  2855.731  9.689822 80.55658
## 4 121_MSM_202311_TR_RHZ_SF_Tree_4  961.2419  8335.104 29.982592 235.12283
## 5      131_MSM_202311_WC_RHZ_LysA 1320.9987 10110.144 41.203952 285.19447
## 6      21_MSM_202311_WC_LysB_45cm  145.5572  8790.947  4.540150 247.98159
##   salinity SO4_mM_spk SO4_spk_Conc Dilution SampleVol SO4_Total_unspkd
## 1 21.977495  51.51342 7.797879e-05       1  1e-06  4.109203e-05
## 2 12.971413  25.06514 7.797879e-05       1  1e-06  2.212330e-05
## 3  5.160331  12.41856 7.797879e-05       1  1e-06  9.689822e-06
## 4 15.061560  31.17632 7.797879e-05       1  1e-06  2.998259e-05
## 5 18.269056  46.20168 7.797879e-05       1  1e-06  4.120395e-05
## 6 15.885268  10.23248 7.797879e-05       1  1e-06  4.540150e-06
##   SO4_Total_spkd SO4_expctd_spkd spk_recovery SO4_spks_flag
## 1  0.0005666476  1.190708e-04    475.8912  NO, rerun
## 2  0.0002757166  1.001021e-04    275.4354  NO, rerun
## 3  0.0001366041  8.766861e-05    155.8188  NO, rerun
## 4  0.0003429395  1.079614e-04    317.6501  NO, rerun
## 5  0.0005082185  1.191827e-04    426.4195  NO, rerun
## 6  0.0001125572  8.251894e-05    136.4017  NO, rerun

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)

spksbar

```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(SO4_spks_flag) %>%
  summarise(no_rows = length(SO4_spks_flag))
Perc_spks$Total <- length(QAspks$SO4_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 2 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>          <int> <int>   <dbl>
## 1 NO, rerun      11     13    84.6
## 2 YES            2      13    15.4
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_', fixed=TRUE)))

## Warning in rbind(c("1", "MSM", "202311", "UP", "LysA", "20cm"), c("10", :
## number of columns of result is not a multiple of vector length (arg 1)
```

```

colnames(IDs) <- c("Analysis_No", "Site", "Date", "Zone", "Replicate", "Depth", "RHZ", "RHZ_Rep")
head(IDs)

##   Analysis_No Site Date Zone Replicate Depth RHZ RHZ_Rep
## 1             1 MSM 202311   UP     LysA 20cm    1    MSM
## 2            10 MSM 202311   TR     LysA 45cm   10    MSM
## 3           100 GWI 202311   WC     LysA 10cm  100   GWI
## 4           101 GWI 202311   WC     LysA 20cm 101   GWI
## 5           102 GWI 202311   WC     LysA 45cm 102   GWI
## 6           103 GWI 202311   WC     LysB 10cm 103   GWI

#rejoin them to the dataframe
alldat <- cbind(IDs, sampledat2)
head(alldat)

##   Analysis_No Site Date Zone Replicate Depth RHZ RHZ_Rep
## 1             1 MSM 202311   UP     LysA 20cm    1    MSM
## 2            10 MSM 202311   TR     LysA 45cm   10    MSM
## 3           100 GWI 202311   WC     LysA 10cm  100   GWI
## 4           101 GWI 202311   WC     LysA 20cm 101   GWI
## 5           102 GWI 202311   WC     LysA 45cm 102   GWI
## 6           103 GWI 202311   WC     LysB 10cm 103   GWI
##               Sample_ID   SO4_ppm     Cl_ppm   SO4_mM     Cl_mM salinity
## 1 1_MSM_202311_UP_LysA_20cm 200.7171 2967.900 6.260671 83.72074 5.363022
## 2 10_MSM_202311_TR_LysA_45cm 532.8335 4351.421 16.619885 122.74811 7.863043
## 3 100_GWI_202311_WC_LysA_10cm 1566.6521 12820.329 48.866254 361.64539 23.166360
## 4 101_GWI_202311_WC_LysA_20cm 1317.4105 12162.407 41.092031 343.08623 21.977495
## 5 102_GWI_202311_WC_LysA_45cm 1199.1944 10824.812 37.404691 305.35435 19.560461
## 6 103_GWI_202311_WC_LysB_10cm 1592.8158 13596.121 49.682339 383.52951 24.568217

```

## Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_SO4_202311.csv")
```

END