

Dionex_COMPASS_November2022

Stephanie J. Wilson

2023-01-23

Daily Set up

Read in and Format the raw data - change wd & file names (Make sure files are UTF-8 encoded)

```
# SULFATE DATA:
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Sdat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202211_S04.txt", sep='\t' , header=T, skip=3)  
head(Sdat)
```

```
##      X      X.1      X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4  
## 1 1 Lab Blank      Unknown  n.a.      n.a.      n.a.      n.a.      n.a.  
## 2 2 Lab Blank      Unknown  n.a.      n.a.      n.a.      n.a.      n.a.  
## 3 3 Standard 1      Unknown  4.397    0.5017    5.85    0.4070    2.53  
## 4 4 Blank          Unknown  n.a.      n.a.      n.a.      n.a.      n.a.  
## 5 5 Standard 1 Calibration Standard 4.393    0.5034    5.87    0.4084    2.54  
## 6 6 Standard 2 Calibration Standard 4.400    1.0993    6.12    0.8917    5.39
```

```
## Only keep the columns that we need
```

```
Sdat <- Sdat[,c(2,5)] # dont need this here  
head(Sdat)
```

```
##      X.1 IC.S04.1  
## 1 Lab Blank      n.a.  
## 2 Lab Blank      n.a.  
## 3 Standard 1 0.5017  
## 4 Blank      n.a.  
## 5 Standard 1 0.5034  
## 6 Standard 2 1.0993
```

```
## Name the columns correctly
```

```
colnames(Sdat) <- c( "Sample_ID", "S04_ppm")  
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)  
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID S04_ppm
## 1 Lab Blank      NA
## 2 Lab Blank      NA
## 3 Standard 1  0.5017
## 4      Blank      NA
## 5 Standard 1  0.5034
## 6 Standard 2  1.0993
```

```
#Chloride data
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Cldat <- read.table("Raw Data/COMPASS_Synoptic_CB_MonMon_202211_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X      X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank      Unknown 3.270  0.0241  100.00  0.0271  0.19
## 2 2 Lab Blank      Unknown 3.270  0.0251  100.00  0.0282  0.20
## 3 3 Standard 1      Unknown 3.273  5.7878  93.58  6.5054  43.94
## 4 4      Blank      Unknown 3.270  0.0260  100.00  0.0292  0.21
## 5 5 Standard 1 Calibration Standard 3.273  5.7959  93.56  6.5145  44.12
## 6 6 Standard 2 Calibration Standard 3.273  12.0880  93.22  13.5866  90.14
```

```
## Only keep the columns that we need
```

```
Cldat <- Cldat[,c(2,5)]
head(Cldat)
```

```
##      X.1 IC.Cl.1
## 1 Lab Blank 0.0241
## 2 Lab Blank 0.0251
## 3 Standard 1 5.7878
## 4      Blank 0.0260
## 5 Standard 1 5.7959
## 6 Standard 2 12.0880
```

```
## Name the columns correctly
```

```
colnames(Cldat) <- c("Sample_ID", "Cl_ppm")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##      Sample_ID Cl_ppm
## 1 Lab Blank 0.0241
## 2 Lab Blank 0.0251
## 3 Standard 1 5.7878
## 4      Blank 0.0260
## 5 Standard 1 5.7959
## 6 Standard 2 12.0880
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##      Sample_ID      S04_ppm      Cl_ppm
## 1           2544.4789 19269.7855
## 2           2544.4789  5007.1144
## 3           2544.4789    0.0179
## 4           2544.4789  5111.7920
## 5           2544.4789         NA
## 6           784.5324 19269.7855
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##              Sample_ID      S04_ppm      Cl_ppm
## 26      1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644
## 27     10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636
## 28 10_202211_MSM_TR_LysA_20cm_spk 1248.1115 6395.9330
## 29 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034
## 30 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425
## 31 100_202211_GCrew_TR_LysC_45cm 19.3955 134.0034
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)
```

```
##              Sample_ID      S04_ppm      Cl_ppm
## 26      1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644
## 27     10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636
## 28 10_202211_MSM_TR_LysA_20cm_spk 1248.1115 6395.9330
## 29 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034
## 30 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425
## 31 100_202211_GCrew_TR_LysC_45cm 19.3955 134.0034
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grep1("Standard", all_dat$Sample_ID),]
head(stds)
```

```
##      Sample_ID S04_ppm Cl_ppm
## 1730 Standard 1  0.5034 5.7959
## 1731 Standard 1  0.5034 5.8500
## 1732 Standard 1  0.5034 5.8551
## 1733 Standard 1  0.5034 5.9283
## 1734 Standard 1  0.5034 5.9000
## 1735 Standard 1  0.5034 5.7635
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(SO4_ppm), sd = sd(SO4_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID    mean      sd    cv flag
##   <fct>      <dbl>  <dbl> <dbl> <chr>
## 1 Standard 1  0.506 0.0109 2.16 NO, rerun
## 2 Standard 2  1.11  0.0162 1.46 YES
## 3 Standard 3  2.22  0.0416 1.88 YES
## 4 Standard 4  8.80  0.100  1.14 YES
## 5 Standard 5 20.4   0.190  0.935 YES
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID    mean      sd    cv flag
##   <fct>      <dbl>  <dbl> <dbl> <chr>
## 1 Standard 1   5.84 0.0570 0.975 YES
## 2 Standard 2  12.2 0.0665 0.546 YES
## 3 Standard 3  23.4 0.186  0.796 YES
## 4 Standard 4  89.3 0.628  0.704 YES
## 5 Standard 5 203.   1.53   0.753 YES
```

Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##           Sample_ID  SO4_ppm  Cl_ppm
## 26  1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644
## 27  10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636
## 28 10_202211_MSM_TR_LysA_20cm_spk 1248.1115 6395.9330
## 29 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034
## 30 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425
## 31 100_202211_GCrew_TR_LysC_45cm 19.3955 134.0034
```

```
# Constants needed for calculations:
clmw <- 35.45    #molecular weight of Chloride: 35.45
smw <- 32.06     #molecular weight of sulfur: 32.06
```

```

# Convert ppm to mmol/L
sampledat$SO4_mM <- (sampledat$SO4_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)

```

```

##           Sample_ID   SO4_ppm   Cl_ppm   SO4_mM   Cl_mM
## 26      1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644 8.254775 92.266979
## 27     10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636 41.341089 193.951018
## 28 10_202211_MSM_TR_LysA_20cm_spk 1248.1115 6395.9330 38.930490 180.421241
## 29 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034 4.951578 3.780068
## 30 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425 4.951578 3.758604
## 31 100_202211_GCrew_TR_LysC_45cm 19.3955 134.0034 0.604975 3.780068
##      salinity
## 26 5.9104780
## 27 12.4241694
## 28 11.5574769
## 29 0.2421701
## 30 0.2407952
## 31 0.2421701

```

Pull out dups and check with percent difference

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##           Sample_ID   SO4_ppm   Cl_ppm   SO4_mM   Cl_mM
## 26      1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644 8.254775 92.266979
## 27     10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636 41.341089 193.951018
## 28 10_202211_MSM_TR_LysA_20cm_spk 1248.1115 6395.9330 38.930490 180.421241
## 29 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034 4.951578 3.780068
## 30 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425 4.951578 3.758604
## 31 100_202211_GCrew_TR_LysC_45cm 19.3955 134.0034 0.604975 3.780068
##      salinity
## 26 5.9104780
## 27 12.4241694
## 28 11.5574769
## 29 0.2421701
## 30 0.2407952
## 31 0.2421701

```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledats %>%
  filter(str_detect(Sample_ID, "dup")) #have to change this to match data
head(dups)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1 109_202211_GCrew_WC_SipC_45cm_dup 490.3857 5847.447 15.29587 164.9491
## 2  129_202211_GWI_WC_SipA_10cm_dup 2192.9807 15519.800 68.40239 437.7941
## 3      139_202211_GWI_SW_C_dup 1972.7418 14364.361 61.53281 405.2006
## 4   19_202211_MSM_WC_SipA_20cm_dup 2527.3082 14813.024 78.83057 417.8568
## 5      29_202211_MSM_SW_C_dup 2245.2770 11988.478 70.03359 338.1799
## 6   39_202211_MSM_RHZ_TR_SF_2_dup 658.3043 4849.288 20.53351 136.7923
##   salinity
## 1 10.56636
## 2 28.04430
## 3 25.95643
## 4 26.76716
## 5 21.66321
## 6 8.76269
```

```
#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledats %>%
  filter(!str_detect(Sample_ID, "dup")) %>%
  filter(!str_detect(Sample_ID, "spk"))
head(sampledat2)
```

```
##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM
## 1   1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644 8.254775 92.266979
## 2  10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636 41.341089 193.951018
## 3 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034 4.951578 3.780068
## 4 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425 4.951578 3.758604
## 5 100_202211_GCrew_TR_LysC_45cm 19.3955 134.0034 0.604975 3.780068
## 6 100_202211_GCrew_TR_LysC_45cm 19.3955 133.2425 0.604975 3.758604
##   salinity
## 1 5.9104780
## 2 12.4241694
## 3 0.2421701
## 4 0.2407952
## 5 0.2421701
## 6 0.2407952
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup", "", as.character(dups$Sample_ID))
dups <- dups[, -c(2,3)]
colnames(dups) <- c('Sample_ID', 'S04_mM_dup', 'Cl_mM_dup', 'salinity_dup')
head(dups)
```

```
##           Sample_ID S04_mM_dup Cl_mM_dup salinity_dup
## 1 109_202211_GCrew_WC_SipC_45cm 15.29587 164.9491 10.56636
## 2  129_202211_GWI_WC_SipA_10cm 68.40239 437.7941 28.04430
## 3      139_202211_GWI_SW_C 61.53281 405.2006 25.95643
## 4   19_202211_MSM_WC_SipA_20cm 78.83057 417.8568 26.76716
```

```
## 5          29_202211_MSM_SW_C    70.03359  338.1799    21.66321
## 6          39_202211_MSM_RHZ_TR_SF_2    20.53351  136.7923    8.76269
```

#put it back together with the old data set and look for duplicates

```
QAdups <- merge(sampledat2, dups)
head(QAdups)
```

```
##          Sample_ID    S04_ppm    Cl_ppm    S04_mM    Cl_mM    salinity
## 1 109_202211_GCrew_WC_SipC_45cm  467.7261  5581.663  14.58909  157.4517  10.086092
## 2          139_202211_GWI_SW_C  2131.0356  15474.590  66.47023  436.5188  27.962610
## 3   19_202211_MSM_WC_SipA_20cm  2365.1949  14059.697  73.77401  396.6064  25.405898
## 4          29_202211_MSM_SW_C  2098.9347  12472.348  65.46896  351.8293  22.537559
## 5   39_202211_MSM_RHZ_TR_SF_2   713.7607   5163.679  22.26328  145.6609   9.330794
## 6   49_202211_MSM_RHZ_WC_Co14  1787.4021  11210.941  55.75178  316.2466  20.258196
##    S04_mM_dup Cl_mM_dup salinity_dup
## 1   15.29587  164.9491    10.56636
## 2   61.53281  405.2006    25.95643
## 3   78.83057  417.8568    26.76716
## 4   70.03359  338.1799    21.66321
## 5   20.53351  136.7923     8.76269
## 6   58.95680  324.4860    20.78600
```

```
QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*100
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')
```

```
QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')
```

```
head(QAdups)
```


```
##          Sample_ID    S04_ppm    Cl_ppm    S04_mM    Cl_mM    salinity
## 1 109_202211_GCrew_WC_SipC_45cm  467.7261  5581.663  14.58909  157.4517  10.086092
## 2          139_202211_GWI_SW_C  2131.0356  15474.590  66.47023  436.5188  27.962610
## 3   19_202211_MSM_WC_SipA_20cm  2365.1949  14059.697  73.77401  396.6064  25.405898
## 4          29_202211_MSM_SW_C  2098.9347  12472.348  65.46896  351.8293  22.537559
## 5   39_202211_MSM_RHZ_TR_SF_2   713.7607   5163.679  22.26328  145.6609   9.330794
## 6   49_202211_MSM_RHZ_WC_Co14  1787.4021  11210.941  55.75178  316.2466  20.258196
##    S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag Cl_dups_chk
## 1   15.29587  164.9491    10.56636      4.730053        YES    4.650997
## 2   61.53281  405.2006    25.95643      7.714541        YES    7.441473
## 3   78.83057  417.8568    26.76716      6.627009        YES    5.218266
## 4   70.03359  338.1799    21.66321      6.737347        YES    3.956287
## 5   20.53351  136.7923     8.76269      8.083640        YES    6.279669
## 6   58.95680  324.4860    20.78600      5.588112        YES    2.571869
##    Cl_dups_flag
## 1          YES
## 2          YES
## 3          YES
## 4          YES
## 5          YES
## 6          YES
```

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar,ncol=2, nrow=1)
```




```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##           Flag S_no_rows      Flag Cl_no_rows Total S_Percent Cl_Percent
## 1 NO, rerun          3 NO, rerun          1     12         25  8.333333
## 2      YES          9      YES          11     12         75 91.666667

```

Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)

```

```

##           Sample_ID  S04_ppm  Cl_ppm  S04_mM  Cl_mM  salinity
## 1 109_202211_GCrew_WC_SipC_45cm 467.7261 5581.663 14.58909 157.4517 10.086092
## 2           139_202211_GWI_SW_C 2131.0356 15474.590 66.47023 436.5188 27.962610
## 3    19_202211_MSM_WC_SipA_20cm 2365.1949 14059.697 73.77401 396.6064 25.405898
## 4           29_202211_MSM_SW_C 2098.9347 12472.348 65.46896 351.8293 22.537559
## 5    39_202211_MSM_RHZ_TR_SF_2  713.7607  5163.679 22.26328 145.6609  9.330794
## 6    49_202211_MSM_RHZ_WC_Col4 1787.4021 11210.941 55.75178 316.2466 20.258196
##   S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag Cl_dups_chk
## 1  15.29587  164.9491   10.56636    4.730053          YES    4.650997
## 2  61.53281  405.2006   25.95643    7.714541          YES    7.441473
## 3  78.83057  417.8568   26.76716    6.627009          YES    5.218266
## 4  70.03359  338.1799   21.66321    6.737347          YES    3.956287
## 5  20.53351  136.7923    8.76269    8.083640          YES    6.279669
## 6  58.95680  324.4860   20.78600    5.588112          YES    2.571869
##   Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1           YES    3.344653          YES

```


## 2	YES	5.455004	YES
## 3	YES	4.686003	YES
## 4	YES	4.764023	YES
## 5	YES	5.715997	YES
## 6	YES	3.951392	YES

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this

```
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```



```

#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(S04_dups_flag) %>%
  summarise(S_no_rows = length(S04_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$S04_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)

```

```

##           Flag S_no_rows           Flag Cl_no_rows Total S_Percent Cl_Percent
## 1 NO, rerun           3 NO, rerun           1      12          25  8.333333
## 2      YES           9      YES           11      12          75 91.666667

```

Pull out spikes and check

```

#Show me the data that we have from the calculations
head(sampledat)

```

```

##           Sample_ID      S04_ppm      Cl_ppm      S04_mM      Cl_mM
## 26  1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644  8.254775 92.266979
## 27  10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636 41.341089 193.951018
## 28 10_202211_MSM_TR_LysA_20cm_spk 1248.1115 6395.9330 38.930490 180.421241
## 29 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034  4.951578  3.780068
## 30 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425  4.951578  3.758604
## 31 100_202211_GCrew_TR_LysC_45cm  19.3955 134.0034  0.604975  3.780068
##           salinity
## 26  5.9104780
## 27 12.4241694
## 28 11.5574769
## 29  0.2421701
## 30  0.2407952
## 31  0.2421701

```

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)

```

```

##           Sample_ID      S04_ppm      Cl_ppm      S04_mM      Cl_mM
## 1  10_202211_MSM_TR_LysA_20cm_spk 1248.1115 6395.933 38.93049 180.4212
## 2    110_202211_GCrew_SW_A_spk  956.5879 5163.239 29.83743 145.6485
## 3 120_202211_GWI_TR_LysA_10cm_spk 1664.3307 12466.940 51.91300 351.6767
## 4 130_202211_GWI_WC_SipA_20cm_spk 1900.5139 13145.045 59.27991 370.8052

```

```
## 5 140_202211_GWI_RHZ_TR_SF_1_spk 1606.2517 11339.052 50.10143 319.8604
## 6 20_202211_MSM_WC_SipA_45cm_spk 326.3741 12616.349 10.18010 355.8914
## salinity
## 1 11.557477
## 2 9.329999
## 3 22.527787
## 4 23.753122
## 5 20.489693
## 6 22.797768
```

```
#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[ ,-c(2,3, 5,6)]
colnames(spks) <- c('Sample_ID', 'SO4_mM_spk')
head(spks)
```

```
## Sample_ID SO4_mM_spk
## 1 10_202211_MSM_TR_LysA_20cm 38.93049
## 2 110_202211_GCrew_SW_A 29.83743
## 3 120_202211_GWI_TR_LysA_10cm 51.91300
## 4 130_202211_GWI_WC_SipA_20cm 59.27991
## 5 140_202211_GWI_RHZ_TR_SF_1 50.10143
## 6 20_202211_MSM_WC_SipA_45cm 10.18010
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampldat, spks)
head(QAspks)
```

```
## Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM salinity
## 1 10_202211_MSM_TR_LysA_20cm 1325.3953 6875.564 41.3410886 193.9510 12.424169
## 2 110_202211_GCrew_SW_A 663.6499 5151.559 20.7002464 145.3190 9.308893
## 3 120_202211_GWI_TR_LysA_10cm 1862.6278 14578.719 58.0981847 411.2474 26.343772
## 4 130_202211_GWI_WC_SipA_20cm 1879.0113 13943.574 58.6092109 393.3307 25.196064
## 5 140_202211_GWI_RHZ_TR_SF_1 1641.0884 12090.940 51.1880349 341.0702 21.848355
## 6 20_202211_MSM_WC_SipA_45cm 29.5728 12145.786 0.9224205 342.6174 21.947462
## SO4_mM_spk
## 1 38.93049
## 2 29.83743
## 3 51.91300
## 4 59.27991
## 5 50.10143
## 6 10.18010
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw) # in mM
spkvol <- 10 # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$SO4_spk_Conc <- (spkconc)*spkvol # mmoles of SO4
head(QAspks)
```

```
## Sample_ID SO4_ppm Cl_ppm SO4_mM Cl_mM salinity
```

```
## 1 10_202211_MSM_TR_LysA_20cm 1325.3953 6875.564 41.3410886 193.9510 12.424169
## 2 110_202211_GCrew_SW_A 663.6499 5151.559 20.7002464 145.3190 9.308893
## 3 120_202211_GWI_TR_LysA_10cm 1862.6278 14578.719 58.0981847 411.2474 26.343772
## 4 130_202211_GWI_WC_SipA_20cm 1879.0113 13943.574 58.6092109 393.3307 25.196064
## 5 140_202211_GWI_RHZ_TR_SF_1 1641.0884 12090.940 51.1880349 341.0702 21.848355
## 6 20_202211_MSM_WC_SipA_45cm 29.5728 12145.786 0.9224205 342.6174 21.947462
## S04_mM_spk S04_spk_Conc
## 1 38.93049 7.797879e-05
## 2 29.83743 7.797879e-05
## 3 51.91300 7.797879e-05
## 4 59.27991 7.797879e-05
## 5 50.10143 7.797879e-05
## 6 10.18010 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
```

```
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)
```

```
#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)
```

```
## Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM salinity
## 1 10_202211_MSM_TR_LysA_20cm 1325.3953 6875.564 41.3410886 193.9510 12.424169
## 2 110_202211_GCrew_SW_A 663.6499 5151.559 20.7002464 145.3190 9.308893
## 3 120_202211_GWI_TR_LysA_10cm 1862.6278 14578.719 58.0981847 411.2474 26.343772
```

```
## 4 130_202211_GWI_WC_SipA_20cm 1879.0113 13943.574 58.6092109 393.3307 25.196064
## 5 140_202211_GWI_RHZ_TR_SF_1 1641.0884 12090.940 51.1880349 341.0702 21.848355
## 6 20_202211_MSM_WC_SipA_45cm 29.5728 12145.786 0.9224205 342.6174 21.947462
## S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 38.93049 7.797879e-05 50 0.001501
## 2 29.83743 7.797879e-05 1 0.000001
## 3 51.91300 7.797879e-05 100 0.001475
## 4 59.27991 7.797879e-05 200 0.001462
## 5 50.10143 7.797879e-05 1 0.000001
## 6 10.18010 7.797879e-05 100 0.001475
```

#gives us the total SO4 in the sample in mmoles

```
QAspks$S04_Total_unspkd <- (QAspks$S04_mM/QAspks$Dilution)*(QAspks$SampleVol)
```

##total SO4 in spiked sample in mmoles

```
QAspks$S04_Total_spkd <- (QAspks$S04_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)
```

```
QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
```

```
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun
```


```
head(QAspks)
```

```
## Sample_ID S04_ppm Cl_ppm S04_mM Cl_mM salinity
## 1 10_202211_MSM_TR_LysA_20cm 1325.3953 6875.564 41.3410886 193.9510 12.424169
## 2 110_202211_GCrew_SW_A 663.6499 5151.559 20.7002464 145.3190 9.308893
## 3 120_202211_GWI_TR_LysA_10cm 1862.6278 14578.719 58.0981847 411.2474 26.343772
## 4 130_202211_GWI_WC_SipA_20cm 1879.0113 13943.574 58.6092109 393.3307 25.196064
## 5 140_202211_GWI_RHZ_TR_SF_1 1641.0884 12090.940 51.1880349 341.0702 21.848355
## 6 20_202211_MSM_WC_SipA_45cm 29.5728 12145.786 0.9224205 342.6174 21.947462
## S04_mM_spk S04_spk_Conc Dilution SampleVol S04_Total_unspkd S04_Total_spkd
## 1 38.93049 7.797879e-05 50 0.001501 1.241059e-03 0.0011764794
## 2 29.83743 7.797879e-05 1 0.000001 2.070025e-05 0.0003282117
## 3 51.91300 7.797879e-05 100 0.001475 8.569482e-04 0.0007709080
## 4 59.27991 7.797879e-05 200 0.001462 4.284333e-04 0.0004363001
## 5 50.10143 7.797879e-05 1 0.000001 5.118803e-05 0.0005511157
## 6 10.18010 7.797879e-05 100 0.001475 1.360570e-05 0.0001511745
## S04_expctd_spkd spk_recovery S04_spks_flag
## 1 1.319038e-03 89.19221 YES
## 2 9.867904e-05 332.60529 NO, rerun
## 3 9.349270e-04 82.45649 YES
## 4 5.064121e-04 86.15515 YES
## 5 1.291668e-04 426.66968 NO, rerun
## 6 9.158449e-05 165.06564 NO, rerun
```

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

spksbar



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 2 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 NO, rerun          10    13   76.9
## 2 YES                 3    13   23.1
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_', fixed=TRUE)))
```

```
## Warning in rbind(c("1", "202211", "MSM", "UP", "LysA", "10cm"), c("10", :
## number of columns of result is not a multiple of vector length (arg 1)
```

```
colnames(IDs) <- c("Analysis_No" , "Date","Site", "Zone", "Replicate", "Depth", "Tree")
head(IDs)
```

```
##   Analysis_No   Date   Site Zone Replicate Depth Tree
## 1           1 202211   MSM  UP      LysA   10cm    1
## 2          10 202211   MSM  TR      LysA   20cm   10
## 3         100 202211 GCrew  TR      LysC   45cm  100
## 4         100 202211 GCrew  TR      LysC   45cm  100
## 5         100 202211 GCrew  TR      LysC   45cm  100
## 6         100 202211 GCrew  TR      LysC   45cm  100
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampled2)
head(alldat)
```

```
##   Analysis_No   Date   Site Zone Replicate Depth Tree
## 1           1 202211   MSM  UP      LysA   10cm    1
## 2          10 202211   MSM  TR      LysA   20cm   10
## 3         100 202211 GCrew  TR      LysC   45cm  100
## 4         100 202211 GCrew  TR      LysC   45cm  100
## 5         100 202211 GCrew  TR      LysC   45cm  100
## 6         100 202211 GCrew  TR      LysC   45cm  100
##           Sample_ID   S04_ppm   Cl_ppm   S04_mM   Cl_mM
## 1 1_202211_MSM_UP_LysA_10cm 264.6481 3270.8644 8.254775 92.266979
## 2 10_202211_MSM_TR_LysA_20cm 1325.3953 6875.5636 41.341089 193.951018
## 3 100_202211_GCrew_TR_LysC_45cm 158.7476 134.0034 4.951578 3.780068
## 4 100_202211_GCrew_TR_LysC_45cm 158.7476 133.2425 4.951578 3.758604
## 5 100_202211_GCrew_TR_LysC_45cm 19.3955 134.0034 0.604975 3.780068
## 6 100_202211_GCrew_TR_LysC_45cm 19.3955 133.2425 0.604975 3.758604
##   salinity
## 1 5.9104780
## 2 12.4241694
## 3 0.2421701
## 4 0.2407952
## 5 0.2421701
## 6 0.2407952
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_S04_202211.csv") #Change fil
```

END