

Synoptic CB: Porewater DIC

November 2025 Samples

2026-05-28

Contents

0.1	Import Data Functions	2
0.2	Import Sample Data	2
0.3	Assessing Standard Curves	3
0.4	CRM Check - No CRMs on this run	4
0.5	Assess Check Standards	4
0.6	Assess Blanks	5
0.7	Assess Duplicates - No duplicates included on this run	6
0.8	Sample Flagging - Are samples Within the range of the curve?	6
0.9	Visualize Data by Plot	7
0.10	Convert data from mg/L to uMoles/L	8
0.11	Check to see if samples run match metadata & merge info	8
0.12	Export Processed Data	8

```
##Setup - Change things here & write any notes
```

```
#identify section  
cat("Setup Information")
```

```
## Setup Information
```

```
##### Run information - PLEASE CHANGE  
Date_Run = "05/26/25" #Date that instrument was run  
Run_by = "Melanie Giessner" #Instrument user  
Script_run_by = "Stephanie J. Wilson" #Code user  
run_notes = "Check Standards are out of range. Standard curve looks good - going to accept run." #an  
samples <- c("GCW", "GWI", "MSM", "SWH") #whatever identifies your samples within the same names  
samples_pattern <- paste(samples, collapse = "|")  
  #samples_pattern <- "GCW" #use this instead of the line above if you have only one site code  
chks_name = "Chk_Std_" #what did you name your check standards?  
crm_name = "CRM|crm" #what did you name your CRMS?  
  
##### File Names - PLEASE CHANGE  
#file path and name for raw summary data file  
raw_file_name = "Raw Data/COMPASS_SynopticCB_PW_DIC_202605.txt"  
  
#file path and name for raw all peaks file  
raw_allpeaks_name = "Raw Data/COMPASS_SynopticCB_PW_DIC_202605_allpeaks.txt"  
  
#file path and name of processed data file  
processed_file_name = "Processed Data/COMPASS_SynopticCB_PW_Processed_DIC_202605.csv"  
  
##### Log Files - PLEASE CHECK  
#downloaded metadata csv - downloaded from Google drive as csv for this year  
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2026.csv"  
  
#qaqc log file path for this year  
Log_path = "Raw Data/COMPASS_Synoptic_DIC_QAQClog_2026.csv"
```

```
##Set Up Code
```

```
##Read in metadata and create similar sample IDs for matching to samples
```

0.1 Import Data Functions

0.2 Import Sample Data

```
## Import Sample Data
```

```
## New names:  
## * `` -> `...14`
```


```
## # A tibble: 6 x 3  
##   sample_name      ic_raw run_datetime  
##   <chr>          <dbl> <chr>
```

```
## 1 202605_GCW_TR_LysB_10cm 5.50 5/26/2026 4:35:34 PM
## 2 202605_GCW_WC_LysA_10cm 60.0 5/26/2026 4:50:50 PM
## 3 202605_GCW_WC_LysA_20cm 90.8 5/26/2026 5:06:51 PM
## 4 202605_GCW_WC_LysA_45cm 198. 5/26/2026 5:23:08 PM
## 5 202605_GCW_WC_LysB_10cm 72.6 5/26/2026 5:38:25 PM
## 6 202605_GCW_WC_LysB_20cm 89.0 5/26/2026 5:54:25 PM
```

0.3 Assessing Standard Curves


```
## Assess the Standard Curves
```

```
## New names:
## `geom_smooth()` using formula = 'y ~ x'
## * `` -> `...18`
```



```
## Warning: Removed 4 rows containing missing values or values outside the scale range
## (`geom_point()`).
```

```
## Warning: Removed 4 rows containing missing values or values outside the scale range
## (`geom_line()`).
```




```
## [1] "IC Curve r2 GOOD"
```

0.4 CRM Check - No CRMs on this run

0.5 Assess Check Standards

```
## Assess the Check Standards
```

```
## New names:
## * `` -> `...14`
```




```
## [1] ">60% of IC Check Standards are within range of expected concentration"
```

0.6 Assess Blanks

```
## Assess Blanks
```

```
## New names:
## * `` -> `...14`
```

```
## [1] ">60% of Carbon Blank concentrations are lower 25% quartile of samples"
```



carbon blanks:


[1] 1.313655

0.7 Assess Duplicates - No duplicates included on this run

0.8 Sample Flagging - Are samples Within the range of the curve?

Sample Flagging


C: Grey = Within Range of Curve



0.9 Visualize Data by Plot

```
## Visualize Data
```

```
## Warning in rbind(c("202605", "GCW", "TR", "LysB", "10cm"), c("202605", "GCW", :
## number of columns of result is not a multiple of vector length (arg 11)
```



0.10 Convert data from mg/L to uMoles/L

0.11 Check to see if samples run match metadata & merge info

```
## Check Sample IDs with Metadata
```

```
## All sample IDs are present in metadata.
```

0.12 Export Processed Data

```
## Export Processed Data
```

```
## # A tibble: 6 x 18
```

```
## Project      Region Site Zone Replicate Depth_cm Sample_ID Year Month Day
## <chr>         <chr> <chr> <fct> <chr>          <int> <chr>    <int> <int> <int>
## 1 COMPASS: Sy~ CB      MSM UP      A              10 202605_M~ 2026      5 14
## 2 COMPASS: Sy~ CB      MSM UP      B              10 202605_M~ 2026      5 14
## 3 COMPASS: Sy~ CB      MSM UP      C              10 202605_M~ 2026      5 14
## 4 COMPASS: Sy~ CB      SWH UP      A              20 202605_S~ 2026      5 19
## 5 COMPASS: Sy~ CB      SWH UP      A              45 202605_S~ 2026      5 19
## 6 COMPASS: Sy~ CB      SWH UP      B              20 202605_S~ 2026      5 19
## # i 8 more variables: Time <chr>, Time_Zone <chr>, ic_mgL <dbl>, ic_uM <dbl>,
## # ic_flag <chr>, Analysis_runtime <chr>, Run_notes <chr>, Field_notes <lg1>
```

```
#end
```