

Synoptic_CB_Nutrients_2023_Sept

September

2025-07-01

Contents

0.1	Import Data & Clean	2
0.2	Assessing standard Curves	3
0.3	TEMPORARY CHUNK FOR SEPT 2023 - Run finished then restarted for dilutions	7
0.4	Dilution Corrections - ensure the latest dilution is kept	7
0.5	Performance Check	7
0.6	Analyze the Check Standards	9
0.7	Analyze Blanks	10
0.8	Analyze Duplicates	11
0.9	Spikes	12
0.10	Matrix Effects	13
0.11	Unit Converted Data Column Added (mg/L to uM)	13
0.12	Sample Flagging - Within range of standard curve	13
0.13	Pull out sample id information	13
0.14	Pulling Rhizon Samples	13
0.15	Check to see if samples run match metadata & merge info	14
0.16	Visualize Data	15
0.17	Export Processed Data	17

##Run Information

```
cat("Run Information: Input by User") #lets you know what section you're in
```

Run Information: Input by User

```
#set the run date & user name
run_date <- "20240125"
sample_year <- 2023
sample_month <- 09
user <- "Stephanie Wilson"
```

```
#identify the files you want to read in
```

```
#read in as a list to accomodate ultiple runs in a month
```

```
NOx_files <- c("Raw Data/SEAL_COMPASS_Synoptic_NOx_Sept2023_1.csv",
              "Raw Data/SEAL_COMPASS_Synoptic_NOx_Sept2023_2.csv",
              "Raw Data/SEAL_COMPASS_Synoptic_NOx_Sept2023_3.csv")
NH3_PO4_files <- c("Raw Data/SEAL_COMPASS_Synoptic_NH3_PO4_Sept2023_1.csv",
                  "Raw Data/SEAL_COMPASS_Synoptic_NH3_PO4_Sept2023_2.csv",
                  "Raw Data/SEAL_COMPASS_Synoptic_NH3_PO4_Sept2023_3.csv")
```

```
# Define the file path for QAQC log file - NO Need to change just check year
```

```
file_path <- "Raw Data/SEAL_COMPASS_Synoptic_QAQC_Log_2023.csv"
final_path <- "Processed Data/COMPASS_Synoptic_Nutrients_202309.xlsx"
```

```
#record any notes about the run or anything other info here:
```

```
run_notes <- "NOx Blanks are out of range because of the slightly negative y intercept of the standard"
```

```
#Set up file path for metadata
```

```
#downloaded metadata csv - downloaded from Google drive as csv for this year
```

```
#https://docs.google.com/spreadsheets/d/1HCAN0\_q6y17x0RUXVzID09hVal-RfwWc/edit?usp=sharing&ouid=10899
```

```
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2023.csv"
```

##Setup

##Read in metadata and create similar sample IDs for matching to samples

0.1 Import Data & Clean


0.2 Assessing standard Curves

```
#Pull out standards data
```

```
## Assess Standard Curves
```


```
#Plot standards data
```

```
## Assess Standard Curves
```




```
## 'geom_smooth()' using formula = 'y ~ x'
```

NH3 Standard Curve



```
## 'geom_smooth()' using formula = 'y ~ x'
```

PO4 Standard Curve



```
## [1] "NOx Curve r2 GOOD - PROCEED"
```

```
## [1] "NH3 Curve r2 GOOD - PROCEED"
```

```
## [1] "PO4 Curve r2 GOOD - PROCEED"
```

```
## [1] "QAQC log file exists and has been read into the code."
```




Table 1: Average Slope by Analyte

Test	avg_slope
NH3	1.899
NOx	1.261
PO4	2.353

0.3 TEMPORARY CHUNK FOR SEPT 2023 - Run finished then restarted for dilutions

0.4 Dilution Corrections - ensure the latest dilution is kept

Dilution Corrections

Duplicated samples: MSM_202309_TR_LysA_20cm, MSM_202309_TR_LysB_10cm, MSM_202309_TR_LysB_20cm, MSM_202309_TR_LysB_30cm

Duplicated samples with NO dilution present (possible input error): MSM_202309_UP_LysB_10cm

0.5 Performance Check

[1] "NOx pe Check has a % Difference <10% - PROCEED"

Run mean = 0.738716

Expected = 0.706

[1] "NH3 pe Check has a % Difference >10% - REASSESS"

Run mean = 1.531811

Expected = 0.948


[1] "PO4 pe Check has a % Difference <10% - PROCEED"

Run mean = 0.864723

Expected = 0.818

#Check NOx Reduction Efficiency

Assess Reduction Efficiency



```
## [1] "Mean NOx Reduction Efficiency >95% - PROCEED"
```

```
## [1] 98.47613
```



0.6 Analyze the Check Standards

```
## Analyze Check Standards
```

```
## [1] "NOx Check Standard RSD within Range - PROCEED"
```

```
## [1] "NH3 Check Standard RSD within Range - PROCEED"
```

```
## [1] "PO4 Check Standard RSD within Range - PROCEED"
```



```
## [1] ">60% of NOx Check Standards are within range of expected concentration - PROCEED"
```

```
## [1] ">60% of NH3 Check Standards are within range of expected concentration - PROCEED"
```

```
## [1] ">60% of PO4 Check Standards are within range of expected concentration - PROCEED"
```

0.7 Analyze Blanks

```
## Assess Blanks
```

```
## [1] "<60% of NOx blaks are lower 25% quartile of samples - REASSESS"
```

```
## [1] ">60% of NH3 Blank concentrations are lower than the lower 25% quartile of samples - PROCEED"
```

```
## [1] ">60% of PO4 Blank concentrations are lower than the lower 25% quartile of samples- PROCEED"
```

```
## [1] "NOx pe Check has a % Difference <10% - PROCEED"
```

```
## Run mean = 0.738716
```

```
## Expected = 0.706
```




Table 2: Mean Concentration of Blanks

Test	Blank_Mean_Conc
NOx	-0.0176
NH3	-0.0022
PO4	0.0023

0.8 Analyze Duplicates


```
## Analyze Duplicates

## [1] "<60% of NOx Duplicates have a CV <10% - REASSESS"

## [1] ">60% of NH3 Duplicates have a CV <10% - PROCEED"

## [1] ">60% of PO4 Duplicates have a CV <10% - PROCEED"

## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```




0.9 Spikes

```
## [1] ">60% of Spikes have a CV <50% - PROCEED"
```

```
## [1] ">60% of Spikes have a CV <50% - PROCEED"
```

```
## [1] "Flag not calculated for this run"
```



0.10 Matrix Effects

```
## [1] "NO NOx Matrix Effect, PROCEED"
```

```
## [1] "NO NH3 Matrix Effect, PROCEED"
```

```
## [1] "NO PO4 Matrix Effect, PROCEED"
```

0.11 Unit Converted Data Column Added (mg/L to uM)

0.12 Sample Flagging - Within range of standard curve

```
## Sample Flagging
```

0.13 Pull out sample id information

```
## Sample Processing
```

```
## Warning: Expected 5 pieces. Missing pieces filled with 'NA' in 36 rows [18, 19, 20, 42,
## 43, 44, 77, 78, 79, 101, 102, 103, 121, 122, 123, 138, 139, 140, 157, 158,
## ...].
```

```
## Warning: There was 1 warning in 'mutate()'.
## i In argument: 'Samp_Time = ym(Samp_Time)'.
## Caused by warning:
## ! 10 failed to parse.
```

0.14 Pulling Rhizon Samples

```
# Filter rhizon and peeper samples
df_rhizon <- df_all %>%
  filter(str_detect(Sample_Name, "RHZ"))
df_peep <- df_all %>%
  filter(str_detect(Sample_Name, "PPR"))

# Timestamp for backups
timestamp <- format(Sys.time(), "%Y-%m-%d_%H%M")

# Paths
folder_path <- file.path("Raw Data", "Rhizon+Peeper")
dir.create(folder_path, recursive = TRUE, showWarnings = FALSE)

rhizon_main <- file.path(folder_path, "rhizon_data.csv")
peeper_main <- file.path(folder_path, "peeper_data.csv")

rhizon_backup <- file.path(folder_path, paste0("rhizon_data_", timestamp, ".csv"))
peeper_backup <- file.path(folder_path, paste0("peeper_data_", timestamp, ".csv"))

# Write timestamped backups
```

```

write.csv(df_rhizon, rhizon_backup, row.names = FALSE)
write.csv(df_peep, peeper_backup, row.names = FALSE)

# Overwrite the main files with latest data
write.csv(df_rhizon, rhizon_main, row.names = FALSE)
write.csv(df_peep, peeper_main, row.names = FALSE)

## ^^ I think there is a cleaner way to write this out, but this should work for now ^^

```

0.15 Check to see if samples run match metadata & merge info

```
## Check Sample IDs with Metadata
```

```
## Some sample IDs are missing from metadata.
```

```

## [1] "GCW_202039_SW_A"          "MSM_202309_TR_LysC_45cm"
## [3] "SWH_202309_UPCON_LysA_10cm" "SWH_202309_UPCON_LysA_20cm"
## [5] "SWH_202309_UPCON_LysA_45cm" "SWH_202309_UPCON_LysB_20cm"
## [7] "SWH_202309_UPCON_LysB_45cm" "SWH_202309_UPCON_LysC_10cm"
## [9] "SWH_202309_UPCON_LysC_20cm" "SWH_202309_UPCON_LysC_45cm"
## [11] "SWH_202039_TR_LysC_45cm"   "SWH_202309_WC_LysA_45cm"
## [13] "SWH_202309_WC_LysC_10cm"   "SWH_202309_WC_LysC_20cm"
## [15] "GWI_202039_WC_LysA_20cm"   "GWI_202039_WC_LysB_45cm"
## [17] "MSM_202039_WC_LysC_10cm"   "MSM_202039_WC_LysC_20cm"
## [19] "SWH_202039_TR_LysB_45cm"

```

```

## Warning: Expected 5 pieces. Missing pieces filled with 'NA' in 36 rows [18, 19, 20, 42,
## 43, 44, 77, 78, 79, 101, 102, 103, 121, 122, 123, 138, 139, 140, 157, 158,
## ...].

```

```


## Warning: There was 1 warning in 'mutate()'.
## i In argument: 'Samp_Time = ym(Samp_Time)'.
## Caused by warning:
## ! 10 failed to parse.

```


Visualize Data



Porewater NH3



Porewater PO4



0.17 Export Processed Data

#end