

COMPASS_SynopticCB_PW_SO4_Cl_202311_Template

Stephanie J. Wilson

2025-08-28

Contents

0.1	Import Sample Data	3
0.2	Assessing Standard Curves	3
0.3	Assess Check Standards	6
0.4	Assess Blanks	7
0.5	Assess Duplicates	8
0.6	Calculate mmol/L concentrations & salinity, add dilutions	10
0.7	Assess Analytical Spikes	11
0.8	Check if samples Within the range of the standard curve	11
0.9	Check to see if samples run match metadata & merge info	12
0.10	Visualize Data by Plot	12
0.11	Export Processed Data	12

```
##Setup - Change things here & write any notes
```

```
#identify section  
cat("Setup Information")
```

```
## Setup Information
```

```
##### Run information - PLEASE CHANGE
```

```
Date_Run = "2023-09-05" #Date that instrument was run  
Run_by = "Unknown" #Instrument user  
Script_run_by = "Zoe Read" #Code user  
run_notes = "Samples missing from metadata: MSM_202311_UP_LYSC_45CM MSM_202311_WC_LYSB_45CM MSM_202311_WC_RHZ_LYSC. 118_MSM_202311_TR_RHZ_SF_Tree_1 value above cal curve for S04, but only slightly below" #Notes on run  
  
samples <- c("GCW", "GWI", "MSM", "SWH") #whatever identifies your samples within the same names  
samples_pattern <- paste(samples, collapse = "|")  
  #samples_pattern <- "GCW" #use this instead of the line above if you have only one site code  
chks_name = "Check Standard" #what did you name your check standards?
```

```
##### File Names - PLEASE CHANGE
```

```
#file path and name for raw summary data file  
  # raw_file_name_cl = "Porewater/Sulfate_Chloride/Synoptic_CB/2023/Raw Data/COMPASS_Synoptic_CB_MonMon_202311_C1.txt"  
  # raw_file_name_so4 = "Porewater/Sulfate_Chloride/Synoptic_CB/2023/Raw Data/COMPASS_Synoptic_CB_MonMon_202311_S04.txt"
```

```
raw_file_name_cl = "Raw Data/COMPASS_Synoptic_CB_MonMon_202311_C1.txt"  
raw_file_name_so4 = "Raw Data/COMPASS_Synoptic_CB_MonMon_202311_S04.txt"
```

```
#file path and name of processed data file
```

```
  # processed_file_name = "Porewater/Sulfate_Chloride/Synoptic_CB/2023/Processed Data/COMPASS_SynopticCB_PW_Processed_C1_S04_202311.csv"  
processed_file_name = "Processed Data/COMPASS_SynopticCB_PW_Processed_C1_S04_202311.csv"
```

```
##### Log Files - PLEASE CHECK
```

```
#downloaded metadata csv - downloaded from Google drive as csv for this year  
  # Raw_Metadata = "Porewater/Sulfate_Chloride/Synoptic_CB/2023/Raw Data/COMPASS_SynopticCB_PW_SampleLog_2023.csv"  
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2023.csv"
```

```
#qac log file path for this year
```

```
  # Log_path = "Porewater/Sulfate_Chloride/Synoptic_CB/2023/Raw Data/COMPASS_Synoptic_C1_S04_QACLog_2023.csv"  
Log_path = "Raw Data/COMPASS_Synoptic_C1_S04_QACLog_2023.csv"
```

```
##Add Required Packages
```

```
##Set Up Code
```

```
#Link to the protocol used for analysis  
  #steph will add this soon
```

```
#Any coefficients / constants that are needed for calculations
```

```
cl_mw <- 35.45 #molecular weight of Chloride: 35.45  
s_mw <- 32.06 #molecular weight of sulfur: 32.06
```

```
#Flags for Dionex
```

```
r2_cutoff = 0.98 #this is the level below which we want to rerun or consider a curve
```

```

chk_flag_std_s = 10 #this is the maximum cv allowed for sulfate check standards
chk_flag_std_cl = 5 #this is the maximum cv allowed for chloride check standards
chk_flag_dups = 10 #this is the maximum cv allowed for duplicates
chks_flag = 0.80 #this is the percent of chks we want to have a CV less than the max allowable cv
#blank_flag - calculated based on samples later in this code as lower 25% quantile of sample concentr

#Standard concentrations - Update if running different standard curve:
standards <- tibble(
  sample_ID = c("Standard 1", "Standard 2", "Standard 3", "Standard 4", "Standard 5"),
  SO4_std_conc = c(0.5, 1.0, 2.0, 10, 20),
  Cl_std_conc = c(5, 10, 20, 100, 200))

#Spike concentration calc
spk_std <- (250/s_mw) # in mM
spkvol <- 10 # in uL
spkvol <- spkvol/1000000 # in L
#spike for these samples was 10uL of the 250 µg/mL standard
spk_Conc <- (spk_std)*spkvol # mmoles of SO4

#Top standard Concentrations- Update if running different standard curve:
top_std_cl = 200
top_std_so4 = 20

#Set time zone
common_tz = "Etc/GMT+5"
Sys.setenv(TZ = "America/New_York")

#plot indicators
site_order <- c('GCW', 'MSM', 'GWI', 'SWH')
plot_order <- c('UP', 'SWAMP', 'TR', 'WC', 'SW')
plot_colors <- c("#20063B", "darkgrey", "#FFBC42", "#419973", "#25ABE6" )

```

```

##Read in metadata and create similar sample IDs for matching to samples

```

0.1 Import Sample Data

0.2 Assessing Standard Curves

```

## Assess the Standard Curves


```

```

## 'geom_smooth()' using formula = 'y ~ x'


```

Chloride Std Curve




```
## 'geom_smooth()' using formula = 'y ~ x'
```

Sulfate Std Curve



```
## 'geom_line()': Each group consists of only one observation.  
## i Do you need to adjust the group aesthetic?
```




```
## [1] "Cl Curve r2 GOOD"
```

```
## [1] "SO4 Curve r2 GOOD"
```

0.3 Assess Check Standards

```
## Assess the Check Standards
```



```
## [1] ">80% of Chloride Check Standards are within range of expected concentration"
```


```
## [1] ">80% of Sulfate Check Standards are within range of expected concentration"
```

0.4 Assess Blanks

```
## Assess Blanks
```

```
## [1] ">80% of Chloride Blank concentrations are lower 25% quartile of samples"
```

```
## [1] ">80% of Sulfate Blank concentrations are lower 25% quartile of samples"
```



```
## Chloride blanks mean ppm:
```



```
## [1] 0.01699524
```

```
## Sulfate blanks mean ppm:
```

```
## [1] 0.004957143
```


0.5 Assess Duplicates

```
## Assess Duplicates
```

```
## [1] ">80% of Chloride Duplicates have a CV <10%"
```

```
## [1] ">80% of Sulfate Duplicates have a CV <10%"
```



```
## [1] ">80% of Chloride Duplicates have a percent difference <10%"
```

```
## [1] ">80% of Sulfate Duplicates have a percent difference <10%"
```

0.6 Calculate mmol/L concentrations & salinity, add dilutions

```
cat("Unit Conversion and Salinity Calculation")

# Convert ppm to mmol/L
all_dat$S04_mM <- (all_dat$S04_ppm / s_mw)
all_dat$Cl_mM <- (all_dat$Cl_ppm / cl_mw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# = ((1.807*Cl_ppm)+0.026)/1000
all_dat$salinity <- ((1.8070 * all_dat$Cl_ppm) + 0.026) / 1000

#Need to determine dilution factors for your samples
#for Steph / COMPASS this depends on the site so...
all_dat$Dilution <- 1
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "UP"),
```

```

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "SW"),

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "UP"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "SW"),

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "UP"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "SW"),

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "SWH"), 50, all_dat$Dilution)


head(all_dat)

```

0.7 Assess Analytical Spikes

0.8 Check if samples Within the range of the standard curve


Sample Flagging



0.9 Check to see if samples run match metadata & merge info

0.10 Visualize Data by Plot

Visualize Data



0.11 Export Processed Data

#end