

COMPASS_Synoptic_SEAL_Data_Analysis_Oct2022

Stephanie J. Wilson

2023-03-03

Information

```
#####
##### COMPASS Synoptic
##### Chesapeake Bay Sites
##### Data Analysis Code: Porewater Nutrients
##### MONTH: August 2022
#####

#####
##### Information #####
#Author: Stephanie J. Wilson
#Edited: 20220727

#Samples taken from Lysimeters & Sippers at CB Synoptic Sites
#Samples were filtered with 0.45 uM, kept on ice, frozen until analysis
#Field Protocol:
#Samples Analyzed on a SEAL discrete auto analyzer
#Lab Protocol: https://docs.google.com/document/d/1VaJT7Wb9AcdmM1tgsR_9ZtQ6kwcaoNmp/edit?usp=sharing&ou
#NOx method = https://drive.google.com/file/d/1sicqBFnzVxmDd5I2_pu8s8pj7iNUAhF/view?usp=sharing
#NH4 method = https://drive.google.com/file/d/1ENGemUEvm_rffZqv3lz9BjD0pAMX5nzu/view?usp=sharing
#PO4 method = https://drive.google.com/file/d/1m3gXDZnJoIo_QmyhvZG4HRgGShCzm9Wq/view?usp=sharing
#Units from SEAL = mg/L and converted to uMoles/L
```

QAQC on Slopes

```
library(ggplot2)

## Warning: package 'ggplot2' was built under R version 4.3.3

library(data.table)

#read in datafile with all the slopes
qlog <- read.csv("Raw Data/SERC_SEAL_STDS_Log.csv")
head(qlog)
```

##	Date	User	Machine	Analysis	Slope	Intercept	R2
----	------	------	---------	----------	-------	-----------	----

```

## 1 20220513 Stephanie Wilson SERC SEAL      NH3 1.9370 -0.030 0.9995
## 2 20220513 Stephanie Wilson SERC SEAL      P04 2.8690  0.002 0.9994
## 3 20220525 Stephanie Wilson SERC SEAL      V-Nox 1.0000  0.000 0.9990
## 4 20220614 Stephanie Wilson SERC SEAL      NH3 1.1487 -0.011 0.9999
## 5 20220614 Stephanie Wilson SERC SEAL      P04 2.3950 -0.001 0.9998
## 6 20220615 Stephanie Wilson SERC SEAL      V-Nox 1.0000  0.000 0.9998
##   Nox_Red_Eff X
## 1      NA
## 2      NA
## 3      NA
## 4      NA
## 5      NA
## 6      NA

```

```

#pull out each method
qNH3 <- qlog[qlog$Analysis %like% "NH3", ]
head(qNH3)

```

```

##           Date       User Machine Analysis Slope Intercept     R2
## 1 20220513 Stephanie Wilson SERC SEAL      NH3 1.9370 -0.030 0.9995
## 4 20220614 Stephanie Wilson SERC SEAL      NH3 1.1487 -0.011 0.9999
## 7 20220615 Stephanie Wilson SERC SEAL      NH3 1.7070 -0.008 1.0000
## 10 20220718 Stephanie Wilson SERC SEAL      NH3 1.7080 -0.017 0.9995
## 12 20220718 Stephanie Wilson SERC SEAL      NH3 1.6590 -0.017 0.9997
## 14 20220719 Stephanie Wilson SERC SEAL      NH3 1.5560  0.000 0.9994
##   Nox_Red_Eff X
## 1      NA
## 4      NA
## 7      NA
## 10     NA
## 12     NA
## 14     NA

```

```

qP04 <- qlog[qlog$Analysis %like% "P04", ]
head(qP04)

```

```

##           Date       User Machine Analysis Slope Intercept     R2
## 2 20220513 Stephanie Wilson SERC SEAL      P04 2.869   0.002 0.9994
## 5 20220614 Stephanie Wilson SERC SEAL      P04 2.395 -0.001 0.9998
## 8 20220615 Stephanie Wilson SERC SEAL      P04 2.349   0.000 0.9999
## 11 20220718 Stephanie Wilson SERC SEAL      P04 2.445 -0.003 0.9993
## 13 20220718 Stephanie Wilson SERC SEAL      P04 2.349 -0.004 0.9992
## 15 20220719 Stephanie Wilson SERC SEAL      P04 2.271 -0.003 0.9993
##   Nox_Red_Eff X
## 2      NA
## 5      NA
## 8      NA
## 11     NA
## 13     NA
## 15     NA

```

```

qN0x <- qlog[qlog$Analysis %like% "V-Nox", ]
head(qN0x)


##           Date      User Machine Analysis Slope Intercept     R2
## 3 20220525 Stephanie Wilson SERC SEAL V-Nox     1        0 0.9990
## 6 20220615 Stephanie Wilson SERC SEAL V-Nox     1        0 0.9998
## 9 20220718 Stephanie Wilson SERC SEAL V-Nox     1        0 0.9999
## 16 20220730 Stephanie Wilson SERC SEAL V-Nox     1        0 1.0000
## 17 20220731 Stephanie Wilson SERC SEAL V-Nox     1        0 1.0000
## 22 20220824 Stephanie Wilson SERC SEAL V-Nox     1        0 0.9988
##   Nox_Red_Eff X
## 3          NA
## 6          NA
## 9          NA
## 16         NA
## 17         NA
## 22         NA

##### NH3
#plot the slopes to make sure there are no crazy outliers
slope1 <- ggplot(data=qNH3, aes(x=Date, y=Slope)) +
  geom_line() +
  geom_point(aes(size=3)) +
  theme_classic() + ylim(0, 5) +
  theme(legend.position="none") +
  ggtitle("NH3 Slopes")

slope1

## Warning: Removed 1 row containing missing values or values outside the scale range
## ('geom_point()').


```



```
#plot the intercepts to make sure there are no crazy outliers
int1 <- ggplot(data=qNH3, aes(x=Date, y=Intercept)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(-0.5,0.5) +
  theme(legend.position="none")+
  ggtitle("NH3 Intercepts")

int1

## Warning: Removed 1 row containing missing values or values outside the scale range
## ('geom_point()').
```




```
#plot the R2s to make sure there are no crazy outliers
Rsq1 <- ggplot(data=qNH3, aes(x=Date, y=R2)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(0.75, 1.25) +
  theme(legend.position="none")+
  ggtitle("NH3 R2s")
```

```
Rsq1
```

```
## Warning: Removed 39 rows containing missing values or values outside the scale range
## ('geom_point()').
```


NH3 R2s



```
##### P04
slope2 <- ggplot(data=qNH3, aes(x=Date, y=Slope)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(0, 5) +
  theme(legend.position="none")+
  ggtitle("P04 Slopes")

slope2
```


```
## Warning: Removed 1 row containing missing values or values outside the scale range
## ('geom_point()').
```



```
#plot the intercepts to make sure there are no crazy outliers
int2 <- ggplot(data=qNH3, aes(x=Date, y=Intercept)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(-0.5,0.5) +
  theme(legend.position="none")+
  ggtitle("P04 Intercepts")

int2

## Warning: Removed 1 row containing missing values or values outside the scale range
## ('geom_point()').
```




```
#plot the R2s to make sure there are no crazy outliers
Rsq2 <- ggplot(data=qNH3, aes(x=Date, y=R2)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(0.75, 1.25) +
  theme(legend.position="none")+
  ggtitle("P04 R2s")
```

```
Rsq2
```

```
## Warning: Removed 39 rows containing missing values or values outside the scale range
## ('geom_point()').
```


PO4 R2s



```
##### NOx
slope3 <- ggplot(data=qNH3, aes(x=Date, y=Slope)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(0, 5) +
  theme(legend.position="none")+
  ggtitle("NOx Slopes")

slope3
```


```
## Warning: Removed 1 row containing missing values or values outside the scale range
## ('geom_point()').
```



```
#plot the intercepts to make sure there are no crazy outliers
int3 <- ggplot(data=qNH3, aes(x=Date, y=Intercept)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(-0.5,0.5) +
  theme(legend.position="none")+
  ggtitle("NOx Intercepts")

int3

## Warning: Removed 1 row containing missing values or values outside the scale range
## ('geom_point()').
```




```
#plot the R2s to make sure there are no crazy outliers
Rsq3 <- ggplot(data=qNH3, aes(x=Date, y=R2)) +
  #geom_line()+
  geom_point(aes(size=3)) +
  theme_classic() + ylim(0.75, 1.25) +
  theme(legend.position="none")+
  ggtitle("NOx R2s")
```

```
Rsq3
```

```
## Warning: Removed 39 rows containing missing values or values outside the scale range
## ('geom_point()').
```

NOx R2s



Code Set up

```
#packages:  
library(ggplot2)  
library(dplyr)  
  
##  
## Attaching package: 'dplyr'  
  
## The following objects are masked from 'package:data.table':  
##  
##     between, first, last  
  
## The following objects are masked from 'package:stats':  
##  
##     filter, lag  
  
## The following objects are masked from 'package:base':  
##  
##     intersect, setdiff, setequal, union
```

```

library(data.table)
library(matrixStats)

## 
## Attaching package: 'matrixStats'

## The following object is masked from 'package:dplyr':
## 
##     count

library(gridExtra)

## 
## Attaching package: 'gridExtra'

## The following object is masked from 'package:dplyr':
## 
##     combine

library(ggpubr)
library(grid)

```

Ammonia & Phosphate

```

#read in data
file1 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NH3_P04_202210_1.csv")
#Quick look at dataframe
head(file1)

##   RUNSTARTED X1674009384 X1.17.2023.21.36          X X.1      X.2      X.3
## 1      RESULT      -1           S1 Standard 1  0 0.013208 0.013208
## 2      RESULT      -2           S90 Standard .0389  1 0.033938 0.033938
## 3      RESULT      -2           S91 Standard .1000  2 0.064862 0.064862
## 4      RESULT      -2           S92 Standard .2000  3 0.118781 0.118781
## 5      RESULT      -2           S93 Standard .5000  4 0.274742 0.274742
## 6      RESULT      -2           S94 Standard 1.0000  5 0.536226 0.536226
##   X.4 X.5 X.6 X.7    X.8      X.9      X.10      X.11
## 1  0  0  0  0 mg N/L Ammonia 2 1674010634 1/17/2023 21:57
## 2  0  0  0  0 mg N/L Ammonia 2 1674010794 1/17/2023 21:59
## 3  0  0  0  0 mg N/L Ammonia 2 1674010954 1/17/2023 22:02
## 4  0  0  0  0 mg N/L Ammonia 2 1674011114 1/17/2023 22:05
## 5  0  0  0  0 mg N/L Ammonia 2 1674011274 1/17/2023 22:07
## 6  0  0  0  0 mg N/L Ammonia 2 1674011434 1/17/2023 22:10

#take out only the columns that we need
dat1 <- file1[,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(dat1) <- c('Run_Info', 'Sample_Name', 'Conc', "Abs", "Units", "Test")
head(dat1)

```

```

##   Run_Info      Sample_Name      Conc      Abs  Units      Test
## 1    RESULT      Standard 1 0.013208 0.013208 mg N/L Ammonia 2
## 2    RESULT      Standard .0389 0.033938 0.033938 mg N/L Ammonia 2
## 3    RESULT      Standard .1000 0.064862 0.064862 mg N/L Ammonia 2
## 4    RESULT      Standard .2000 0.118781 0.118781 mg N/L Ammonia 2
## 5    RESULT      Standard .5000 0.274742 0.274742 mg N/L Ammonia 2
## 6    RESULT      Standard 1.0000 0.536226 0.536226 mg N/L Ammonia 2

#read in data
file2 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NH3_P04_202210_2.csv")
#Quick look at dataframe
head(file2)

##   RUNSTARTED X1674070636 X1.18.2023.14.37          X X.1      X.2      X.3
## 1    RESULT      -1           S1      Standard 1 0 0.015211 0.015211
## 2    RESULT      -2           S90     Standard .0389 1 0.034194 0.034194
## 3    RESULT      -2           S91     Standard .1000 2 0.064507 0.064507
## 4    RESULT      -2           S92     Standard .2000 3 0.117400 0.117400
## 5    RESULT      -2           S93     Standard .5000 4 0.265003 0.265003
## 6    RESULT      -2           S94     Standard 1.0000 5 0.536287 0.536287
##   X.4 X.5 X.6 X.7      X.8      X.9      X.10          X.11
## 1  0   0   0   0      0 mg N/L Ammonia 2 1674071922 1/18/2023 14:58
## 2  0   0   0   0      0 mg N/L Ammonia 2 1674072082 1/18/2023 15:01
## 3  0   0   0   0      0 mg N/L Ammonia 2 1674072242 1/18/2023 15:04
## 4  0   0   0   0      0 mg N/L Ammonia 2 1674072402 1/18/2023 15:06
## 5  0   0   0   0      0 mg N/L Ammonia 2 1674072562 1/18/2023 15:09
## 6  0   0   0   0      0 mg N/L Ammonia 2 1674072722 1/18/2023 15:12

#take out only the columns that we need
dat2 <- file2[ ,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(dat2) <- c('Run_Info','Sample_Name','Conc', "Abs", "Units", "Test")
head(dat2)

##   Run_Info      Sample_Name      Conc      Abs  Units      Test
## 1    RESULT      Standard 1 0.015211 0.015211 mg N/L Ammonia 2
## 2    RESULT      Standard .0389 0.034194 0.034194 mg N/L Ammonia 2
## 3    RESULT      Standard .1000 0.064507 0.064507 mg N/L Ammonia 2
## 4    RESULT      Standard .2000 0.117400 0.117400 mg N/L Ammonia 2
## 5    RESULT      Standard .5000 0.265003 0.265003 mg N/L Ammonia 2
## 6    RESULT      Standard 1.0000 0.536287 0.536287 mg N/L Ammonia 2

#read in data
file3 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NH3_P04_202210_3.csv")
#Quick look at dataframe
head(file3)

##   RUNSTARTED X1674261032 X1.20.2023.19.30          X X.1      X.2      X.3
## 1 RUNSTARTED 1674261064 1/20/2023 19:31          NA      NA      NA
## 2 RUNSTARTED 1674270611 1/20/2023 22:10          NA      NA      NA
## 3    RESULT      -1           S1      Standard 1 0 0.015424 0.015424
## 4    RESULT      -2           S90     Standard .0389 1 0.037041 0.037041

```

```

## 5      RESULT          -2           S91 Standard .1000   2 0.065641 0.065641
## 6      RESULT          -2           S92 Standard .2000   3 0.120285 0.120285
## X.4 X.5 X.6 X.7     X.8       X.9       X.10        X.11
## 1  NA  NA  NA  NA
## 2  NA  NA  NA  NA
## 3  0   0   0   0 mg N/L Ammonia 2 1674271855 1/20/2023 22:30
## 4  0   0   0   0 mg N/L Ammonia 2 1674272015 1/20/2023 22:33
## 5  0   0   0   0 mg N/L Ammonia 2 1674272175 1/20/2023 22:36
## 6  0   0   0   0 mg N/L Ammonia 2 1674272335 1/20/2023 22:38

```

```

#take out only the columns that we need
dat3 <- file3[ ,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(dat3) <- c('Run_Info', 'Sample_Name', 'Conc', "Abs", "Units", "Test")
head(dat3)

```

```

##      Run_Info    Sample_Name    Conc     Abs   Units     Test
## 1 RUNSTARTED             NA        NA
## 2 RUNSTARTED             NA        NA
## 3      RESULT Standard 1 0.015424 0.015424 mg N/L Ammonia 2
## 4      RESULT Standard .0389 0.037041 0.037041 mg N/L Ammonia 2
## 5      RESULT Standard .1000 0.065641 0.065641 mg N/L Ammonia 2
## 6      RESULT Standard .2000 0.120285 0.120285 mg N/L Ammonia 2

```

```

#read in data
file4 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NH3_P04_202210_4.csv")
#Quick look at dataframe
head(file4)

```

```

##  RUNSTARTED X1674261032 X1.20.2023.19.30          X X.1      X.2      X.3
## 1 RUNSTARTED 1674261064 1/20/2023 19:31           NA      NA      NA
## 2 RUNSTARTED 1674270611 1/20/2023 22:10           NA      NA      NA
## 3      RESULT      -1           S1 Standard 1 0 0.015424 0.015424
## 4      RESULT      -2           S90 Standard .0389 1 0.037041 0.037041
## 5      RESULT      -2           S91 Standard .1000 2 0.065641 0.065641
## 6      RESULT      -2           S92 Standard .2000 3 0.120285 0.120285
## X.4 X.5 X.6 X.7     X.8       X.9       X.10        X.11
## 1  NA  NA  NA  NA
## 2  NA  NA  NA  NA
## 3  0   0   0   0 mg N/L Ammonia 2 1674271855 1/20/2023 22:30
## 4  0   0   0   0 mg N/L Ammonia 2 1674272015 1/20/2023 22:33
## 5  0   0   0   0 mg N/L Ammonia 2 1674272175 1/20/2023 22:36
## 6  0   0   0   0 mg N/L Ammonia 2 1674272335 1/20/2023 22:38

```

```

#take out only the columns that we need
dat4 <- file4[ ,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(dat4) <- c('Run_Info', 'Sample_Name', 'Conc', "Abs", "Units", "Test")
head(dat4)

```

```

##      Run_Info    Sample_Name    Conc     Abs   Units     Test
## 1 RUNSTARTED             NA        NA

```

```

## 2 RUNSTARTED NA NA
## 3 RESULT Standard 1 0.015424 0.015424 mg N/L Ammonia 2
## 4 RESULT Standard .0389 0.037041 0.037041 mg N/L Ammonia 2
## 5 RESULT Standard .1000 0.065641 0.065641 mg N/L Ammonia 2
## 6 RESULT Standard .2000 0.120285 0.120285 mg N/L Ammonia 2

#read in data
file5 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NH3_P04_202210-11_4.csv")
#Quick look at dataframe
head(file5)

## RUNSTARTED X1674434960 X1.22.2023.19.49 X X.1 X.2
## 1 RUNENDED 1674434979 STOP IMMEDIATE PRESSED NA NA
## 2 RUNSTARTED 1674435072 1/22/2023 19:51 NA NA
## 3 RESULT -1 S1 Standard 1 0 0.015639
## 4 RESULT -2 S90 Standard .0389 1 0.035279
## 5 RESULT -2 S91 Standard .1000 2 0.064059
## 6 RESULT -2 S92 Standard .2000 3 0.117862
## X.3 X.4 X.5 X.6 X.7 X.8 X.9 X.10 X.11
## 1 NA NA NA NA NA NA NA
## 2 NA NA NA NA NA NA NA
## 3 0.015639 0 0 0 512 mg N/L Ammonia 2 1674436305 1/22/2023 20:11
## 4 0.035279 0 0 0 512 mg N/L Ammonia 2 1674436465 1/22/2023 20:14
## 5 0.064059 0 0 0 512 mg N/L Ammonia 2 1674436625 1/22/2023 20:17
## 6 0.117862 0 0 0 512 mg N/L Ammonia 2 1674436785 1/22/2023 20:19

#take out only the columns that we need
dat5 <- file5[,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(dat5) <- c('Run_Info', 'Sample_Name', 'Conc', "Abs", "Units", "Test")
head(dat5)

## Run_Info Sample_Name Conc Abs Units Test
## 1 RUNENDED NA NA
## 2 RUNSTARTED NA NA
## 3 RESULT Standard 1 0.015639 0.015639 mg N/L Ammonia 2
## 4 RESULT Standard .0389 0.035279 0.035279 mg N/L Ammonia 2
## 5 RESULT Standard .1000 0.064059 0.064059 mg N/L Ammonia 2
## 6 RESULT Standard .2000 0.117862 0.117862 mg N/L Ammonia 2

alldat <- rbind(dat1, dat2, dat3, dat4, dat5)

#Pull out standards
stds <- alldat[alldat$Sample_Name %like% "Standard", ]
head(stds)

## Run_Info Sample_Name Conc Abs Units Test
## 1 RESULT Standard 1 0.013208 0.013208 mg N/L Ammonia 2
## 2 RESULT Standard .0389 0.033938 0.033938 mg N/L Ammonia 2
## 3 RESULT Standard .1000 0.064862 0.064862 mg N/L Ammonia 2
## 4 RESULT Standard .2000 0.118781 0.118781 mg N/L Ammonia 2
## 5 RESULT Standard .5000 0.274742 0.274742 mg N/L Ammonia 2
## 6 RESULT Standard 1.0000 0.536226 0.536226 mg N/L Ammonia 2

```

```
#Pull out samples
alldat2 <- alldat[alldat$Sample_Name %like% "MSM_", ]
alldat2 <- rbind(alldat2, (alldat[alldat$Sample_Name %like% "GWI_", ]))
alldat2 <- rbind(alldat2, (alldat[alldat$Sample_Name %like% "GCrew_", ]))
head(alldat)
```

```
##      Run_Info      Sample_Name    Conc     Abs   Units     Test
## 489  RESULT MSM_202210_UP_PPR_1 0.313359 0.178272 mg N/L Ammonia 2
## 490  RESULT MSM_202210_UP_PPR_2 0.343022 0.193470 mg N/L Ammonia 2
## 491  RESULT MSM_202210_UP_PPR_3 0.467758 0.257378 mg N/L Ammonia 2
## 492  RESULT MSM_202210_UP_PPR_4 0.645589 0.348490 mg N/L Ammonia 2
## 493  RESULT MSM_202210_UP_PPR_5 0.794490 0.424779 mg N/L Ammonia 2
## 494  RESULT MSM_202210_UP_PPR_6 0.450865 0.248723 mg N/L Ammonia 2
```

NOx

```
#read in data
Nfile1 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NOx_202210_1.csv")
#Quick look at dataframe
head(Nfile1)
```

```
##  RUNSTARTED X1674058557 X1.18.2023.11.15          X X.1      X.2      X.3 X.4
## 1  RESULT      -1           S1 Standard 1  0 0.138425 0.138425 0
## 2  RESULT      -2           S90 Standard 90  1 0.156085 0.156085 0
## 3  RESULT      -2           S91 Standard 91  2 0.181254 0.181254 0
## 4  RESULT      -2           S92 Standard 92  3 0.224592 0.224592 0
## 5  RESULT      -2           S93 Standard 93  4 0.338805 0.338805 0
## 6  RESULT      -2           S94 Standard 94  5 0.546026 0.546026 0
##  X.5 X.6 X.7   X.8       X.9      X.10      X.11
## 1  0  0  0 mg N/L Vanadium NOx 1674060737 1/18/2023 11:52
## 2  0  0  0 mg N/L Vanadium NOx 1674060816 1/18/2023 11:53
## 3  0  0  0 mg N/L Vanadium NOx 1674060896 1/18/2023 11:54
## 4  0  0  0 mg N/L Vanadium NOx 1674060977 1/18/2023 11:56
## 5  0  0  0 mg N/L Vanadium NOx 1674061057 1/18/2023 11:57
## 6  0  0  0 mg N/L Vanadium NOx 1674061137 1/18/2023 11:58
```

```
#take out only the columns that we need
Ndat1 <- Nfile1[, c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(Ndat1) <- c('Run_Info', 'Sample_Name', 'Conc', "Abs", "Units", "Test")
head(Ndat1)
```

```
##      Run_Info Sample_Name    Conc     Abs   Units     Test
## 1  RESULT Standard 1 0.138425 0.138425 mg N/L Vanadium NOx
## 2  RESULT Standard 90 0.156085 0.156085 mg N/L Vanadium NOx
## 3  RESULT Standard 91 0.181254 0.181254 mg N/L Vanadium NOx
## 4  RESULT Standard 92 0.224592 0.224592 mg N/L Vanadium NOx
## 5  RESULT Standard 93 0.338805 0.338805 mg N/L Vanadium NOx
## 6  RESULT Standard 94 0.546026 0.546026 mg N/L Vanadium NOx
```

```

#read in data
Nfile2 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NOx_202210_2.csv")
#Quick look at dataframe
head(Nfile2)

## RUNSTARTED X1674088994 X1.18.2023.19.43          X X.1      X.2      X.3 X.4
## 1   RESULT      -1           S1 Standard 1  0 0.142252 0.142252 0
## 2   RESULT      -2           S90 Standard 90 1 0.161788 0.161788 0
## 3   RESULT      -2           S91 Standard 91 2 0.187268 0.187268 0
## 4   RESULT      -2           S92 Standard 92 3 0.229701 0.229701 0
## 5   RESULT      -2           S93 Standard 93 4 0.350446 0.350446 0
## 6   RESULT      -2           S94 Standard 94 5 0.552649 0.552649 0
##   X.5 X.6 X.7    X.8          X.9      X.10      X.11
## 1  0  0  0 mg N/L Vanadium NOx 1674091173 1/18/2023 20:19
## 2  0  0  0 mg N/L Vanadium NOx 1674091253 1/18/2023 20:20
## 3  0  0  0 mg N/L Vanadium NOx 1674091334 1/18/2023 20:22
## 4  0  0  0 mg N/L Vanadium NOx 1674091413 1/18/2023 20:23
## 5  0  0  0 mg N/L Vanadium NOx 1674091493 1/18/2023 20:24
## 6  0  0  0 mg N/L Vanadium NOx 1674091573 1/18/2023 20:26

```

```

#take out only the columns that we need
Ndat2 <- Nfile2[ ,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(Ndat2) <- c('Run_Info', 'Sample_Name', 'Conc', "Abs", "Units", "Test")
head(Ndat2)

```

```

## Run_Info Sample_Name     Conc     Abs   Units      Test
## 1   RESULT Standard 1 0.142252 0.142252 mg N/L Vanadium NOx
## 2   RESULT Standard 90 0.161788 0.161788 mg N/L Vanadium NOx
## 3   RESULT Standard 91 0.187268 0.187268 mg N/L Vanadium NOx
## 4   RESULT Standard 92 0.229701 0.229701 mg N/L Vanadium NOx
## 5   RESULT Standard 93 0.350446 0.350446 mg N/L Vanadium NOx
## 6   RESULT Standard 94 0.552649 0.552649 mg N/L Vanadium NOx

```

```

#read in data
Nfile3 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NOx_202210_3.csv")
#Quick look at dataframe
head(Nfile3)

```

```

## RUNSTARTED X1674251527 X1.20.2023.16.52          X X.1      X.2      X.3 X.4
## 1   RESULT      -1           S1 Standard 1  0 0.181183 0.181183 0
## 2   RESULT      -2           S90 Standard 90 1 0.198214 0.198214 0
## 3   RESULT      -2           S91 Standard 91 2 0.223312 0.223312 0
## 4   RESULT      -2           S92 Standard 92 3 0.265110 0.265110 0
## 5   RESULT      -2           S93 Standard 93 4 0.383191 0.383191 0
## 6   RESULT      -2           S94 Standard 94 5 0.579907 0.579907 0
##   X.5 X.6 X.7    X.8          X.9      X.10      X.11
## 1  0  0  0 mg N/L Vanadium NOx 1674253722 1/20/2023 17:28
## 2  0  0  0 mg N/L Vanadium NOx 1674253804 1/20/2023 17:30
## 3  0  0  0 mg N/L Vanadium NOx 1674253884 1/20/2023 17:31
## 4  0  0  0 mg N/L Vanadium NOx 1674253966 1/20/2023 17:32
## 5  0  0  0 mg N/L Vanadium NOx 1674254046 1/20/2023 17:34
## 6  0  0  0 mg N/L Vanadium NOx 1674254127 1/20/2023 17:35

```

```

#take out only the columns that we need
Ndat3 <- Nfile3[ ,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(Ndat3) <- c('Run_Info','Sample_Name','Conc', "Abs", "Units", "Test")
head(Ndat3)

```

```

##   Run_Info Sample_Name      Conc      Abs  Units      Test
## 1  RESULT  Standard 1 0.181183 0.181183 mg N/L Vanadium NOx
## 2  RESULT  Standard 90 0.198214 0.198214 mg N/L Vanadium NOx
## 3  RESULT  Standard 91 0.223312 0.223312 mg N/L Vanadium NOx
## 4  RESULT  Standard 92 0.265110 0.265110 mg N/L Vanadium NOx
## 5  RESULT  Standard 93 0.383191 0.383191 mg N/L Vanadium NOx
## 6  RESULT  Standard 94 0.579907 0.579907 mg N/L Vanadium NOx

```

```

#read in data
Nfile4 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NOx_202210_4.csv")
#Quick look at dataframe
head(Nfile4)

```

```

##   RUNSTARTED X1674330201 X1.21.2023.14.43          X X.1      X.2      X.3 X.4
## 1  RESULT      -1           S1 Standard 1  0 0.182628 0.182628 0
## 2  RESULT      -2           S90 Standard 90  1 0.201200 0.201200 0
## 3  RESULT      -2           S91 Standard 91  2 0.224336 0.224336 0
## 4  RESULT      -2           S92 Standard 92  3 0.264526 0.264526 0
## 5  RESULT      -2           S93 Standard 93  4 0.384563 0.384563 0
## 6  RESULT      -2           S94 Standard 94  5 0.587757 0.587757 0
##   X.5 X.6 X.7     X.8      X.9      X.10      X.11
## 1  0  0  0 mg N/L Vanadium NOx 1674332394 1/21/2023 15:19
## 2  0  0  0 mg N/L Vanadium NOx 1674332475 1/21/2023 15:21
## 3  0  0  0 mg N/L Vanadium NOx 1674332556 1/21/2023 15:22
## 4  0  0  0 mg N/L Vanadium NOx 1674332637 1/21/2023 15:23
## 5  0  0  0 mg N/L Vanadium NOx 1674332718 1/21/2023 15:25
## 6  0  0  0 mg N/L Vanadium NOx 1674332799 1/21/2023 15:26

```

```

#take out only the columns that we need
Ndat4 <- Nfile4[ ,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(Ndat4) <- c('Run_Info','Sample_Name','Conc', "Abs", "Units", "Test")
head(Ndat4)

```

```

##   Run_Info Sample_Name      Conc      Abs  Units      Test
## 1  RESULT  Standard 1 0.182628 0.182628 mg N/L Vanadium NOx
## 2  RESULT  Standard 90 0.201200 0.201200 mg N/L Vanadium NOx
## 3  RESULT  Standard 91 0.224336 0.224336 mg N/L Vanadium NOx
## 4  RESULT  Standard 92 0.264526 0.264526 mg N/L Vanadium NOx
## 5  RESULT  Standard 93 0.384563 0.384563 mg N/L Vanadium NOx
## 6  RESULT  Standard 94 0.587757 0.587757 mg N/L Vanadium NOx

```

```

#read in data
Nfile5 <- read.csv("Raw Data/SEAL_COMPASS_Synoptic_NOx_202210-11_4.csv")
#Quick look at dataframe
head(Nfile5)

```

```

## RUNSTARTED X1674567913      X1.24.2023.8.45          X.X.1      X.2      X.3
## 1  RUNENDED  1674567954 INSUFFICIENT REAGENT          NA       NA       NA
## 2 RUNSTARTED  1674569191      1/24/2023 9:06          NA       NA       NA
## 3   RESULT      -1           S1 Standard 1  0 0.178421 0.178421
## 4   RESULT      -2           S90 Standard 90  1 0.198991 0.198991
## 5   RESULT      -2           S91 Standard 91  2 0.218318 0.218318
## 6   RESULT      -2           S92 Standard 92  3 0.260783 0.260783
##   X.4 X.5 X.6 X.7     X.8      X.9      X.10      X.11
## 1  NA  NA  NA  NA          NA
## 2  NA  NA  NA  NA          NA
## 3  0  0  0 512 mg N/L Vanadium NOx 1674571379 1/24/2023 9:42
## 4  0  0  0 512 mg N/L Vanadium NOx 1674571460 1/24/2023 9:44
## 5  0  0  0 512 mg N/L Vanadium NOx 1674571540 1/24/2023 9:45
## 6  0  0  0 512 mg N/L Vanadium NOx 1674571621 1/24/2023 9:47

```

```

#take out only the columns that we need
Ndat5 <- Nfile5[ ,c(1,4,6,7, 12, 13)]
# assigning new names to the columns of the data frame
colnames(Ndat5) <- c('Run_Info', 'Sample_Name', 'Conc', "Abs", "Units", "Test")
head(Ndat5)

```

```

##   Run_Info Sample_Name      Conc      Abs    Units      Test
## 1  RUNENDED          NA       NA
## 2 RUNSTARTED          NA       NA
## 3   RESULT Standard 1 0.178421 0.178421 mg N/L Vanadium NOx
## 4   RESULT Standard 90 0.198991 0.198991 mg N/L Vanadium NOx
## 5   RESULT Standard 91 0.218318 0.218318 mg N/L Vanadium NOx
## 6   RESULT Standard 92 0.260783 0.260783 mg N/L Vanadium NOx

```

```
Nalldat <- rbind(Ndat1, Ndat2, Ndat3, Ndat4, Ndat5)
```

```

#Pull out standards
Nstds <- Nalldat[Nalldat$Sample_Name %like% "Standard", ]
head(Nstds)

```

```

##   Run_Info Sample_Name      Conc      Abs    Units      Test
## 1   RESULT Standard 1 0.138425 0.138425 mg N/L Vanadium NOx
## 2   RESULT Standard 90 0.156085 0.156085 mg N/L Vanadium NOx
## 3   RESULT Standard 91 0.181254 0.181254 mg N/L Vanadium NOx
## 4   RESULT Standard 92 0.224592 0.224592 mg N/L Vanadium NOx
## 5   RESULT Standard 93 0.338805 0.338805 mg N/L Vanadium NOx
## 6   RESULT Standard 94 0.546026 0.546026 mg N/L Vanadium NOx

```

```

#Pull out samples
Nalldat2 <- Nalldat[Nalldat$Sample_Name %like% "MSM_", ]
Nalldat2 <- rbind(Nalldat2, (Nalldat[Nalldat$Sample_Name %like% "GWI_", ]))
Nalldat2 <- rbind(Nalldat2, (Nalldat[Nalldat$Sample_Name %like% "GCrew_", ]))
head(Nalldat2)

```

```

##   Run_Info          Sample_Name      Conc      Abs    Units      Test
## 231   RESULT MSM_202210_UP_LysA_20cm 0.000456 0.182796 mg N/L Vanadium NOx
## 232   RESULT MSM_202210_UP_LysA_45cm 0.000638 0.182941 mg N/L Vanadium NOx

```

```

## 233  RESULT MSM_202210_UP_LysB_10cm -0.002491 0.180455 mg N/L Vanadium NOx
## 234  RESULT MSM_202210_UP_LysB_20cm  0.001671 0.183762 mg N/L Vanadium NOx
## 235  RESULT MSM_202210_UP_LysB_45cm -0.001690 0.181091 mg N/L Vanadium NOx
## 236  RESULT MSM_202210_UP_LysC_10cm -0.005525 0.178046 mg N/L Vanadium NOx

```

Constants

```

N_mw <- 14.0067    # molecular weight of N

P_mw <- 30.973762  # molecular weight of P

Con1 <- 1000        # conversion factor value

Con2 <- 1000000     # conversion factor value

```

Convert Data from mg/L to uM

```
head(alldat2)
```

```

##      Run_Info      Sample_Name      Conc      Abs      Units      Test
## 489  RESULT MSM_202210_UP_PPR_1 0.313359 0.178272 mg N/L Ammonia 2
## 490  RESULT MSM_202210_UP_PPR_2 0.343022 0.193470 mg N/L Ammonia 2
## 491  RESULT MSM_202210_UP_PPR_3 0.467758 0.257378 mg N/L Ammonia 2
## 492  RESULT MSM_202210_UP_PPR_4 0.645589 0.348490 mg N/L Ammonia 2
## 493  RESULT MSM_202210_UP_PPR_5 0.794490 0.424779 mg N/L Ammonia 2
## 494  RESULT MSM_202210_UP_PPR_6 0.450865 0.248723 mg N/L Ammonia 2

```

```
head(Nalldat2)
```

```

##      Run_Info      Sample_Name      Conc      Abs      Units      Test
## 231  RESULT MSM_202210_UP_LysA_20cm 0.000456 0.182796 mg N/L Vanadium NOx
## 232  RESULT MSM_202210_UP_LysA_45cm 0.000638 0.182941 mg N/L Vanadium NOx
## 233  RESULT MSM_202210_UP_LysB_10cm -0.002491 0.180455 mg N/L Vanadium NOx
## 234  RESULT MSM_202210_UP_LysB_20cm  0.001671 0.183762 mg N/L Vanadium NOx
## 235  RESULT MSM_202210_UP_LysB_45cm -0.001690 0.181091 mg N/L Vanadium NOx
## 236  RESULT MSM_202210_UP_LysC_10cm -0.005525 0.178046 mg N/L Vanadium NOx

```

```
#subset by test
NH4samples <- subset(alldat2, Test == "Ammonia 2")
head(NH4samples)
```

```

##      Run_Info      Sample_Name      Conc      Abs      Units      Test
## 489  RESULT MSM_202210_UP_PPR_1 0.313359 0.178272 mg N/L Ammonia 2
## 490  RESULT MSM_202210_UP_PPR_2 0.343022 0.193470 mg N/L Ammonia 2
## 491  RESULT MSM_202210_UP_PPR_3 0.467758 0.257378 mg N/L Ammonia 2
## 492  RESULT MSM_202210_UP_PPR_4 0.645589 0.348490 mg N/L Ammonia 2
## 493  RESULT MSM_202210_UP_PPR_5 0.794490 0.424779 mg N/L Ammonia 2
## 494  RESULT MSM_202210_UP_PPR_6 0.450865 0.248723 mg N/L Ammonia 2

```

```
P04samples <- subset(alldat2, Test == "o-PHOS 0.3")
head(P04samples)
```

##	Run_Info	Sample_Name	Conc	Abs	Units	Test
##	631	RESULT MSM_202210_UP_PPR_1	0.003792	0.003058	mg P/L	o-PHOS 0.3
##	632	RESULT MSM_202210_UP_PPR_2	0.003947	0.003122	mg P/L	o-PHOS 0.3
##	633	RESULT MSM_202210_UP_PPR_3	0.007060	0.004405	mg P/L	o-PHOS 0.3
##	634	RESULT MSM_202210_UP_PPR_4	0.009802	0.005535	mg P/L	o-PHOS 0.3
##	635	RESULT MSM_202210_UP_PPR_5	0.009451	0.005390	mg P/L	o-PHOS 0.3
##	636	RESULT MSM_202210_UP_PPR_6	0.005077	0.003588	mg P/L	o-PHOS 0.3

```
NOXsamples <- subset(Nalldat2, Test == "Vanadium NOx")
head(NOXsamples)
```

##	Run_Info	Sample_Name	Conc	Abs	Units	Test
##	231	RESULT MSM_202210_UP_LysA_20cm	0.000456	0.182796	mg N/L	Vanadium NOx
##	232	RESULT MSM_202210_UP_LysA_45cm	0.000638	0.182941	mg N/L	Vanadium NOx
##	233	RESULT MSM_202210_UP_LysB_10cm	-0.002491	0.180455	mg N/L	Vanadium NOx
##	234	RESULT MSM_202210_UP_LysB_20cm	0.001671	0.183762	mg N/L	Vanadium NOx
##	235	RESULT MSM_202210_UP_LysB_45cm	-0.001690	0.181091	mg N/L	Vanadium NOx
##	236	RESULT MSM_202210_UP_LysC_10cm	-0.005525	0.178046	mg N/L	Vanadium NOx

```
NH4samples$Conc_uM <- (((as.numeric(NH4samples$Conc))/Con1)/N_mw)*Con2
head(NH4samples)
```

##	Run_Info	Sample_Name	Conc	Abs	Units	Test	Conc_uM
##	489	RESULT MSM_202210_UP_PPR_1	0.313359	0.178272	mg N/L	Ammonia 2	22.37208
##	490	RESULT MSM_202210_UP_PPR_2	0.343022	0.193470	mg N/L	Ammonia 2	24.48985
##	491	RESULT MSM_202210_UP_PPR_3	0.467758	0.257378	mg N/L	Ammonia 2	33.39530
##	492	RESULT MSM_202210_UP_PPR_4	0.645589	0.348490	mg N/L	Ammonia 2	46.09144
##	493	RESULT MSM_202210_UP_PPR_5	0.794490	0.424779	mg N/L	Ammonia 2	56.72214
##	494	RESULT MSM_202210_UP_PPR_6	0.450865	0.248723	mg N/L	Ammonia 2	32.18924

```
P04samples$Conc_uM <- (((as.numeric(P04samples$Conc))/Con1)/N_mw)*Con2
head(P04samples)
```

##	Run_Info	Sample_Name	Conc	Abs	Units	Test	Conc_uM
##	631	RESULT MSM_202210_UP_PPR_1	0.003792	0.003058	mg P/L	o-PHOS 0.3	0.2707276
##	632	RESULT MSM_202210_UP_PPR_2	0.003947	0.003122	mg P/L	o-PHOS 0.3	0.2817937
##	633	RESULT MSM_202210_UP_PPR_3	0.007060	0.004405	mg P/L	o-PHOS 0.3	0.5040445
##	634	RESULT MSM_202210_UP_PPR_4	0.009802	0.005535	mg P/L	o-PHOS 0.3	0.6998079
##	635	RESULT MSM_202210_UP_PPR_5	0.009451	0.005390	mg P/L	o-PHOS 0.3	0.6747485
##	636	RESULT MSM_202210_UP_PPR_6	0.005077	0.003588	mg P/L	o-PHOS 0.3	0.3624694

```
NOXsamples$Conc_uM_raw <- (((as.numeric(NOXsamples$Conc))/Con1)/N_mw)*Con2
head(NOXsamples)
```

##	Run_Info	Sample_Name	Conc	Abs	Units	Test
##	231	RESULT MSM_202210_UP_LysA_20cm	0.000456	0.182796	mg N/L	Vanadium NOx
##	232	RESULT MSM_202210_UP_LysA_45cm	0.000638	0.182941	mg N/L	Vanadium NOx

```

## 233 RESULT MSM_202210_UP_LysB_10cm -0.002491 0.180455 mg N/L Vanadium NOx
## 234 RESULT MSM_202210_UP_LysB_20cm 0.001671 0.183762 mg N/L Vanadium NOx
## 235 RESULT MSM_202210_UP_LysB_45cm -0.001690 0.181091 mg N/L Vanadium NOx
## 236 RESULT MSM_202210_UP_LysC_10cm -0.005525 0.178046 mg N/L Vanadium NOx
## Conc_uM_raw
## 231 0.03255585
## 232 0.04554963
## 233 -0.17784346
## 234 0.11930005
## 235 -0.12065654
## 236 -0.39445408

#add step to make negative values equal to bd (below detection) and replace with zeros
NOXsamples$Conc_uM <- ifelse(NOXsamples$Conc_uM_raw<0, 0, (NOXsamples$Conc_uM_raw) )
head(NOXsamples)

```

```

## Run_Info Sample_Name Conc Abs Units Test
## 231 RESULT MSM_202210_UP_LysA_20cm 0.000456 0.182796 mg N/L Vanadium NOx
## 232 RESULT MSM_202210_UP_LysA_45cm 0.000638 0.182941 mg N/L Vanadium NOx
## 233 RESULT MSM_202210_UP_LysB_10cm -0.002491 0.180455 mg N/L Vanadium NOx
## 234 RESULT MSM_202210_UP_LysB_20cm 0.001671 0.183762 mg N/L Vanadium NOx
## 235 RESULT MSM_202210_UP_LysB_45cm -0.001690 0.181091 mg N/L Vanadium NOx
## 236 RESULT MSM_202210_UP_LysC_10cm -0.005525 0.178046 mg N/L Vanadium NOx
## Conc_uM_raw Conc_uM
## 231 0.03255585 0.03255585
## 232 0.04554963 0.04554963
## 233 -0.17784346 0.00000000
## 234 0.11930005 0.11930005
## 235 -0.12065654 0.00000000
## 236 -0.39445408 0.00000000

```

Pull all data back together and add flags

```

#pull out the columns we want from each dataframe
NH4_pull <- NH4samples[ ,c(2,3,7) ]
head(NH4_pull)

```

```

## Sample_Name Conc Conc_uM
## 489 MSM_202210_UP_PPR_1 0.313359 22.37208
## 490 MSM_202210_UP_PPR_2 0.343022 24.48985
## 491 MSM_202210_UP_PPR_3 0.467758 33.39530
## 492 MSM_202210_UP_PPR_4 0.645589 46.09144
## 493 MSM_202210_UP_PPR_5 0.794490 56.72214
## 494 MSM_202210_UP_PPR_6 0.450865 32.18924

```

```

P04_pull <- P04samples[ ,c(2,3,7) ]
head(P04_pull)

```

```

## Sample_Name Conc Conc_uM
## 631 MSM_202210_UP_PPR_1 0.003792 0.2707276

```

```

## 632 MSM_202210_UP_PPR_2 0.003947 0.2817937
## 633 MSM_202210_UP_PPR_3 0.007060 0.5040445
## 634 MSM_202210_UP_PPR_4 0.009802 0.6998079
## 635 MSM_202210_UP_PPR_5 0.009451 0.6747485
## 636 MSM_202210_UP_PPR_6 0.005077 0.3624694

```

```

NOX_pull <- NOXsamples[ ,c(2,3,8) ]
head(NOX_pull)

```

```

##             Sample_Name      Conc      Conc_uM
## 231 MSM_202210_UP_LysA_20cm 0.000456 0.03255585
## 232 MSM_202210_UP_LysA_45cm 0.000638 0.04554963
## 233 MSM_202210_UP_LysB_10cm -0.002491 0.00000000
## 234 MSM_202210_UP_LysB_20cm 0.001671 0.11930005
## 235 MSM_202210_UP_LysB_45cm -0.001690 0.00000000
## 236 MSM_202210_UP_LysC_10cm -0.005525 0.00000000

```

```

#Bring all this data back together:
all_data <- merge(NH4_pull, PO4_pull, by="Sample_Name", all.x=TRUE)
all_data <- merge(all_data, NOX_pull, by="Sample_Name", all.x=TRUE)
head(all_data)

```

```

##             Sample_Name  Conc.x Conc_uM.x  Conc.y Conc_uM.y      Conc      Conc_uM
## 1 GCrew_202210_SW_A 0.249386 17.804765 0.125976 8.9939815 0.007561 0.5398131
## 2 GCrew_202210_SW_B 0.305534 21.813418 0.137691 9.8303669 0.008117 0.5795084
## 3 GCrew_202210_SW_C 0.124176  8.865472 0.002949 0.2105421 -0.069967 0.0000000
## 4 GCrew_202210_SW_C 0.124176  8.865472 0.002949 0.2105421  0.004603 0.3286284
## 5 GCrew_202210_SW_C 0.124176  8.865472 0.069649 4.9725489 -0.069967 0.0000000
## 6 GCrew_202210_SW_C 0.124176  8.865472 0.069649 4.9725489  0.004603 0.3286284

```

```

colnames(all_data) <- c("Sample_Name", "NH3_mgL", "NH3_uM", "PO4_mgL", "PO4_uM", "NOx_mgL", "NOx_uM")
head(all_data)

```

```

##             Sample_Name  NH3_mgL    NH3_uM  PO4_mgL     PO4_uM  NOx_mgL    NOx_uM
## 1 GCrew_202210_SW_A 0.249386 17.804765 0.125976 8.9939815 0.007561 0.5398131
## 2 GCrew_202210_SW_B 0.305534 21.813418 0.137691 9.8303669 0.008117 0.5795084
## 3 GCrew_202210_SW_C 0.124176  8.865472 0.002949 0.2105421 -0.069967 0.0000000
## 4 GCrew_202210_SW_C 0.124176  8.865472 0.002949 0.2105421  0.004603 0.3286284
## 5 GCrew_202210_SW_C 0.124176  8.865472 0.069649 4.9725489 -0.069967 0.0000000
## 6 GCrew_202210_SW_C 0.124176  8.865472 0.069649 4.9725489  0.004603 0.3286284

```

```

#add in an if then statement that tells us if they are within the range of the test - check this after
all_data$NH3_range <- ifelse(all_data$NH3_mgL<0.02, "bdl", ifelse(all_data$NH3_mgL>2, "adl", "Within_R"))
all_data$PO4_range <- ifelse(all_data$PO4_mgL<0.003, "bdl", ifelse(all_data$PO4_mgL>3, "adl", "Within_R"))
all_data$NOx_range <- ifelse(all_data$NOx_mgL<0.025, "bdl", ifelse(all_data$NOx_mgL>1, "adl", "Within_R"))
head(all_data)

```

```

##             Sample_Name  NH3_mgL    NH3_uM  PO4_mgL     PO4_uM  NOx_mgL    NOx_uM
## 1 GCrew_202210_SW_A 0.249386 17.804765 0.125976 8.9939815 0.007561 0.5398131
## 2 GCrew_202210_SW_B 0.305534 21.813418 0.137691 9.8303669 0.008117 0.5795084
## 3 GCrew_202210_SW_C 0.124176  8.865472 0.002949 0.2105421 -0.069967 0.0000000

```

```

## 4 GCrew_202210_SW_C 0.124176 8.865472 0.002949 0.2105421 0.004603 0.3286284
## 5 GCrew_202210_SW_C 0.124176 8.865472 0.069649 4.9725489 -0.069967 0.0000000
## 6 GCrew_202210_SW_C 0.124176 8.865472 0.069649 4.9725489 0.004603 0.3286284
##      NH3_range    PO4_range NOx_range
## 1 Within_Range Within_Range      bdl
## 2 Within_Range Within_Range      bdl
## 3 Within_Range          bdl      bdl
## 4 Within_Range          bdl      bdl
## 5 Within_Range Within_Range      bdl
## 6 Within_Range Within_Range      bdl

```

Take an initial look at concentrations

```


#plot data to get a sense of any outliers
NH3look <- ggplot(data=all_data, aes(x=Sample_Name, y=NH3_uM)) +
  geom_bar(stat="identity") +
  theme_classic() + ylim(-10, 200) +
  theme(legend.position="none") +
  ggtitle("Sample NH3 Concentrations")
NH3look

```

```


## Warning: Removed 28 rows containing missing values or values outside the scale range
## ('geom_bar()').

```



```
P04look <- ggplot(data=all_data, aes(x=Sample_Name, y=P04_uM)) +
  geom_bar(stat="identity") +
  theme_classic() + ylim(-2, 50) +
  theme(legend.position="none") +
  ggtitle("Sample P04 Concentrations")
P04look
```


```
## Warning: Removed 5 rows containing missing values or values outside the scale range
## ('geom_bar()'').
```



```
NOXlook <- ggplot(data=all_data, aes(x=Sample_Name, y=NOx_uM)) +
  geom_bar(stat="identity") +
  theme_classic() + ylim(-2, 2) +
  theme(legend.position="none") +
  ggtitle("Sample NOx Concentrations")
NOXlook
```

```
## Warning: Removed 80 rows containing missing values or values outside the scale range
## ('geom_bar()'').
```

Sample NOx Concentrations



Pull out data you need, make IDs

```
head(all_data)
```

```
##           Sample_Name  NH3_mgL    NH3_uM  PO4_mgL    PO4_uM  NOx_mgL    NOx_uM
## 1 GCrew_202210_SW_A 0.249386 17.804765 0.125976 8.9939815 0.007561 0.5398131
## 2 GCrew_202210_SW_B 0.305534 21.813418 0.137691 9.8303669 0.008117 0.5795084
## 3 GCrew_202210_SW_C 0.124176  8.865472 0.002949 0.2105421 -0.069967 0.0000000
## 4 GCrew_202210_SW_C 0.124176  8.865472 0.002949 0.2105421  0.004603 0.3286284
## 5 GCrew_202210_SW_C 0.124176  8.865472 0.069649 4.9725489 -0.069967 0.0000000
## 6 GCrew_202210_SW_C 0.124176  8.865472 0.069649 4.9725489  0.004603 0.3286284
##      NH3_range    PO4_range NOx_range
## 1 Within_Range Within_Range      bdl
## 2 Within_Range Within_Range      bdl
## 3 Within_Range          bdl      bdl
## 4 Within_Range          bdl      bdl
## 5 Within_Range Within_Range      bdl
## 6 Within_Range Within_Range      bdl
```

```
out <- all_data[, c(1,3,5,7,8,9,10)]
head(out)
```

	Sample_Name	NH3_uM	PO4_uM	NOx_uM	NH3_range	PO4_range
1	GCrew_202210_SW_A	0.249386	17.804765	0.007561	Within_Range	Within_Range
2	GCrew_202210_SW_B	0.305534	21.813418	0.008117	Within_Range	Within_Range
3	GCrew_202210_SW_C	0.124176	8.865472	-0.069967	Within_Range	Within_Range
4	GCrew_202210_SW_C	0.124176	8.865472	0.004603	Within_Range	Within_Range
5	GCrew_202210_SW_C	0.124176	8.865472	-0.069967	Within_Range	Within_Range
6	GCrew_202210_SW_C	0.124176	8.865472	0.004603	Within_Range	Within_Range

```

## 1 GCrew_202210_SW_A 17.804765 8.9939815 0.5398131 Within_Range Within_Range
## 2 GCrew_202210_SW_B 21.813418 9.8303669 0.5795084 Within_Range Within_Range
## 3 GCrew_202210_SW_C 8.865472 0.2105421 0.0000000 Within_Range bdl
## 4 GCrew_202210_SW_C 8.865472 0.2105421 0.3286284 Within_Range bdl
## 5 GCrew_202210_SW_C 8.865472 4.9725489 0.0000000 Within_Range Within_Range
## 6 GCrew_202210_SW_C 8.865472 4.9725489 0.3286284 Within_Range Within_Range
##   NOx_range
## 1      bdl
## 2      bdl
## 3      bdl
## 4      bdl
## 5      bdl
## 6      bdl

#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(out$Sample_Name), '_', fixed=TRUE)))

## Warning in rbind(c("GCrew", "202210", "SW", "A"), c("GCrew", "202210", "SW", :
## number of columns of result is not a multiple of vector length (arg 1)

colnames(IDs) <- c("Site" , "Date", "Zone", "Replicate", "Depth", "Tree")
IDs$Month <- ifelse(IDs>Date == 202210, "October", "November")
#IDs$Month <- "September"
head(IDs)

##     Site Date Zone Replicate Depth Tree Month
## 1 GCrew 202210 SW      A GCrew 202210 October
## 2 GCrew 202210 SW      B GCrew 202210 October
## 3 GCrew 202210 SW      C GCrew 202210 October
## 4 GCrew 202210 SW      C GCrew 202210 October
## 5 GCrew 202210 SW      C GCrew 202210 October
## 6 GCrew 202210 SW      C GCrew 202210 October

#rejoin them to the dataframe
alldat <- cbind(IDs, out)
head(alldat)

##     Site Date Zone Replicate Depth Tree Month Sample_Name NH3_uM
## 1 GCrew 202210 SW      A GCrew 202210 October GCrew_202210_SW_A 17.804765
## 2 GCrew 202210 SW      B GCrew 202210 October GCrew_202210_SW_B 21.813418
## 3 GCrew 202210 SW      C GCrew 202210 October GCrew_202210_SW_C 8.865472
## 4 GCrew 202210 SW      C GCrew 202210 October GCrew_202210_SW_C 8.865472
## 5 GCrew 202210 SW      C GCrew 202210 October GCrew_202210_SW_C 8.865472
## 6 GCrew 202210 SW      C GCrew 202210 October GCrew_202210_SW_C 8.865472
##     P04_uM NOx_uM NH3_range P04_range NOx_range
## 1 8.9939815 0.5398131 Within_Range Within_Range bdl
## 2 9.8303669 0.5795084 Within_Range Within_Range bdl
## 3 0.2105421 0.0000000 Within_Range bdl bdl
## 4 0.2105421 0.3286284 Within_Range bdl bdl
## 5 4.9725489 0.0000000 Within_Range Within_Range bdl
## 6 4.9725489 0.3286284 Within_Range Within_Range bdl

```

Export final data with flags

```
#Export Data  
  
write.csv(alldat, file="Processed Data/COMPASS_Synoptic_CB_SEAL_NUTR_202210-11.csv")
```