

# COMPASS\_SynopticCB\_PW\_SO4\_Cl\_202304\_Template

Stephanie J. Wilson

2025-08-20

## Contents

0.1	Import Sample Data . . . . .	3
0.2	Assessing Standard Curves . . . . .	3
0.3	Assess Check Standards . . . . .	6
0.4	Assess Blanks . . . . .	7
0.5	Assess Duplicates . . . . .	8
0.6	Calculate mmol/L concentrations & salinity, add dilutions . . . . .	9
0.7	Assess Analytical Spikes . . . . .	10
0.8	Sample Flagging - Are samples Within the range of the curve? . . . . .	10
0.9	Visualize Data by Plot . . . . .	11
0.10	Check to see if samples run match metadata & merge info . . . . .	11
0.11	Export Processed Data . . . . .	11

```
##Setup - Change things here & write any notes
```

```
#identify section  
cat("Setup Information")
```

```
## Setup Information
```

```
##### Run information - PLEASE CHANGE  
Date_Run = "09/05/23" #Date that instrument was run  
Run_by = "Stephanie J. Wilson" #Instrument user  
Script_run_by = "Stephanie J. Wilson" #Code user  
run_notes = "Samples missing from metadata: MSM_202311_UP_LYSC_45CM MSM_202311_WC_LYSB_45CM MSM_202311_WC_RHZ_LYSC, 118_MSM_202311_TR_RHZ_SF_Tree_1 value above cal curve for S04, but only slightly below" #Notes on samples  
  
samples <- c("GCW", "GWI", "MSM", "SWH") #whatever identifies your samples within the same names  
samples_pattern <- paste(samples, collapse = "|")  
  #samples_pattern <- "GCW" #use this instead of the line above if you have only one site code  
chks_name = "Check Standard" #what did you name your check standards?  
  
##### File Names - PLEASE CHANGE  
#file path and name for raw summary data file  
raw_file_name_cl = "Raw Data/COMPASS_Synoptic_CB_MonMon_202311_Cl.txt" #example  
raw_file_name_so4 = "Raw Data/COMPASS_Synoptic_CB_MonMon_202311_S04.txt"  
  
#file path and name of processed data file  
processed_file_name = "Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_S04_202311.csv" #example  
  
##### Log Files - PLEASE CHECK  
#downloaded metadata csv - downloaded from Google drive as csv for this year  
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2023.csv"  
  
#qaqc log file path for this year  
Log_path = "Raw Data/COMPASS_Synoptic_Cl_S04_QAQClog_2023.csv"
```

```
##Required Packages
```

```
##Set Up Code
```

```
#Link to the protocol used for analysis  
  #steph will add this soon  
  
#any coefficients / constants that are needed for calculations  
cl_mw <- 35.45 #molecular weight of Chloride: 35.45  
s_mw <- 32.06 #molecular weight of sulfur: 32.06  
Con1 <- 1000 # conversion factor value  
Con2 <- 1000000 # conversion factor value  
  
#Flag for Dionex ##Need to edit this!  
r2_cutoff = 0.98 #this is the level below which we want to rerun or consider a curve  
chk_flag = 5 #for the cv  
chk_conc_flag = 5 #this is the level cutoff for percent difference of check standards vs. the concentration  
chks_flag = 0.50 #this is the percent of chks we want to have a CV less than 10, usually 60  
rep_flag = 25 #this is a 25% error between samples
```

```

#blank_flag - calculated based on samples later in this code as lower 25% quantile of sample concentr

#standard concentrations - Update if running different checks:
standards <- tibble(
  sample_ID = c("Standard 1", "Standard 2", "Standard 3", "Standard 4", "Standard 5"),
  SO4_std_conc = c(0.5, 1.0, 2.0, 10, 20),
  Cl_std_conc = c(5, 10, 20, 100, 200)
)

#Spike concentration calc:
spk_std <- (250/s_mw) # in mM
spkvol <- 10 # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
spk_Conc <- (spk_std)*spkvol

#Top standard Concentrations- Update if running different standard curve:
top_std_cl = 200
top_std_so4 = 20

#Set time zone
common_tz = "Etc/GMT+5"
Sys.setenv(TZ = "America/New_York")

#plot indicators
site_order <- c('GCW', 'MSM', 'GWI', 'SWH')
plot_order <- c('UP', 'SWAMP', 'TR', 'WC', 'SW')
plot_colors <- c("#20063B", "darkgrey", "#FFBC42", "#419973", "#25ABE6" )

```

```
##Read in metadata and create similar sample IDs for matching to samples
```


## 0.1 Import Sample Data

## 0.2 Assessing Standard Curves

```
## Assess the Standard Curves
```


```
## 'geom_smooth()' using formula = 'y ~ x'
```

Chloride Std Curve by Date




```
## 'geom_smooth()' using formula = 'y ~ x'
```

Sulfate Std Curve by Date



```
## 'geom_line()': Each group consists of only one observation.  
## i Do you need to adjust the group aesthetic?
```



```
## [1] "Cl Curve r2 GOOD"
```

```
## [1] "SO4 Curve r2 GOOD"
```

### 0.3 Assess Check Standards


```
## Assess the Check Standards
```

```
## # A tibble: 5 x 5
```

```
##   sample_ID mean_Cl sd_Cl cv_Cl flag_Cl
##   <chr>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1    5.20 0.0382 0.00735 Chloride Check Standard RSD within Range
## 2 Standard 2   10.5 0.0743 0.00711 Chloride Check Standard RSD within Range
## 3 Standard 3   20.7 0.200 0.00965 Chloride Check Standard RSD within Range
## 4 Standard 4  103. 0.819 0.00796 Chloride Check Standard RSD within Range
## 5 Standard 5  204. 1.24 0.00606 Chloride Check Standard RSD within Range
```

```
## # A tibble: 5 x 5
```

```
##   sample_ID mean_S04 sd_S04 cv_S04 flag_S04
##   <chr>      <dbl> <dbl> <dbl> <chr>
## 1 Standard 1    0.504 0.0111 0.0220 Sulfate Check Standard RSD within Range
## 2 Standard 2    1.03 0.0144 0.0140 Sulfate Check Standard RSD within Range
## 3 Standard 3    2.05 0.0215 0.0105 Sulfate Check Standard RSD within Range
## 4 Standard 4   10.3 0.0655 0.00639 Sulfate Check Standard RSD within Range
## 5 Standard 5   20.4 0.132 0.00646 Sulfate Check Standard RSD within Range
```



```
## [1] ">60% of Chloride Check Standards are within range of expected concentration"
```


```
## [1] ">60% of Sulfate Check Standards are within range of expected concentration"
```

#### 0.4 Assess Blanks

```
## Assess Blanks
```

```
## [1] ">60% of Carbon Blank concentrations are lower 25% quartile of samples"
```

```
## [1] ">60% of Nitrogen Blank concentrations are lower 25% quartile of samples"
```



Blank Conc <25% Quartile Samples ■ YI

Blank Conc <25% Quartile Samples ■ Y

```
## Chloride blanks:
```

```
## [1] 0.01699524
```

```
## Sulfate blanks:
```

```
## [1] 0.004957143
```

## 0.5 Assess Duplicates

```
## Assess Duplicates
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
```


```
## i Please use 'linewidth' instead.
```

```
## This warning is displayed once every 8 hours.
```

```
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
```

```
## generated.
```





```
## [1] ">60% of Chloride Duplicates have a CV <10%"
```

```
## [1] ">60% of Sulfate Duplicates have a CV <10%"
```

## 0.6 Calculate mmol/L concentrations & salinity, add dilutions

```
cat("Unit Conversion and Salinity Calculation")

# Convert ppm to mmol/L
all_dat$S04_mM <- (all_dat$S04_ppm / s_mw)
all_dat$Cl_mM <- (all_dat$Cl_ppm / cl_mw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
all_dat$salinity <- ((1.8070 * all_dat$Cl_ppm) + 0.026) / 1000

#Need to determine dilution factors for your samples
#for Steph / COMPASS this depends on the site so...
all_dat$Dilution <- 1
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "UP"),
```

```

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "SW"),

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "UP"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "SW"),

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "UP"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "SW"),

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "SWH"), 50, all_dat$Dilution)


head(all_dat)

```

## 0.7 Assess Analytical Spikes

## 0.8 Sample Flagging - Are samples Within the range of the curve?

## Sample Flagging



**0.9 Visualize Data by Plot**

**0.10 Check to see if samples run match metadata & merge info**

**0.11 Export Processed Data**

*#end*