

Synoptic CB: Porewater SO₄/Cl

July 2025 Samples

2025-10-03

Contents

0.1	Run Information	2
0.2	Assess Standard Curves	3
0.3	Assess Check Standards	4
0.4	Assess Blanks	5
0.5	Assess Duplicates	7
0.6	Calculate mmol/L concentrations & salinity, add dilutions	8
0.7	Assess Analytical Spikes	9
0.8	Check if samples within the range of the standard curve	10
0.9	Check to see if samples run match metadata & merge info	10
0.10	Visualize Data by Plot	12
0.11	Export Processed Data	12

##Add Required Packages

0.1 Run Information

```
##### Run information - PLEASE CHANGE
Date_Run = "2025-10-02" #Date that instrument was run
Run_by = "Zoe Read" #Instrument user
Script_run_by = "Zoe Read" #Code user
run_notes = "Need to rerun once the 2025 sample log is updated.
Std 1 is lower than the expected concentration.
" #any notes from the run
samples <- c("GCW", "GWI", "MSM", "SWH") #whatever identifies your samples within the same names
samples_pattern <- paste(samples, collapse = "|")
#samples_pattern <- "GCW" #use this instead of the line above if you have only one site code

##### File Names - PLEASE CHANGE
#file path and name for raw summary data file
raw_file_name_cl = "Raw Data/COMPASS_Synoptic_CB_MonMon_202507_Cl.txt"
raw_file_name_so4 = "Raw Data/COMPASS_Synoptic_CB_MonMon_202507_SO4.txt"

#file path and name of processed data file
processed_file_name = "Processed Data/COMPASS_SynopticCB_PW_Processed_Cl_SO4_202507.csv"

##### Log Files - PLEASE CHECK
#downloaded metadata csv - downloaded from Google drive as csv for this year
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2025.csv"

#qaqc log file path for this year
Log_path = "Raw Data/COMPASS_Synoptic_Cl_SO4_QAQClog_2024.csv"
```


##Set Up Code - constants and QAQC cutoffs

##Read in metadata and create similar sample IDs for matching to samples


##Import Sample Data

0.2 Assess Standard Curves


Chloride Std Curve



Sulfate Std Curve



```
## [1] "QAQC log file exists and has been read into the code."
```



```
## [1] "Cl Curve r2 GOOD"
```

```
## [1] "SO4 Curve r2 GOOD"
```


0.3 Assess Check Standards

```
## # A tibble: 5 x 5
##   sample_ID mean_Cl sd_Cl   cv_Cl flag_Cl
##   <chr>      <dbl> <dbl>   <dbl> <chr>
## 1 Standard 1    4.66 0.133 0.0284 Chloride Check Standard RSD within Range - P~
## 2 Standard 2    9.63 0.218 0.0226 Chloride Check Standard RSD within Range - P~
## 3 Standard 3   19.3 0.206 0.0107 Chloride Check Standard RSD within Range - P~
## 4 Standard 4   99.8 0.526 0.00527 Chloride Check Standard RSD within Range - P~
## 5 Standard 5  201. 0.353 0.00176 Chloride Check Standard RSD within Range - P~
```

```
## # A tibble: 5 x 5
##   sample_ID mean_S04 sd_S04   cv_S04 flag_S04
##   <chr>      <dbl> <dbl>   <dbl> <chr>
## 1 Standard 1    0.419 0.0113 0.0269 Sulfate Check Standard RSD within Range - ~
## 2 Standard 2    0.914 0.0166 0.0182 Sulfate Check Standard RSD within Range - ~
## 3 Standard 3    1.91 0.0238 0.0125 Sulfate Check Standard RSD within Range - ~
## 4 Standard 4    9.98 0.0516 0.00517 Sulfate Check Standard RSD within Range - ~
## 5 Standard 5   20.1 0.0485 0.00241 Sulfate Check Standard RSD within Range - ~
```

```
## [1] ">80% of Chloride Check Standards have RSD within range - PROCEED"
```

```
## [1] ">80% of Sulfate Check Standards have RSD within range - PROCEED"
```





```
## [1] ">80% of Chloride Check Standards are within range of expected concentration - PROCEED"
```

```
## [1] "<80% of Sulfate Check Standards are within range of expected concentration - REASSESS"
```

0.4 Assess Blanks

```
## [1] ">80% of Chloride Blank concentrations are lower 25% quartile of samples"
```

```
## [1] ">80% of Sulfate Blank concentrations are lower 25% quartile of samples"
```



Blank Conc <25% Quartile Samples ■ YE

Blank Conc <25% Quartile Samples ■ YE


```
## Chloride blanks mean ppm:
```

```
## [1] 0
```

```
## Sulfate blanks mean ppm:
```


```
## [1] 0
```

0.5 Assess Duplicates



```
## [1] ">80% of Chloride Duplicates have a CV <10% - PROCEED"
```

```
## [1] ">80% of Sulfate Duplicates have a CV <10% - PROCEED"
```



```
## [1] ">80% of Chloride Duplicates have a percent difference <10% - PROCEED"
```

```
## [1] ">80% of Sulfate Duplicates have a percent difference <10% - PROCEED"
```

0.6 Calculate mmol/L concentrations & salinity, add dilutions

```
# Convert ppm to mmol/L
all_dat$SO4_Conc_mM <- (all_dat$SO4_ppm / s_mw)
all_dat$Cl_Conc_mM <- (all_dat$Cl_ppm / cl_mw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# = ((1.807 * Cl_ppm) + 0.026) / 1000
all_dat$salinity <- ((1.8070 * all_dat$Cl_ppm) + 0.026) / 1000

#Need to determine dilution factors for your samples
#for Steph / COMPASS this depends on the site so...
all_dat$Dilution <- 1
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "UP"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "WC"),
```



```

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "MSM") & str_detect(all_dat$sample_ID, "SW"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "UP"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GCW") & str_detect(all_dat$sample_ID, "SW"),

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "UP"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "TR"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "WC"),
all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "GWI") & str_detect(all_dat$sample_ID, "SW"),


all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "SWH"), 50, all_dat$Dilution)

all_dat$Dilution <- ifelse(str_detect(all_dat$sample_ID, "Lucie"), 50, all_dat$Dilution)

# head(all_dat)

```

0.7 Assess Analytical Spikes



```
## [1] ">80% of S04 spikes have a recovery between the high and low cutoff - PROCEED"
```

0.8 Check if samples within the range of the standard curve

Sample Flagging




Table 1: SO4 samples

SO4_Conc_flag	Percent_samples
Within_Range	89.20863
bdl	10.79137

Table 2: Cl samples

Cl_Conc_flag	Percent_samples
Within_Range	98.561151
bdl	1.438849

0.9 Check to see if samples run match metadata & merge info


Some sample IDs are missing from metadata.

```
## [1] "GCW_202507_UP_LYSA_20CM"    "GCW_202507_UP_LYSB_10CM"
## [3] "GCW_202507_UP_LYSB_20CM"    "GCW_202507_TR_LYSA_20CM"
```

##	[5]	"GCW_202507_TR_LYSA_45CM"	"GCW_202507_TR_LYSB_10CM"
##	[7]	"GCW_202507_TR_LYSB_20CM"	"GCW_202507_TR_LYSB_45CM"
##	[9]	"GCW_202507_TR_LYSC_10CM"	"GCW_202507_TR_LYSC_20CM"
##	[11]	"GCW_202507_TR_LYSC_45CM"	"GCW_202507_WC_LYSA_10CM"
##	[13]	"GCW_202507_WC_LYSA_20CM"	"GCW_202507_WC_LYSA_45CM"
##	[15]	"GCW_202507_WC_LYSB_10CM"	"GCW_202507_WC_LYSB_20CM"
##	[17]	"GCW_202507_WC_LYSB_45CM"	"GCW_202507_WC_LYSC_10CM"
##	[19]	"GCW_202507_WC_LYSC_20CM"	"GCW_202507_WC_LYSC_45CM"
##	[21]	"GCW_202507_SW_A"	"GCW_202507_SW_B"
##	[23]	"GCW_202507_SW_C"	"MSM_202507_UP_LYSA_10CM"
##	[25]	"MSM_202507_UP_LYSA_20CM"	"MSM_202507_UP_LYSA_45CM"
##	[27]	"MSM_202507_UP_LYSB_10CM"	"MSM_202507_UP_LYSB_20CM"
##	[29]	"MSM_202507_UP_LYSB_45CM"	"MSM_202507_UP_LYSC_10CM"
##	[31]	"MSM_202507_UP_LYSC_20CM"	"MSM_202507_UP_LYSC_45CM"
##	[33]	"MSM_202507_TR_LYSA_10CM"	"MSM_202507_TR_LYSA_20CM"
##	[35]	"MSM_202507_TR_LYSA_45CM"	"MSM_202507_TR_LYSB_10CM"
##	[37]	"MSM_202507_TR_LYSB_20CM"	"MSM_202507_TR_LYSB_45CM"
##	[39]	"MSM_202507_TR_LYSC_20CM"	"MSM_202507_TR_LYSC_45CM"
##	[41]	"MSM_202507_WC_LYSA_10CM"	"MSM_202507_WC_LYSA_20CM"
##	[43]	"MSM_202507_WC_LYSA_45CM"	"MSM_202507_WC_LYSB_10CM"
##	[45]	"MSM_202507_WC_LYSB_20CM"	"MSM_202507_WC_LYSC_10CM"
##	[47]	"MSM_202507_WC_LYSC_20CM"	"MSM_202507_SW_A"
##	[49]	"MSM_202507_SW_B"	"MSM_202507_SW_C"
##	[51]	"GWI_202507_UP_LYSA_20CM"	"GWI_202507_UP_LYSA_45CM"
##	[53]	"GWI_202507_UP_LYSB_10CM"	"GWI_202507_UP_LYSB_20CM"
##	[55]	"GWI_202507_UP_LYSB_45CM"	"GWI_202507_UP_LYSC_20CM"
##	[57]	"GWI_202507_UP_LYSC_45CM"	"GWI_202507_TR_LYSA_10CM"
##	[59]	"GWI_202507_TR_LYSA_20CM"	"GWI_202507_TR_LYSB_10CM"
##	[61]	"GWI_202507_TR_LYSB_20CM"	"GWI_202507_TR_LYSB_45CM"
##	[63]	"GWI_202507_TR_LYSC_10CM"	"GWI_202507_TR_LYSC_20CM"
##	[65]	"GWI_202507_WC_LYSA_10CM"	"GWI_202507_WC_LYSB_10CM"
##	[67]	"GWI_202507_WC_LYSB_20CM"	"GWI_202507_WC_LYSB_45CM"
##	[69]	"GWI_202507_WC_LYSC_10CM"	"GWI_202507_WC_LYSA_45CM"
##	[71]	"GWI_202507_SW_A"	"GWI_202507_SW_B"
##	[73]	"GWI_202507_SW_C"	"SWH_202507_UP_LYSA_10CM"
##	[75]	"SWH_202507_UP_LYSA_20CM"	"SWH_202507_UP_LYSA_45CM"
##	[77]	"SWH_202507_UP_LYSB_10CM"	"SWH_202507_UP_LYSB_20CM"
##	[79]	"SWH_202507_UP_LYSB_45CM"	"SWH_202507_UP_LYSC_10CM"
##	[81]	"SWH_202507_UP_LYSC_20CM"	"SWH_202507_UP_LYSC_45CM"
##	[83]	"SWH_202507_SWAMP_LYSA_10CM"	"SWH_202507_SWAMP_LYSA_20CM"
##	[85]	"SWH_202507_SWAMP_LYSA_45CM"	"SWH_202507_SWAMP_LYSB_10CM"
##	[87]	"SWH_202507_SWAMP_LYSB_20CM"	"SWH_202507_SWAMP_LYSB_45CM"
##	[89]	"SWH_202507_SWAMP_LYSC_10CM"	"SWH_202507_SWAMP_LYSC_20CM"
##	[91]	"SWH_202507_SWAMP_LYSC_45CM"	"SWH_202507_TR_LYSA_10CM"
##	[93]	"SWH_202507_TR_LYSA_20CM"	"SWH_202507_TR_LYSA_45CM"
##	[95]	"SWH_202507_TR_LYSB_10CM"	"SWH_202507_TR_LYSB_20CM"
##	[97]	"SWH_202507_TR_LYSB_45CM"	"SWH_202507_TR_LYSC_10CM"
##	[99]	"SWH_202507_TR_LYSC_20CM"	"SWH_202507_TR_LYSC_45CM"
##	[101]	"SWH_202507_WC_LYSA_10CM"	"SWH_202507_WC_LYSA_20CM"
##	[103]	"SWH_202507_WC_LYSA_45CM"	"SWH_202507_WC_LYSB_10CM"
##	[105]	"SWH_202507_WC_LYSB_20CM"	"SWH_202507_WC_LYSB_45CM"
##	[107]	"SWH_202507_WC_LYSC_10CM"	"SWH_202507_WC_LYSC_20CM"
##	[109]	"SWH_202507_WC_LYSC_45CM"	"SWH_202507_SW_A"
##	[111]	"SWH_202507_SW_B"	"SWH_202507_SW_C"

```
## [113] "GCW_202507_TR_LYSA_10CM"    "9"
## [115] "10"                          "88"
## [117] "89"
```

0.10 Visualize Data by Plot



0.11 Export Processed Data

```
#end
```