

Synoptic CB: Porewater Sulfide

May 2026 Samples Ran on 20260526

2026-05-29

Run Information

###things that need to be changed

Date_Run = "20260526-20260527"

Month = "May"

Year = "2026"

Sample_Month_Year = "202605"

Run_by = "Zoe Read" #Instrument user

Script_run_by = "Zoe Read" #Code user

Project = "COMPASS"


Run_notes="Used 20260527_Plate2 std curve. 3/16 Matrix checks were bad. 3/8 dups were bad, probably because"

#Stds that should be excluded

stds_to_remove<-data.frame(Plate=c("Plate6"),IDs=c("Std 3"))

stds_to_remove<-NA

STD Curves



Checking STD Data against QAQC file








Table 1: Best std curve to use:

Date	Project	R2	Slope	Intercept	Top_STD	Plate
20260526-20260527	COMPASS	0.9895755	0.0125559	0.0792854	100	20260527_Plate2

Matrix Check QAQC




Sample triplicate means and sd dev before bad reps removed




Remove bad reps

Sample triplicate means and sd dev after bad reps removed




Check the dups for QAQC



```
## [1] ">60% of Duplicates are within <10%"
```

Check the spks for QAQC




```
## [1] "<60% of Spikes are out of range - REASSESS"
```

Check to see if samples run match metadata & merge info

```
## ***All sample IDs are present in metadata.***
```

Visualize Data by Plot



###END