

COMPASS Synoptic CB Porewater: DOC

September 2023 Samples

2026-06-03

##Setup - Change things here & write any notes

```
#identify section  
cat("Setup Information")
```

Setup Information

```
##### Run information - PLEASE CHANGE  
Date_Run = "09/16/23" #Date that instrument was run  
Run_by = "Stephanie J. Wilson" #Instrument user  
Script_run_by = "Stephanie J. Wilson" #Code user  
run_notes = "NA" #any notes from the run  
samples <- c("GCW", "GWI", "MSM", "SWH") #whatever identifies your samples within the same names  
samples_pattern <- paste(samples, collapse = "|")  
  #samples_pattern <- "GCW" #use this instead of the line above if you have only one site code  
chks_name = "Chk_Std_50C_2N" #what did you name your check standards?  
  
##### File Names - PLEASE CHANGE  
#file path and name for raw summary data file  
  #raw_file_name = "Raw Data/COMPASS_SynopticCB_PW_DOC_202505.txt" #example  
raw_file_name = c("Raw Data/TOCTN_COMPASS_Synoptic_TCTN_202309_1.txt",  
  "Raw Data/TOCTN_COMPASS_Synoptic_TCTN_202309_2.txt")  
  
#file path and name for raw all peaks file  
  #raw_allpeaks_name = "Raw Data/COMPASS_SynopticCB_PW_DOC_202505_allpeaks.txt" #example  
# raw_allpeaks_name = "Raw Data/COMPASS_SynopticCB_PW_DOC_202411_allpeaks.txt"  
  
#file path and name of processed data file  
  #processed_file_name = "Processed Data/COMPASS_SynopticCB_PW_Processed_DOC_202505.csv" #example  
processed_file_name = "Processed Data/COMPASS_SynopticCB_PW_Processed_DOC_202309.csv" #example  
  
##### Log Files - PLEASE CHECK  
#downloaded metadata csv - downloaded from Google drive as csv for this year  
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2023.csv"  
  
#qac log file path for this year  
Log_path = "Raw Data/COMPASS_Synoptic_TOCTN_QACLog_2024.csv"
```

##Set Up Code

##Read in metadata and create similar sample IDs for matching to samples

Import Data Functions

Import Sample Data

```
## Import Sample Data

## New names:
## New names:
## * ' ' -> '...14'

## # A tibble: 6 x 4
##   sample_name      npoc_raw tdn_raw run_datetime
##   <chr>          <dbl>   <dbl> <chr>
## 1 202309_GCW_TR_LysA_20cm 35.8    1.42  9/16/2023 7:07:41 PM
## 2 202309_GCW_TR_LysA_45cm 13.9    0.527 9/16/2023 7:32:14 PM
## 3 202309_GCW_TR_LysB_10cm  4.85   0.393 9/16/2023 7:55:43 PM
## 4 202309_GCW_TR_LysB_20cm  3.60   2.85  9/16/2023 8:15:13 PM
## 5 202309_GCW_TR_LysB_45cm 21.1    0.776 9/16/2023 8:40:45 PM
## 6 202309_GCW_TR_LysC_10cm 13.0    0.336 9/16/2023 9:01:01 PM
```

Assessing Standard Curves - assessed manually on the instrument


Assess Check Standards

```
## Assess the Check Standards

## New names:
## New names:
## * ' ' -> '...14'

## [1] "Carbon Check Standard RSD within Range"

## [1] "Nitrogen CHECK STANDARD RSD TOO HIGH - REASSESS"
```



```
## [1] ">60% of Carbon Check Standards are within range of expected concentration"
```

```
## [1] ">60% of Nitrogen Check Standards are within range of expected concentration"
```



Assess Blanks

```
## Assess Blanks
```

```
## New names:
## New names:
## * ' ' -> '...14'
```

```
## [1] ">60% of Carbon Blank concentrations are lower 25% quartile of samples"
```

```
## [1] ">60% of Nitrogen Blank concentrations are lower 25% quartile of samples"
```



Blank Conc <25% Quartile Samples ■ YE

Blank Conc <25% Quartile Samples ■ YE

```
## carbon blanks:
```

```
## [1] 0.2844268
```


```
## nitrogen blanks:
```

```
## [1] 0.01043415
```

Assess Duplicates

```
## Assess Duplicates
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once per session.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```




```
## [1] ">60% of Carbon Duplicates have a CV <10%"
```

```
## [1] ">60% of Nitrogen Duplicates have a CV <10%"
```

Sample Flagging - Are samples Within the range of the curve?

```
## Sample Flagging
```




Visualize Data by Plot


```
## Visualize Data
```

```
## Warning in rbind(c("202309", "GCW", "TR", "LysA", "20cm"), c("202309", "GCW", :
## number of columns of result is not a multiple of vector length (arg 17)
```

```
## [1] Date          Site_Code    Zone         Lysimeter    Depth
## [6] sample_name  npoc_raw    tdn_raw      run_datetime npoc_flag
## [11] tdn_flag
## <0 rows> (or 0-length row.names)
```

Samples: NPOC





Convert data from mg/L to uMoles/L

Check to see if samples run match metadata & merge info

```
## Check Sample IDs with Metadata

## Some sample IDs are missing from metadata.

## [1] "202309_MSM_UP_LysB_10cm" "202309_GWI_TR_LysC_45cm"
## [3] "202309_MSM_WC_LysA_45cm" "202309_MSM_WC_LysC_20cm"
## [5] "202309_SWH_WC_LysB_10cm" "202309_SWH_WC_LysC_10cm"
## [7] "202309_SWH_WC_LysC_20cm"
```

Export Processed Data

```
## Export Processed Data

## # A tibble: 6 x 21
##   Project      Region Site  Zone  Replicate Depth_cm Sample_ID  Year Month  Day
##   <chr>        <chr> <chr> <fct> <chr>      <int> <chr>    <int> <int> <int>
## 1 COMPASS: Sy~ CB    MSM  UP    <NA>      NA  202309_M~    NA    NA    NA
## 2 COMPASS: Sy~ CB    MSM  UP    C         10  202309_M~   2023    9    14
## 3 COMPASS: Sy~ CB    MSM  UP    C         20  202309_M~   2023    9    14
```



```

## 4 COMPASS: Sy~ CB      SWH  UP   A           10 202309_S~ 2023      9    20
## 5 COMPASS: Sy~ CB      SWH  UP   A           20 202309_S~ 2023      9    20
## 6 COMPASS: Sy~ CB      SWH  UP   A           45 202309_S~ 2023      9    20
## # i 11 more variables: Time <lgl>, Time_Zone <lgl>, npoc_mgL <dbl>,
## #   npoc_uM <dbl>, npoc_flag <chr>, tdn_mgL <dbl>, tdn_uM <dbl>,
## #   tdn_flag <chr>, Analysis_runtime <chr>, Run_notes <chr>, Field_notes <chr>

```

```

#end

```