

GENX Organs/Compass FeII/III Microplate Template

Name: Melanie
Date: 5/15/24

Plate ID: STD

	1	2	3	4	5	6	7	8	9	10	11	12
A	0.5 M HCl	0.6 uM Fe(II) H	3 uM Fe(II) G	6 uM Fe(II) F	9 uM Fe(II) E	18 uM Fe(II) D	30 uM Fe(II) C	45 uM Fe(II) B	60 uM Fe(II) A	100 (")	250 (")	500 (")
B	0.5 M HCl	0.6 uM Fe(II) H	3 uM Fe(II) G	6 uM Fe(II) F	9 uM Fe(II) E	18 uM Fe(II) D	30 uM Fe(II) C	45 uM Fe(II) B	60 uM Fe(II) A			
C	0.5 M HCl	0.6 uM Fe(II) H	3 uM Fe(II) G	6 uM Fe(II) F	9 uM Fe(II) E	18 uM Fe(II) D	30 uM Fe(II) C	45 uM Fe(II) B	60 uM Fe(II) A			
D	0.6 uM Fe(II) H	0.6 uM Fe(II) H	0.6 uM Fe(II) H	0.6 uM Fe(II) H	0.6 uM Fe(II) H	0.6 uM Fe(II) H	0.6 uM Fe(II) H	0.6 uM Fe(II) H	0.6 uM Fe(II) H	745 1:50		0 C
E	0.6 uM Fe(III) H	0.6 uM Fe(III) H	0.6 uM Fe(III) H	0.6 uM Fe(III) H	0.6 uM Fe(III) H	0.6 uM Fe(III) H	0.6 uM Fe(III) H	0.6 uM Fe(III) H	0.6 uM Fe(III) H			
F	0.5 M HCl	0.6 uM Fe(III) H	3 uM Fe(III) G	6 uM Fe(III) F	9 uM Fe(III) E	18 uM Fe(III) D	30 uM Fe(III) C	45 uM Fe(III) B	60 uM Fe(III) A	100 (")	250 (")	500 (")
G	0.5 M HCl	0.6 uM Fe(III) H	3 uM Fe(III) G	6 uM Fe(III) F	9 uM Fe(III) E	18 uM Fe(III) D	30 uM Fe(III) C	45 uM Fe(III) B	60 uM Fe(III) A			
H	0.5 M HCl	0.6 uM Fe(III) H	3 uM Fe(III) G	6 uM Fe(III) F	9 uM Fe(III) E	18 uM Fe(III) D	30 uM Fe(III) C	45 uM Fe(III) B	60 uM Fe(III) A			

Already run at 16 & 17 min

Ferrozine 1st Read Reducing Agent Buffer 2nd Read

GENX Organs/Compass FeI/III Microplate Template

57

Name: Melanie

Plate ID: Plate 1

Date: 5/15/24

	1	2	3	4	5	6	7	8	9	10	11	12
A	Fe(II) A	Fe(II) A	Fe(II) A	49		*	15 1:50			382 1:50		
B	Fe(II) E	Fe(II) E	Fe(II) E	864			81 1:50			596 1:50		
C	Fe(II) H	Fe(II) H	Fe(II) H	865			83 1:50			701 1:50		
D	0.5 M HCl	0.5 M HCl	0.5 M HCl	867			83 dup 1:50	dup	dup	707 1:50		
E	Fe(III) A	Fe(III) A	Fe(III) A	869			89 1:50			708 1:50		
F	Fe(III) E	Fe(III) E	Fe(III) E	880			89 1:50 spk	spk	spk	715 1:50		
G	Fe(III) H	Fe(III) H	Fe(III) H	5 1:50			276 1:50		*	716 1:50		
H	0.5 M HCl	0.5 M HCl	0.5 M HCl	12 1:50			182 1:50			804 1:50		

Ferrozine 1st Read Reducing Agent Buffer 2nd Read

1:16

10 mL sample: 150 mL DI

2 mL 1:50

40 mL sample 1960 mL DI

Spike?
Yes

58

GENX Organs/Compass FeII/III Microplate Template

Name: Melanie
 Date: 5/15/24

Plate ID: plate 2

	1	2	3	4	5	6	7	8	9	10	11	12
A	Fe(II) A	Fe(II) A	Fe(II) A	2			86			205		
B	Fe(II) E	Fe(II) E	Fe(II) E	10			87			208		
C	Fe(II) H	Fe(II) H	Fe(II) H	16		*	90			268		
D	0.5 M HCl	0.5 M HCl	0.5 M HCl	38			90 dup	dup	dup	269		
E	Fe(III) A	Fe(III) A	Fe(III) A	41		*	93		*	270		
F	Fe(III) E	Fe(III) E	Fe(III) E	44			93 spk	spk	spk	273		
G	Fe(III) H	Fe(III) H	Fe(III) H	78			190			317		0
H	0.5 M HCl	0.5 M HCl	0.5 M HCl	84		*	200		*	322		

☐ Ferrozine ☐ 1st Read ☐ Reducing Agent ☐ Buffer ☐ 2nd Read

1:50 Dilution

Spike?
Yes

2mL in vial 40mL sample: 1960mL DI

59

GENX Organs/Compass FeII/III Microplate Template

Name: MelaniePlate ID: Plate 3Date: 5/15/24

	1	2	3	4	5	6	7	8	9	10	11	12
A	Fe(II) A	Fe(II) A	Fe(II) A	324			586			697		
B	Fe(II) E	Fe(II) E	Fe(II) E	325			587			700		A
C	Fe(II) H	Fe(II) H	Fe(II) H	326			591			703		
D	0.5 M HCl	0.5 M HCl	0.5 M HCl	373			591 dup	dup	dup	718		
E	Fe(III) A	Fe(III) A	Fe(III) A	377			600			796		
F	Fe(III) E	Fe(III) E	Fe(III) E	388			600 spk	spk	spk	801		
G	Fe(III) H	Fe(III) H	Fe(III) H	393			601			806		
H	0.5 M HCl	0.5 M HCl	0.5 M HCl	490			692			809		

 Ferrozine 1st Read Reducing Agent Buffer 2nd Read

1:50

 spike?
 yes

GENX Organs/Compass FeII/III Microplate Template

60

Name: Melanie

Plate ID: Plate 4

Date: 5/15/24

	1	2	3	4	5	6	7	8	9	10	11	12
A	Fe(II) A	Fe(II) A	Fe(II) A	811			978			990		
B	Fe(II) E	Fe(II) E	Fe(II) E	816			985			991		
C	Fe(II) H	Fe(II) H	Fe(II) H	817			986			992		
D	0.5 M HCl	0.5 M HCl	0.5 M HCl	938			986 dup	dup	dup	1015		
E	Fe(III) A	Fe(III) A	Fe(III) A	939			987			1020		
F	Fe(III) E	Fe(III) E	Fe(III) E	974			987 spk	spk	spk	1021		
G	Fe(III) H	Fe(III) H	Fe(III) H	975			988			1027		
H	0.5 M HCl	0.5 M HCl	0.5 M HCl	976			989			1028		

Ferrozine 1st Read Reducing Agent Buffer 2nd Read

1:50

spike?
yes

GENX Organs/Compass FeII/III Microplate Template

Name: me lonie

Plate ID: plate 5

Date: 5/15/24

61

	1	2	3	4	5	6	7	8	9	10	11	12
A	Fe(II) A	Fe(II) A	Fe(II) A	425			433			533		
B	Fe(II) E	Fe(II) E	Fe(II) E	426		★	528			534		
C	Fe(II) H	Fe(II) H	Fe(II) H	427			529		★	535		
D	0.5 M HCl	0.5 M HCl	0.5 M HCl	428			529 dup	dup	dup	631		
E	Fe(III) A	Fe(III) A	Fe(III) A	429			530			632		★
F	Fe(III) E	Fe(III) E	Fe(III) E	430		★	530 spk	spk	spk	633		★
G	Fe(III) H	Fe(III) H	Fe(III) H	431			531			634		
H	0.5 M HCl	0.5 M HCl	0.5 M HCl	432			532			635		

☐ Ferrozine ☐ 1st Read ☐ Reducing Agent ☐ Buffer ☐ 2nd Read

spike?
yes

MSM Tr

1:50